

May 02, 2025

NU iPSC general culture

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgb98znlpk/v1>

De Xing¹, Ozlem Neyisci¹, Tao Bai¹, Mazhar Adli¹

¹Northwestern University

MorPhiC Consortium



Gina Kirsammer

Northwestern University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgb98znlpk/v1>

Protocol Citation: De Xing, Ozlem Neyisci, Tao Bai, Mazhar Adli 2025. NU iPSC general culture. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.81wgb98znlpk/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working.

Created: April 30, 2025

Last Modified: May 02, 2025

Protocol Integer ID: 204081

Keywords: substrate in stemflex media, stemflex media, cell

Funders Acknowledgements:

Molecular and Cellular Characterization of Essential Human Genes

Grant ID: HG012649

Abstract

KOLF2.2J iPS cells are grown on substrate in StemFlex media

Coating Plates with Substrate

1 **Option A:** Matrigel coating

- 1.1 Reconstitute aliquot of Matrigel (200ul) with cold medium (20mL, 1:100) and mix well. (Cold tips)
- 1.2 Immediately use Matrigel solution to coat 6-well plates (1mL for each well, 19 wells) & gently rock to spread evenly.
- 1.3 Place in 37 °C incubator for 30 minutes or room temperature 1 hour.

2 **Option B:** Synthemax coating

- 2.1 Prepare the required volume of a 1:40 dilution of Synthemax to coat the plates (2 ml water+50 µl 1 mg/ml Synthemax for each well of a 6-well plate; 8 ml water+200 µl 1 mg/ml Synthemax for each 10 cm dish).
- 2.2 Leave plates in hood for at least 2 h (or at 37 °C for at least 1 h).
- 2.3 Plates can be prepared ahead of time and stored dry at 4 °C for up to 3 months.

Thawing Cells

- 3 We generally thaw a vial of one million cells KOLF2.2J cells into one well of a 6-well plate.
- 3.1 Treat one well of a 6-well plate with Matrigel and stand at room temperature for 1 h.
- 3.2 Prepare 6 ml of StemFlex media containing ClonR2 (10x, 1:10 dilution) or rock inhibitor y-27632 (10 nM, 1:1000 dilution) in Falcon tube. Aspirate Matrigel and add 2 mL StemFlex media containing ClonR2 or rock inhibitor.



- 3.3 Remove 1 vial of cells in the 96-well Matrix plate or cryovial from the liquid nitrogen tank.
- 3.4 Hold the tube between your index finger and thumb to thaw quickly. The layer of mineral oil will thaw first. For cryovial, hold the top, put the bottom into 37 °C waterbath to thaw until a tiny ice left.
- 3.5 Once thawed, with a pipettor, transfer the cell suspension under the oil layer or the cell suspension from cryovial to a 15 mL Faclon tube containing 3 ml StemFlex+ ClonR2 or rock inhibitor media.
- 3.6 Centrifuge for 5 min at 300 x g to pellet the cells
- 3.7 Aspirate the media, taking care not to disturb the cell pellet and resuspend the cell pellet in 0.5 ml StemFlex + ClonR2 or rock inhibitor by gentle trituration.
- 3.8 Transfer the cell suspension to one well of a Matrigel-treated 6-well plate containing 2 ml StemFlex + ClonR2 or rock inhibitor and distribute evenly
- 3.9 Place 6-well plate in 37°C/5% CO2 incubator and rock plate gently to make sure cells are well distributed.
- 3.10 Next day, change the media to StemFlex (without ClonR2 or rock inhibitor) and then every 2 days thereafter.
- 3.11 When the plate is 80% confluent (3-4 days), freeze into 3 × 1 cryovials OR pass 1:15 to expand the cell line stock (See Passage and Expansion). (Note: We recommend passaging 1:15 to 3 × 10cm plates to allow the subsequent freezing of many cryovials to ensure future experiments can access an inventory of low passage stocks.)
- 3.12 The next day, replace the media in the well with StemFlex without ClonR2 or rock inhibitor and change the media every two days thereafter.

*Note: To ensure the cells are growing optimally, we recommend passing the cells at a dilution of 1:12 with ReLeSR at least once before initiating genome editing experiments.

Passaging Cells with ReLeSR



- 4 On reaching 80% confluence, passage the iPS cells at a dilution of 1:12 using ReLeSR as follows:
 - 4.1 Aspirate the culture media and wash the cells with 3 mL room temperature D-PBS (Life Technologies, 14190–094).
 - 4.2 Add 1 ml room temperature ReLeSR (Stem Cell Technologies, 05872) and incubate at room temperature in the hood for 5 min then aspirate ReLeSR. On Metrigel substrate, the iPS cells should not detach from the plate.
 - 4.3 Add 2 ml room temperature StemFlex media (3 ml per 6-well). Use a cell lifter (Corning, 3008) to gently scrape the cells from the surface of the plate. Break up the cells to small clumps by pipetting up and down 2–3 times using a 5ml pipette.
 - 4.4 Transfer 0.2 ml to one well of a new Metrigel -coated well of a 6-well plate containing 2 ml StemFlex media. Check under the microscope to make sure the cell clumps evenly distributed across the well. Place cells in a 37 °C/5% CO₂ incubator. Change the media every 2 days until the cells reach confluence, usually in 3 to 4 days.

*Note: When using ReLeSR, ROCK inhibitor is not required for passaging of small clumps of cells.

Freezing Cells

- 5 Generally, one well of 6-well plate will yield 3 vials of KOLF2.2J cells.
 - 5.1 Pre-warm ACCUTASE at room temperature.
 - 5.2 Pipette old media from the 6-well plate into sterile conical tube.
 - 5.3 Wash each well of a 6-well with 1 ml of PBS and then aspirate PBS.
 - 5.4 Add 1ml of ACCUTASE to each well of 6-well plate.

- 5.5 Incubate plate(s) for 7 minutes in 37°C incubator.
- 5.6 Cells on Metrigel-coated plate will detach, add 2 mL medium. Disperse the medium by pipetting over the cell layer surface several times. and transfer to sterile conical tube.
- Cells on Synthemax-coated plate will still adhere to plate. Aspirate ACCUTASE, add 2 mL medium and Using a cell lifter, thoroughly scrape plate to dislodge all cells. Gently pipette media up and down 3 times while washing plate surface with the media/cell mixture, and transfer to sterile conical tube.
- 5.7 Count cells using a haemocytometer (include a vital dye such as Trypan blue)
- 5.8 Centrifuge for 5 min at 300 x g to pellet the cells, and aspirate media.
- 5.9 Resuspend cells in appropriate volume of freezing media (90% KOSR + 10% DMSO) to a concentration of 1×10^6 live cells/ml of freezing media. Aliquot 1 mL freezing media/cell mixture per cryovial.
- 5.10 After capping tubes, put the cryovials into Cell Freezing Container (Corning, 432001 or Thermo Fisher Scientific, 15-350-50), and transfer to -80°C for overnight.
- 5.11 Transfer cryovials to liquid N₂ for long term storage.
- *Note: One cryovial should be thawed to test the viability of the frozen Lot. Expand it for 2-3 days and then subject it to pathogen testing. This will validate the quality of your frozen Lot for future experiments.