



Jul 08, 2023

Nuclear cytoplasmic fractionation

DOI

dx.doi.org/10.17504/protocols.io.ewov1oejylr2/v1

Narayana Yadavalli^{1,2,3,4,5}, Shawn M. Ferguson^{1,2,3,4,6,5}

¹Department of Cell Biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

³Program in Cellular Neuroscience, Neurodegeneration and Repair;

⁴Wu Tsai Institute Yale University School of Medicine, New Haven, Connecticut 06510, USA;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA;

⁶Kavli Institute for Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA



Berrak Ugur

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1oejylr2/v1>

Protocol Citation: Narayana Yadavalli, Shawn M. Ferguson 2023. Nuclear cytoplasmic fractionation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ewov1oejylr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 24, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 82369

Keywords: nuclear cytoplasmic fractionation, ASAPCRN, nuclear cytoplasmic fractionation this protocol, nuclear cytoplasmic fractionation

Funders Acknowledgements:

ASAP

Grant ID: 000580

Abstract

This protocol describes nuclear cytoplasmic fractionation.

Attachments



[724-1816.docx](#)

15KB



Materials

Reagents required














- PBS
- SDS
- Benzonase
- Hypotonic buffer
- 10 mM Hepes (pH 7.9)
- 10 mM KCl
- 0.1 mM EDTA
- 0.1 mM EGTA
- 1 mM dithiothreitol (DTT)
- high-salt buffer
- 20 mM HEPES
- 400 mM NaCl
- 1 mM EDTA
- 1 mM EGTA
- 1 mM DTT
- 0.5% NP-40

Troubleshooting



Nuclear cytoplasmic fractionation

6h 5m

- 1 Treat 500,000 BMDM's with DMSO or  50 nanomolar (nM) MLi-2 for  06:00:00 . 
- 2 After the treatment, wash the cells 3X in PBS. 
- 3 Harvest the cells in  500 μ L of ice-cold hypotonic buffer.
- 4 Homogenize the cells with 20 strokes of a Dounce homogenizer.
- 5 To  100 μ L of this homogenate, add SDS (1% final) and  25 U of Benzonase  (Novagen). Use this as total lysate.
- 6 Spin the rest of the homogenate at  14000 rpm, 4°C, 00:05:00 . 

- 7 Collect the supernatant (cytoplasmic fraction) into a new tube.
- 8 Resuspend the pellet (nuclear fraction) in  200 μ L of high-salt buffer and solubilize with SDS (1% final) in the presence of  25 U of Benzonase.
- 9 Measure protein concentration with the BCA reagent (Thermo Scientific), and subsequently analyze the samples by immunoblotting. 