

Jan 15, 2022

## Nuclear extraction from endometrial tumors for single nuclei sequencing

DOI

<https://dx.doi.org/10.17504/protocols.io.b3vfqn3n>


gracefoley<sup>1</sup>

<sup>1</sup>Northwestern University

Northwestern University



gracefoley



### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.b3vfqn3n>

**External link:** <https://www.nature.com/articles/s41591-020-0844-1#Sec13>

**Protocol Citation:** gracefoley 2022. Nuclear extraction from endometrial tumors for single nuclei sequencing . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.b3vfqn3n>



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** January 15, 2022

**Last Modified:** January 15, 2022

**Protocol Integer ID:** 56967

**Keywords:** nuclear extraction from endometrial tumor, endometrial tumor, endometrial adenocarcinoma, protocol for nuclear extraction, nuclear extraction, endometrioid type, single nuclei, nucleus rna, rna, frozen human tumor, sequencing

## Abstract

Modified from: Slyper, M. *et al.* A single-cell and single-nucleus RNA-Seq toolbox for fresh and frozen human tumors. *Nat. Med.* 2020 26**526**, 792–802 (2020).

A optimized protocol for nuclear extraction from endometrial tumors. Performed with endometrial adenocarcinoma, endometrioid type, FIGO grade 1.



Nuclear extraction optimization for ...



## Materials

Buffers:

### 2X ST Buffer

|  |                   | Stock  | Final  | Volume for 10 mL |
|--|-------------------|--------|--------|------------------|
|  | NaCl              | 5 M    | 292 mM | 584 uL           |
|  | Tris-HCl (pH 7.5) | 1 M    | 20 mM  | 200 uL           |
|  | CaCl <sub>2</sub> | 1M     | 2 mM   | 20 uL            |
|  | MgCl <sub>2</sub> | 100 mM | 42 mM  | 4.2 mL           |
|  | Ultrapure Water   | -      | -      | 5 mL             |

### NST Buffer

|  |                        | Stock | Final | Volume for 10 mL |
|--|------------------------|-------|-------|------------------|
|  | Nonidet P40 Substitute | 10%   |       | 200 ul           |
|  | BSA                    | 5%    |       | 20 ul            |
|  | ST Buffer              | 2X    | 1X    | 5 ml             |
|  | Nuclease Free Water    | -     | -     | 4.75 mL          |

### Nuclei Resuspension Buffer

|  |                 | Stock concentration | Final    | Volume Needed for 50 mL |
|--|-----------------|---------------------|----------|-------------------------|
|  | BSA solution    | 5%                  | 1%       | 10 ml                   |
|  | RNase Inhibitor | 40 U/ul             | 0.2 U/ul | 0.25 ml                 |
|  | 1X PBS          | -                   | -        | 39.75 ml                |





- 1 Place frozen tissue in a well of a 6 well plate containing 1 mL of cold NST disassociation buffer.
- 1.1 Ensure plate is on ice
- 2 Chop tissue using a scalpel until mostly homogenized 2m
- 3 Incubate solution on ice for 5 minutes 5m
- 4 Filter the homogenized solution through a 40 uM falcon cell strainer into a 50 mL falcon tube
- 5 Use an additional 1 mL of NST buffer to rinse the well and filter
- 6 Count nuclei
- 7 Centrifuge nuclei at 300 g for 10 minutes and aspirate supernatant completely 10m
- 8 Resuspend pellet and bring the volume up to 5 mL using ST buffer
- 9 Transfer solution to 15 mL conical tube and centrifuge for 10 minutes at 300 g at 4oC 10m
- 10 10. Resuspend pellets on ice in 1 ml of the ST buffer
- 11 11. Filter through 35 um falcon strainer
- 12 12. Count nuclei



13 Centrifuge nuclei at 300 g for 10 minutes and aspirate supernatant completely

10m

14 13. Resuspend nuclei pellet in resuspension buffer to get a concentration of 1,000 nuclei per  $\mu\text{L}$

14.1 Start with .5 mL per 0.2 grams of endometrial tumor tissue

14.2 Look at nuclear membranes at 40x-60x to ensure high nuclear membrane quality with minimal blebbing