



Mar 15, 2024

## Nuclear Extraction from Tissue for FACS

DOI

[dx.doi.org/10.17504/protocols.io.5jyl8j678g2w/v1](https://dx.doi.org/10.17504/protocols.io.5jyl8j678g2w/v1)

anita.adami<sup>1</sup>

<sup>1</sup>Laboratory of Molecular Neurogenetics, Department of Experimental Medical Science, Wallenberg Neuroscience Center and Lund Stem Cell Center, BMC A11, Lund University, 221 84 Lund, Sweden.



**Raquel Garza**

Lund University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8j678g2w/v1>

**Protocol Citation:** anita.adami 2024. Nuclear Extraction from Tissue for FACS. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.5jyl8j678g2w/v1>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working



**Created:** February 07, 2023

**Last Modified:** May 31, 2024

**Protocol Integer ID:** 76528

**Keywords:** ASAPCRN, tissue preparation for nuclear extraction, nuclear extraction from tissue, nuclear extraction, tissue preparation, facs this protocol, fac, tissue, preparation

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research**

Grant ID: ASAP-000520

**Swedish Research Council**

Grant ID: 2018-02694

**Swedish Brain Foundation**

Grant ID: FO2019-0098

**Cancerfonden**

Grant ID: 190326

**Barncancerfonden**

Grant ID: PR2017-0053

**NIHR Cambridge Biomedical Research Centre**

Grant ID: NIHR203312

**Swedish Society for Medical Research**

Grant ID: S19-0100

**National Institutes of Health**

Grant ID: HG002385

**Swedish Research Council**

Grant ID: 2021-03494

**Swedish Research Council**

Grant ID: 2020-01660

## Abstract

This protocol details the process for pre-tissue preparation and tissue preparation for nuclear extraction for FACS.

## Troubleshooting

## Before start

Keep everything On ice.



## Pre-Tissue Preparation

- 1 Coat all tubes (centrifugation, sample and collection tubes) in coating medium

Overnight ( 1 % (v/v) BSA , filtered).



- 2 Cool ultracentrifuge to 4 °C prior to centrifugation

- 3 Pre-weigh the DTT needed.

- 4 Keep douncers On ice .

## Tissue Preparation

**10m**

- 5 Put 2.8 mL sucrose solution in centrifuge tube.

- 6 Mix 6 mL sucrose + 2 mL lysis solution in 15 ml falcon tube.



- 7 If the tissue is large (>2 mm<sup>3</sup>), chop the tissue with a razor.

- 8 Collect chopped tissue with 0.5 mL lysis buffer and transfer to douncer. Repeat with another 0.5 mL lysis buffer to get all tissue.


- 9 Dounce: 5-7x loose, 5-7x tight.


- 10 Pour douncer content into falcon tube with sucrose/lysis buffer mixture and flip upside down a few times. Pour 1 mL of tube content back and forth from douncer to collect remaining nuclei.



- 11 Carefully layer the 8.5 mL mix of lysis/sucrose/nuclei to the centrifuge tube, on top of sucrose.



12 Balance the tubes ( $\leq$   5 mg difference).


13 Place tubes in ultracentrifuge, spin at  15500 rpm , 02:15:00, (rotor dependent, use 30.000xg) .




14 Carefully aspirate off supernatant without disturbing the pellet.




#### Note

Due to the viscosity, try to aspirate the inside of the tube to get rid of all supernatant.



15 Soften pellet for  00:10:00 in  50  $\mu$ L medium A .

10m

16 Add  100  $\mu$ L dilution buffer (containing a nuclei marker, such as Draq7) and filter using 70um pipette tip filters into coated tubes.

17 FACS sort, low flowrate, largest nozzle, cooling system on. The sucrose will cause problems at FACS. Avoid this by adding dilution buffer (with your nuclei marker),  100  $\mu$ L at the time (usually need  300  $\mu$ L -  400  $\mu$ L for dilution to avoid problems – add that even before starting to sort).

18 For single nuclei sequence using 10x: Pre-coat the pipette tips needed to load the 10x machine according to step 1.

19 For RNA extraction: Pellet nuclei at  1300 x g , 4°C, 00:15:00 . Remove supernatant and snap freeze for later RNA extraction (or resuspend in  350  $\mu$ L RLT including beta-mercapto ethanol directly. Use RNeasy micro columns. For very small volumes, resuspend in RLT without pelleting).

