


Jul 07, 2023

Version 2

# Nuclei Isolation and Immunoprecipitation for 10X Sequencing V.2

DOI

[dx.doi.org/10.17504/protocols.io.kxygx9jk4g8j/v2](https://dx.doi.org/10.17504/protocols.io.kxygx9jk4g8j/v2)



Lakme Caceres<sup>1</sup>

<sup>1</sup>Princeton Neuroscience Institute



Lakme Caceres

Princeton University

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx9jk4g8j/v2>

**Protocol Citation:** Lakme Caceres 2023. Nuclei Isolation and Immunoprecipitation for 10X Sequencing . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx9jk4g8j/v2> Version created by **Lakme Caceres**

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** July 06, 2023

**Last Modified:** July 07, 2023

**Protocol Integer ID:** 84603

**Keywords:** isolating nuclei, nuclei isolation, sequencing application, immunoprecipitation for 10x, immunoprecipitation, isolation

## Abstract

This protocol is for isolating nuclei for downstream sequencing applications.

## Guidelines

Keep tissue/nuclei on ice as much as possible.

## Materials

Dynabeads - 4C

anti-GFP - 4C (roche)

Triton X 100 - 4C

DTT - 4C

RNAsin -20C

## Troubleshooting



## Prepare Stock Solutions

- 1 Make 20 mL **10% BSA** by combining 2 mL of BSA with 18 mL of MilliQ water in a 50 mL falcon tube. (4°C - 2 weeks)
- 2 Make 20 mL **10% Triton X-100** by combining 18 mL MilliQ water with 2 mL Triton X-100 in a 50 mL tube. Vortex and then incubate at room temperature for 20 minutes. Filter it through a 0.22 µm filter with a syringe into a clean 50 mL tube. (4°C - 1 month)
- 3 Make 250 mL **Nuclear Isolation Media** by filling a 250 mL bottle with 200 mL of MilliQ water and then adding 2.5 mL 1M Tris, 6.26 mL 1M KCl, 1.25 mL 1M MgCl<sub>2</sub>, and 21.45 g Sucrose. Shake until sucrose is dissolved then fill to 250 mL with MilliQ water. (4°C - 2 weeks)
- 4 Make 50 mL **Citrate-Phosphate Buffer** (pH 5) by adding 0.48 g Citric Acid and 0.91 g Dibasic Sodium Phosphate Dihydrate to a falcon tube and then filling it to 50 mL with MilliQ water. Titrate the pH with NaOH using the pH meter. (4°C - 1 month)

## Immunoprecipitation Prep

- 5 Vortex stock of Dynabeads Protein G and then add 5 µL to a 1.5 mL eppendorf. Place it on the mag rack for a minute and discard the supernatant remaining at the bottom of the tube without disturbing the beads.
- 6 Remove the tube from the magnet and resuspend the beads with 500 µL of Citrate Phosphate Buffer. Place the tube back on the mag rack for a minute and discard the supernatant. Repeat this step once more.
- 7 Add 100 µL of 1X PBS and 2 µL anti-GFP to the washed Dynabeads. Add this to the rotator in the 4°C fridge and let it incubate for 1 hour. Proceed with nuclei isolation while the beads incubate.

## Prepare Fresh Solutions

- 8 Make 3 mL **Homogenization Buffer** per sample by adding 2.895 mL Nuclear Isolation Media (filtered via syringe) to a 5 mL eppendorf. Then add 3 µL 100 mM DTT and 30 µL 10% Triton X-100. Add 15 µL RNasin and invert to mix. Store on ice.
- 9 Make 200 µL **Blocking Buffer** per sample by dividing your total desired volume of blocking buffer by 10 to get the amount of 10% BSA in µL. Add this amount to a tube and then fill the remainder with 1X PBS.

## Homogenization



- 10 Clean dounce, scalpel, and forceps using MilliQ water, ethanol, RNase Zap, then MilliQ again. The red-tape forceps are for unfixed tissue.
- 11 Get tissue sample from -80°C freezer and place on dry ice. Weigh it on a sterile, tared weigh boat.
- 12 Add tissue to dounce and push it to the bottom using 1 mL of Homogenization Buffer and the pestle. Homogenize the tissue without creating bubbles. Then add the remaining 2 mL of the Homogenization Buffer and continue to dounce until homogenized.
- 13 Pass all of the nuclei suspension through three FlowMi filters, 1000 µL at a time into a new 5 mL eppendorf.
- 14 Centrifuge at 900 g/rcf for 10 minutes at 4°C.
- 15 Discard the supernatant and resuspend the pellet in 210 µL Blocking Buffer. Incubate for 10 minutes on ice.

## Cell Count

- 16 Add 9 uL of sample to a PCR tube and then add 1 uL of acridine orange.  
If sample is very concentrated, instead add 2 uL sample to 2 uL of acridine orange and 16 uL 1X PBS.
- 17 Pipette mix and then add the total volume to a three-chamber cell counting chip and make note of the channels used (A, B, and/or C).
- 18 On the cell counter, select Fluorescence Cell Counting → Cell Lines → Advanced → Protocol → and then choose "NUCLEI" from the list of protocols. Load the protocol.  
  
Then go to Settings and choose the appropriate number of channels.  
  
Then hit "Count" and then "Start Count."
- 19 When the cell count is complete, you will get a reading in cells/mL. Convert this to cells/uL by dividing this number by 1,000.
- 20 Save approximately 40,000 nuclei from the original sample to use as our unpurified population for sequencing. (If the concentration was 1,000 cells/uL, then save 40 uL).



- 21 If there are channels left on the cell counter chip, mark the used channels on the back and place it back in the drawer for future use.

## Immunoprecipitation

- 22 After the one hour incubation of the beads and the aliquoting of nuclei for the unpurified population, place the incubated dynabeads tube on the 1.5mL magnet for 2 minutes and discard the supernatant without disturbing the beads.
- 23 Wash with 100 uL 1X PBS Buffer, enough to submerge the beads while they are on the rack. Discard the supernatant.
- 24 Add the total volume of nuclei sample to the beads. Incubate on the rotator in the 4°C fridge for 1 hour.
- 25 After incubation, place the samples on the mag rack for 2 minutes. Collect the supernatant in appropriately labeled eppendorf tubes. This will be our supernatant population for sequencing.
- 26 Wash the bead tubes with enough 1X PBS to submerge the beads (don't resuspend) and then discard the supernatant.
- 27 Remove the tubes from the rack and then resuspend the beads in 50 uL 1X PBS.



## Cell Count

- 28 Take another cell count of the supernatant and bead-bound samples. Reverse pipette to record their exact volumes.
- 29 Dilute all of the samples (unpurified, supernatant, and bead-bound), if necessary, to 1,000 cells/uL. The bead-bound population likely won't need dilution.
- 30 Proceed with 10X RNAseq or ATACseq protocols.

## Protocol references

Adapted from the Allen and Thermofisher protocols.