



Nov 12, 2024

Version 2

Nuclei Isolation and Sorting from Frozen Human Temporal Cortex_V2 V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk4k55g5r/v2>

Nuclei isolation and
sorting from
frozen human
temporal cortex

Clemens Scherzer¹, Beatrice Weykopf¹

¹Stephen & Denise Adams Center for Parkinson's Disease Research at Yale School of Medicine



Beatrice Weykopf

Yale

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk4k55g5r/v2>

Protocol Citation: Clemens Scherzer, Beatrice Weykopf 2024. Nuclei Isolation and Sorting from Frozen Human Temporal Cortex_V2. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlk4k55g5r/v2> Version created by **Beatrice Weykopf**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 31, 2024

Last Modified: November 12, 2024

Protocol Integer ID: 106795

Keywords: nuclei isolation, frozen human temporal cortex, ASAPCRN, sorting from frozen human temporal cortex-v2, frozen human temporal cortex-v2, frozen human temporal cortex, nuclei isolation, nuclei, sorting, isolation, protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000301

Abstract

This protocol is about nuclei isolation and sorting from frozen human temporal cortex.

Guidelines

This protocol need prior approval by the user's institutional review board (IRB) or equivalent ethics committee.



Materials

1.  Dounce tissue grinder set **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8938-1SET**
2.  Nuclei Isolation Kit: Nuclei PURE Prep **Merck MilliporeSigma (Sigma-Aldrich) Catalog #NUC201-1KT**
3.  DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #43816-10ML**
4.  DPBS with no calcium and magnesium **Thermo Fisher Scientific Catalog #14190-144**
5.  UltraPure™ BSA (50 mg/mL) **Thermo Fisher Scientific Catalog # AM2616**
6.  Recombinant RNase Inhibitor **Takara Bio Inc. Catalog #2313B**
7.  70 µm Sterile Cell Strainer **Fisher Scientific Catalog #22363548**
8.  Corning™ Falcon™ Test Tube with 35µm Cell Strainer Snap Cap **Corning Catalog #352235**
9.  ART™ Wide Bore Filtered Pipette Tips **Thermo Fisher Scientific Catalog #2079G**
10.  Flowmi™ Cell Strainer 40 µm **Bel-Art Catalog #H13680-0040**
11.  DAPI **Thermo Fisher Scientific Catalog #D1306**

Buffer Recipes:

1. Lysis buffer (LB):
 - Nuclei PURE lysis Buffer
 - 1:1000 DTT
 - 0.01% Triton-x-100
2. Nuclei wash and resuspension buffer (NWRB):
 - 1x DPBS
 - 0.01% BSA
 - freshly added RNase Inhibitor (40 U/ml)
3. Nuclei wash and resuspension buffer with DAPI (NWRBD):
 - 1x NWRB
 - 0.01 mg/ml DAPI

Safety warnings

 For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

Before start

Prepare Buffers as described in section '[Materials](#)'.



Method

15m

- 1 Cut at total of 10 sections of human temporal cortex (50 μ m thickness) store on dry ice.
- 2 Transfer tissue to a glass dounce tissue grinder on regular ice.
- 3 Add  2 mL LB to the tissue grinder.
- 4 Homogenize tissue with pastel A *10 times* and pastel B *10 times*.
- 5 Transfer the homogenate into a 2 ml microcentrifuge tube.
- 6 Incubate  On ice for  00:05:00 , mix with Wide bore tips once during the incubation.
- 7 Filter the homogenate with a 70 μ m strainer.
- 8 Centrifuge at  500 x g, 4°C, 00:05:00 .
- 9 Remove the supernatant and resuspend the pellet in  1.5 mL LB with Wide bore tips.
- 10 Incubate  On ice for  00:05:00 .
- 11 Centrifuge at  500 x g, 4°C, 00:05:00 .
- 12 Remove the supernatant and add  1500 μ L NWRB to the pellet without resuspension

5m



5m





- 13 Centrifuge at  500 x g, 4°C, 00:05:00 . 5m
- 14 Remove the supernatant and add  500 µL NWRB to the pellet without resuspension.
- 15 Add  1.0 mL additional NWRB to the tube and resuspend the pellet.
- 16 Centrifuge at  500 x g, 4°C, 00:05:00 . 
- 17 Remove the supernatant and resuspend the pellet in  1.8 mL NWRBD . and filter through 40µm FLOWMI cell strainer (BIO-RAD)
- 18 Next, DAPI-stained nuclei were sorted using FACSAria III cell sorter (BD Bioscience)
- 19 Gating strategy: 
 - use 100 µm nozzle
 - FSC-A at 5V / SSC-A at 180V gating
 - second gating is DAPI-positive nuclei using violet laser (405nm) with 450/50 bandpass filter
 - Collect approximately 90,000 events in 1.5ml eppendorf tube and count on hemocytometer for quality control
- 20 Immediately after collection centrifuge at  500 x g, 4°C, 00:05:00 . 
- 21 Remove supernatant and resuspend pellet according to the input amount in NWRB, For 90,000 events you resuspend in 15µL NWRB if you have less number of nuclei adjust the total volume accordingly.
- 22 Use  8.4 µL of nuclei (50,400 events) for 10x GEM generation and barcoding following manufacturer's protocol. Continue with GEM generation as per manufacturer's instructions.

Protocol references

Habib et al 2017; pmid-28846088