

Jan 07, 2020

## Nuclei Isolation for SNARE-seq2

DOI

[dx.doi.org/10.17504/protocols.io.8tvhwn6](https://dx.doi.org/10.17504/protocols.io.8tvhwn6)

Carter Palmer<sup>1</sup>, Christine Liu<sup>1</sup>, Jerold Chun<sup>2</sup>

<sup>1</sup>University of California, San Diego, and Sanford Burnham Prebys Medical Discovery Institute;

<sup>2</sup>Sanford Burnham Prebys Medical Discovery Institute

BICCN / BICAN



Carter Palmer

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.8tvhwn6](https://dx.doi.org/10.17504/protocols.io.8tvhwn6)

**Protocol Citation:** Carter Palmer, Christine Liu, Jerold Chun 2020. Nuclei Isolation for SNARE-seq2 . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8tvhwn6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 28, 2019

**Last Modified:** January 07, 2020

**Protocol Integer ID:** 29269

**Keywords:** Brain, Nuclei, Isolation, Single-Cell



## Abstract

This protocol is intended to be used for the isolation of nuclei from fresh-frozen brain tissue in preparation for analysis by Single-Nucleus Chromatin Accessibility and mRNA Expression sequencing (SNARE-seq). It has been applied to tissues from mouse, marmoset, and human.

## Guidelines

This protocol is designed specifically for isolating tissue for SNARE-seq. RNA stability is considered at every step. Samples are kept on ice throughout the process, all centrifuges are pre-chilled to 4 °C before use, and RNase Away is used to spray down all surfaces and pipets before use.

## Materials

### MATERIALS

- ☒ Sucrose **Fisher Scientific Catalog #S3-212**
- ☒ Protease Inhibitor Tablets cOmplete Mini EDTA free **Roche Catalog #11836170001**
- ☒ RNase Inhibitor **Takara Bio Inc. Catalog #2313A**
- ☒ Bovine Serum Albumin **Gemini Bio-Products Catalog #700-107P**
- ☒ Magnesium acetate tetrahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M5661**
- ☒ Calcium chloride dihydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C5080**
- ☒ Ethylenediaminetetraacetic acid (EDTA) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #EDS**
- ☒ Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787**
- ☒ 16% Formaldehyde (w/v) Methanol-free **Thermo Fisher Scientific Catalog #28906**
- ☒ DAPI **Merck MilliporeSigma (Sigma-Aldrich) Catalog #10236276001**
- ☒ 50um filters **Sysmex Catalog #04-0042-2317**
- ☒ Tissue Homogenizer **Catalog #358005**

