



Jun 22, 2023

🌐 Nuclei isolation from fresh frozen human colon tissue for 10X Genomics Multi-omics (ATAC + GEX) assay

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8pn17g2w/v1>

Antonina Mikorska¹, Koen Theunis^{2,3}, Suresh Poovathingal^{2,3}, sam kint¹, Sarah Geurs¹, Thierry Voet¹

¹Laboratory of Reproductive Genomics, KU Leuven; ²VIB-KU Leuven Center for Brain & Disease Research;

³Laboratory of Computational Biology, KU Leuven



Antonina Mikorska

Laboratory of Reproductive Genomics, KU Leuven

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8pn17g2w/v1>

Protocol Citation: Antonina Mikorska, Koen Theunis, Suresh Poovathingal, sam kint, Sarah Geurs, Thierry Voet 2023. Nuclei isolation from fresh frozen human colon tissue for 10X Genomics Multi-omics (ATAC + GEX) assay. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5jyl8pn17g2w/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: In development

It works, but we are still developing and optimizing this protocol

Created: June 22, 2023

Last Modified: June 01, 2024

Protocol Integer ID: 83804

Keywords: ASAPCRN, fresh frozen human colon tissue, nuclei isolation protocol, 10x genomic, atac

Funders Acknowledgements:

Aligning Science Across Parkinson's

Disclaimer

This protocol is still being optimized! At the moment, we are testing if we can achieve better results without using the collagenase/dispase as dispase is a proteinase.

The buffers are based on the user guides from 10X Genomics.

Abstract

Nuclei isolation protocol used for sn multi-omics assay on fresh frozen human colon tissue (full thickness and muscularis layer only).

Materials

BUFFERS

Buffers based on 10X Genomics user guides.

Collagenase/Dispase 10269638001 from Sigma

	A	B	C	D	E	F
	Collagenase Lysis Buffer	Stoc k con c		Fin al con c		V (µl)
	Tris-HCl pH 7.4	100 0	mM	10	mM	10,5
	NaCl	500 0	mM	10, 0	mM	2,1
	MgCl ₂	100 0	mM	3	mM	3,2
	DTT	100 0	mM	1	mM	1,1
	RNAseIn	40	U/uL	1	U/uL	26,3
	Dispase/collagen ase	0,1	mg/m L	0,0 03	mg/m L	31,5
	NF H ₂ O					975,5
	Total Volume					1050,0
	Wash Buffer 1	Stoc k con c		Fin al con c		V (µl)
	Tris-HCl (pH 7.4)	100 0	mM	10	mM	28,4
	NaCl	500 0	mM	10	mM	5,7
	MgCl ₂	100 0	mM	3	mM	8,5
	BSA	10	%	1	%	283,5
	Tween-20	10	%	0,1	%	28,4



	A	B	C	D	E	F
	DTT	100 0	mM	1	mM	2,8
	RNase in plus (promega)	40	U/ul	1	U/ul	70,9
	Water					2406,9
	Total Volume					2835,0
	0.1x Lysis	Stoc k con c		Fin al con c		V (µl)
	Tris-HCl (pH 7.4)	100 0	mM	10	mM	5,3
	NaCl	500 0	mM	10	mM	1,1
	MgCl ₂	100 0	mM	3	mM	1,6
	Tween-20	10	%	0,0 1	%	0,5
	NP-40	10	%	0,0 1	%	0,5
	Digitonin	0,5	%	0,0 02	%	2,1
	BSA	10	%	1	%	52,5
	DTT	100 0	mM	1	mM	0,5
	Rnase Inhibitor Protector (SIGMA)	40	U/ul	1	U/ul	13,1
	Water					447,8
	Total Volume					525,0
	Wash Buffer 2	Stoc k con c		Fin al con c		V (µl)
	Tris-HCl (pH 7.4)	100 0	mM	10	mM	10,5

	A	B	C	D	E	F
	NaCl	500 0	mM	10	mM	2,1
	MgCl ₂	100 0	mM	3	mM	3,2
	BSA	10	%	1	%	105,0
	Tween-20	10	%	0,1	%	10,5
	DTT	100 0	mM	1	mM	1,1
	Rnase Inhibitor Protector (SIGMA)	40	U/ul	1	U/ul	26,3
	Water					891,5
	Total Volume					1050,0
	1X Diluted Nuclei buffer	Stoc k con c		Fin al con c		V(μl)
	Nuclei buffer	20	X	1	X	10
	DTT	100	mM	1	mM	2
	Rnase Inhibitor Protector (SIGMA)	40	U/ul	1	U/ul	5
	Water					183,0
	Total Volume					200

Instruments and equipment:

- Kai Medical 2 mm Biopsy punch BP-20F
- WPI Noyes scissors 12 cm S/S 15 mm blades
- VWR®, Disposable Pestles and Cordless Motor for Pellet Mixing
- PluriSelect 40 μm, 20 μm filters
- Flowmi 40 μm Cell strainer for 1000P
- Luna FL Cell Counter
- Swing out rotor centrifuge



Sample preparation

- 1 Place the sample preparation instruments on dry ice: biopsy punch, 1.5 mL EP tube, tweezers, petri dish.
- 2 Using 2 mm biopsy punch, aliquot 1.5-2 tissue pieces per patient. Place the pieces in a fresh 1.5 mL EP tube placed on dry ice. Repeat with a fresh biopsy punch for each multiplexed patient.

Experiment preparation

- 3 Prepare all buffers fresh and on wet ice. Add the the detergents, DTT and RNase inhibitor just before use.
- 4 Place all the plastics and filters on wet ice.

Tissue homogenization

- 5 **Perform all steps on wet ice.** Add 500 μ L Collagenase Lysis Buffer to the tube with multiplexed patient samples. Let thaw slightly.
- 6 Cut the tissue with scissors until there are no more pieces visible. **Time** = around 4 min (depending on the amount and the characteristics of the tissue).
 - 6.1 After cutting, place the scissors in a fresh 1.5 mL EP tube placed on ice.
- 7 Further homogenize the tissue using an automatic plastic pellet pestle, for 45 s.
 - 7.1 Place the pellet pestle in the EP tube instead of the scissors.
- 8 Mince the tissue for 1 min.
- 9 Wash the scissors using 250 μ L Collagenase Lysis Buffer (while adding it to the tissue).



- 10 Homogenize the tissue with the used pellet pestle x15 (by hand).
- 11 Wash the pestle with 250µL Collagenase Lysis Buffer and add it to the tissue.
- 12 Gently pipet up and down 10X with P1000.
- 13 Centrifugate for 5 min at 500g at 4°C, remove the supernatant and resuspend in 1 mL Wash Buffer 1.
- 14 With P1000, transfer the sample onto 40 and 20 µm stacked strainers placed on 25 mL EP tube.
- 15 Wash the 1.5 mL EP tube with additional 1 mL Wash buffer 1, transfer onto the top filter.
- 16 Centrifugate the 25 mL EP tube for 5 min at 500g, 4°C.
 - 16.1 Gently remove the supernatant.
- 17 Resuspend the pellet in 700 µL of Wash Buffer 1 and filtrate through Flowmi filter into a fresh DNA LoBind 1.5 mL EP tube.
- 18 Centrifugate for 5 min at 500g at 4°C
 - 18.1 Gently remove the supernatant.

Nuclei permeabilization

- 19 Resuspend the pellet in 500 µL 0.1x Lysis Buffer, pipette mix X5 with P1000.



- 20 Incubate on ice for 1 min (tissue quality dependent).
- 21 Add 1000 μ L of Wash Buffer 2 and pipette mix 5X with P1000.
- 22 Centrifugate for 5 min at 500 g at 5°C.
- 22.1 Very gently remove the supernatant.

Nuclei resuspension and counting

- 23 Resuspend the pellet in 25-50 μ L 1X Diluted Nuclei Buffer.
- 24 Cell count the nuclei and proceed to the 10X Genomics Chromium Next GEM Single Cell Multiome ATAC + Gene Expression assay.