

Sep 14, 2020

Version 2

# Nuclei isolation from human intestinal biopsic tissue for single-cell genomic applications V.2

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.bma2k2ge](https://dx.doi.org/10.17504/protocols.io.bma2k2ge)

Ran RZ Zhou<sup>1</sup>, Oni Basu<sup>1</sup>

<sup>1</sup>University of Chicago


Human Cell Atlas Metho...

Helmsley project\_Basu lab



Ran RZ Zhou

University of Chicago



## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bma2k2ge>

**Protocol Citation:** Ran RZ Zhou, Oni Basu 2020. Nuclei isolation from human intestinal biopsic tissue for single-cell genomic applications. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bma2k2ge> Version created by **Ran RZ Zhou**



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** September 11, 2020

**Last Modified:** September 14, 2020

**Protocol Integer ID:** 42042

**Keywords:** gut, intestine, human, nuclei, single cell, nuclei isolation from human intestinal biopsic tissue, nuclei from human intestinal biopsy sample, human intestinal biopsy sample, human intestinal biopsic tissue, nuclei isolation, cell genomic applications this protocol, cell genomic application, single cell application, atac seq, rna, genomic application,

## Disclaimer

The lysis buffer is formulated from the recipe in:

Drokhlyansky E, Smillie CS, Van Wittenberghe N, et al. The Human and Mouse Enteric Nervous System at Single-Cell Resolution [published online ahead of print, 2020 Aug 21]. *Cell*. 2020;S0092-8674(20)30994-6. doi:10.1016/j.cell.2020.08.003

## Abstract

This protocol provides an efficient method to isolate nuclei from human intestinal biopsy samples for single cell applications (RNA-seq or ATAC-seq).

## Guidelines

The human intestinal tissue were obtained with patient consent and approval by the Institutional Review Board at the University of Chicago (IRB Number: 15573A). All the samples are processed for research use only.



## Materials

### MATERIALS

☒ 5M Sodium Chloride, 1000ml **Promega Catalog #V4221**

☒ BSA **Merck MilliporeSigma (Sigma-Aldrich) Catalog ##A8806**

☒ RiboLock RNase Inhibitor (40 U/μL) **Thermo Fisher Catalog #EO0381**

☒ 0.5M EDTA **Fisher Scientific Catalog #2482-500**

☒ 10 x PBS no calcium no magnesium **Fisher Scientific Catalog #BP399500**

☒ UltraPure™ DNase/RNase-Free Distilled Water **ThermoFisher Catalog #10977023**

☒ Red blood cell lysis buffer 10x **Miltenyi Biotec Catalog #130-094-183**

☒ Tween 20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7949**

☒ 1M Tris-HCl pH 7.5 **Thermo Fisher Scientific Catalog #15567027**

☒ 1M CaCl<sub>2</sub> **Merck MilliporeSigma (Sigma-Aldrich) Catalog #21115**

☒ 1M MgCl<sub>2</sub> **Merck MilliporeSigma (Sigma-Aldrich) Catalog #63069**

#### **Lysis buffer** 10 ml (make fresh)

5 ml 2x ST buffer

300 μl 1% Tween-20

50 μl 2% BSA

10 μl RNase Inhibitor stock

4.64 ml UltraPure water

#### **2x ST buffer** 10 ml (Store at 4 Celsius up to 1 month)

292 mM NaCl

20 mM Tris-HCl pH 7.5

2 mM CaCl<sub>2</sub>

42 mM MgCl<sub>2</sub>

Bring up to volume with UltraPure water

#### **RBC lysis buffer** 10 ml

1 ml Red blood cells lysis buffer 10x

9 ml ultra pure water

#### **2% BSA** 10 ml (Store at 4 Celsius up to 1 month)

0.2 g BSA

10 ml UltraPure water

#### **1% Tween-20** 10 ml (Store up to 1 month)

1ml 10% Tween-20

9 ml UltraPure water



**Nuclei suspension buffer** 10 ml (make fresh)

10 ul RNase Inhibitor stock

50 ul 2% BSA

9.94 ml 1x PBS

**1x PBS** 500 ml (filter through 0.2 uM filter top)

50 ml 10x PBS

450 ml UltraPure water

## Troubleshooting



## Sample preparation

- 1 Rinse fresh samples in ice-cold PBS twice.

### Note

Biopsy tissue can be store up to 3 days in liquid nitrogen/at - 80 Celsius following the steps below:

Tissues are rinsed in ice-cold PBS twice

Flash freeze the tissue in 1.7 ml Eppendorf tube in liquid nitrogen

Store frozen tissue in liquid nitrogen (preferred) or at -80C


Start from step 2 if working with frozen tissues.

## Tissue lysis

- 2 Mince the tissue, with 200 ul lysis buffer added, in a 1.7 ml Eppendorf tube by Iris Scissors on ice x 1 mins.
- 3 Add 1-1.5 ml ice-cold lysis buffer to the tube and incubate on ice x 5 mins. Invert the tube 3 times in the middle of the incubation to mix.
- 4 Wet a 40 micron cell strainer with 1 ml lysis buffer.
- 5 Filter the lysis through the strainer. Wash the strainer by 3 ml lysis buffer and 4 ml nuclei suspension buffer (NSB). Keep the flow through as this is where your nuclei are.

## Nuclei collection

- 6 Spin down the flow-through, at 600 g x 5 mins, at 4 celsius in a 15 ml conical tube.

 600 x g, 4°C, 00:05:00

- 7 Suspend nuclei in 100 ul NSB using gentle pipetting.




### Note

If working with a single-cell platform, e.g. 10x Genomics, the nuclei should be suspended in PBS + 1% BSA + 0.2 U/ul RNase Inhibitor

- 7.1 If red blood cells are present in fresh tissue nuclei suspension, dilute the suspension with NSB to 1 ml. Add 2 ml RBC lysis buffer and incubate on ice for 5 minutes. Pellet nuclei by centrifugation 600 g x 5 mins, 4 Celsius. Suspend nuclei in 100 ul NSB with gentle pipetting.

 600 x g, 4°C, 00:05:00

- 8 Take 10 ul nuclei suspension and mix with 10 ul DAPI or Hoechst dye at 10 ug/ml and 10 ul WGA dye at 1 ug/ml. Count the nuclei.



DAPI/Hoechst staining-blue; WGA staining-orange Intact nuclei are co-stained by DAPI/Hoechst and WGA.

## Nuclei preparation

- 9 Dilute nuclei to the desired density using NSB.