

Jan 13, 2026

Nuclei isolation of fresh frozen human brain samples for for scRNA-seq

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvo1rddg4o/v1>

Koen Theunis^{1,2,3,4}, Suresh Poovathingal^{5,6}, Sarah Geurs⁷

¹Laboratory of Computational Biology, VIB Center for AI & Computational Biology, Leuven, Belgium;

²VIB-KU Leuven Center for Brain & Disease Research, Leuven, Belgium;

³Department of Human Genetics, KU Leuven, Leuven, Belgium;

⁴Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, United States;

⁵VIB KU Leuven Center for Brain & Disease Research;

⁶VIB Center for AI & Computational Biology (VIB.AI), Leuven, Belgium;

⁷VIB Tech Watch, VIB Technologies, Ghent, Belgium

Aertslab



Koen Theunis

Laboratory of Computational Biology, VIB Center for AI & Com...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvo1rddg4o/v1>



Protocol Citation: Koen Theunis, Suresh Poovathingal, Sarah Geurs 2026. Nuclei isolation of fresh frozen human brain samples for for scRNA-seq. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvo1rddg4o/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 09, 2025

Last Modified: January 13, 2026

Protocol Integer ID: 238288

Keywords: ASAPCRN, scRNA-seq, nuclei isolation, nuclei isolation for scrna, nuclei isolation of fresh frozen human brain sample, fresh frozen human brain sample, scrna, nuclei isolation, seq, nuclei

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000430

Aligning Science Across Parkinson's

Grant ID: ASAP-025179

Abstract

This protocol describes nuclei isolation for scRNA-seq.

Attachments



[SN-Isolation-forRNA_...](#)

18KB

Materials

	A	B
	HOW MANY SAMPLES TO PROCESS?	2
	Concentration of Sucrose?	2 M

2X salt-Tris solution

	A	B	C	D
	2X salt-Tris solution	Stock conc	Final conc	V (μl)
	NaCl	5000 mM	292 mM	1168.0
	Tris pH7.5	1000 mM	20 mM	400.0
	CaCl	1000 mM	2 mM	40.0
	MgCl ₂	1000 mM	42 mM	840.0
	Water			17552.0
	Total Volume			20000.0

TW20+ Homogenization buffer

	A	B	C	D
	TW20+ Homogenization buffer	Stock conc	Final conc	V (μl)
	2X salt-Tris solution	2 X	1 X	1100.0
	Tween-20	10%	0.03%	6.6
	BSA	7.5%	0.01%	2.9
	KCl	3000 mM	25 mM	18.3
	Sucrose	2000 mM	250 mM	275.0
	DTT*	1000 mM	1 mM	2.2
	cOmplete protease inhibitor (Roche EDTA free)*	50 X	0.5 X	22.0
	RNasin Plus (Promega)*	40 U/ul	0.5 U/ul	27.5
	Water			745.4

A	B	C	D
Total Volume			2200.0

TW20- Homogenization buffer

A	B	C	D
TW20- Homogenization buffer	Stock conc	Final conc	V (μl)
2X salt-Tris solution	2 X	1 X	572.0
BSA	7.5%	0.01%	1.5
KCl	3000 mM	25 mM	9.5
Sucrose	2000 mM	250 mM	143.0
DTT*	1000 mM	1 mM	1.1
cOmplete protease inhibitor (Roche EDTA free)*	50 X	0.2 X	4.6
RNasin Plus (Promega)*	40 U/ul	0.2 U/ul	5.7
Water			406.5
Total Volume			1144.0

GM (Gradient Medium)

A	B	C	D
GM (Gradient Medium)	Stock conc	Final conc	V(μl)
CaCl	1000 mM	1 mM	1.1
Optiprep*	60%	50%	953.3
MgCl ₂	1000 mM	5 mM	5.7
Tris pH 7.5	1000 mM	10 mM	11.4
Sucrose	2000 mM	75 mM	42.9
DTT*	1000 mM	1 mM	1.1
cOmplete protease inhibitor*	50 X	0.5 X	11.4
RNasin Plus (Promega)*	40 U/ul	0.5 U/ul	14.3
Water			102.6

	A	B	C	D
	Total Volume			1144

ODM (Optiprep Diluent Medium)

	A	B	C	D
	ODM (Optiprep Diluent Medium)	Stock conc	Final conc	V(μl)
	KCl	3000 mM	150 mM	48.1
	MgCl ₂	1000 mM	30 mM	28.9
	Tris pH 8	1000 mM	60 mM	57.8
	Sucrose	2000 mM	250 mM	120.3
	Water			707.6
	Total Volume			963

29% Cushion

	A	B	C	D
	29% Cushion	Stock conc	Final conc	V(μl)
	Optiprep	60%	29%	819
	ODM			875
	Total Volume			1694

Resuspension Buffer

	A	B	C	D
	Resuspension Buffer	Stock conc	Final conc	V(μl)
	PBS	1 X	1 X	192.5
	BSA	7.5%	0.75%	22
	RNasin Plus (Promega)*	40 U/uL	1 U/uL	5.5
				220

*Add these components just before use.



A. Preparation

5m

- 1 Prepare the 2X salt-Tris solution and aliquot to 1 ml tube and freeze at -20 °C .
- 2 Put the centrifuge for 4 °C .
- 3 Prepare the tween-20 lysis buffer fresh for each preparations.
- 4 Prepare cOmplete protease inhibitor (50X) by dissolving a tablet in 1 mL of water.
- 5 Filter sucrose solution and Protease inhibitor through 0.2 um/0.45 um filter.
- 6 Thaw all the reagents and prepare for the buffer reagents.
- 7 Place the homogenizer at -70 °C at the start of the experiment and place On ice 00:05:00 before homogenization.

5m

B. Homogenization

8m

- 8 Cut brain piece on dry ice (if large; > 20 mg - 50 mg tissue).
- 9 Transfer 750 µL of ice-cold TW20+ HBuffer and tissue pieces to a homogenizer and leave it in +HB for 00:03:00 before starting homogenization.
- 10 Homogenize tissue with **10** gentle strokes with **pestle A**, and **10** gentle strokes with **pestle B**.
- 11 Put homogenate through 70 µm cell strainer, wash homogenizer and filter with 250 µL **TW+HB** and leave the filtrate On ice . **Start the timer for 5 mins.**

3m



12 Transfer the filtrate to a 2 mL protein low bind eppendorf and centrifuge @  400 x g -  500 x g for  00:05:00 . Discard the supernatant.

5m



13 Resuspend the pellet gently in  200 μL of **TW20- HB**, and transfer to 2 mL low-bind tube.



14 Add additional **TW20- HB** to reach a final volume of  520 μL .



15 Add  520 μL of Gradient Medium (V_f =  1040 μL).



C. Isolation by centrifugation

20m

16 Layer  770 μL of 29% Cushion in the ultracentrifuge tube.

17 Layer the sample on top of cushion using a P1000, without disturbing the cushion.

18 Check tubes weight, adjust weight with HB/GM.

19 Centrifuge at least  3000 rcf with a swinging bucket rotor, at  4 $^{\circ}\text{C}$ for  00:20:00 with brake off.

20m



20 Remove supernatant, remove the lower supernatant with P200 or less, leaving around  50 μL -  100 μL .

D. Resuspension and counting

5m

21 Disrupt nuclei on bottom by pipetting with P200. And bring total volume to a fresh 1.5 mL LoBind tube.



22 Add  100 μL of 1X Resuspension buffer to ultracentrifuge tube to retain extra nuclei, and add them to the 1.5 mL tube.





23 Spin ⌚ 00:05:00 @ 🌀 350 x g - 🌀 450 x g followed by resuspension in Resuspension buffer. Optional filter and count.

5m



Protocol references

This protocol was based on earlier work from others but with in-house modifications:

Marine Denechaud, Sarah Geurs, Thomas Comptdaer, Séverine Bégard, Alejandro Garcia-Núñez, Louis-Adrien Pechereau, Thomas Bouillet, Yannick Vermeiren, Peter P. De Deyn, Romain Perbet, Vincent Deramecourt, Claude-Alain Maurage, Michiel Vanderhaegen, Sebastiaan Vanuytven, Bruno Lefebvre, Elke Bogaert, Nicole Déglon, Thierry Voet, Morvane Colin, Luc Buée, Bart Dermaut, Marie-Christine Galas.

Tau promotes oxidative stress-associated cycling neurons in S phase as a pro-survival mechanism: Possible implication for Alzheimer's disease.

Progress in Neurobiology, Volume 223, 2023, 102386, ISSN 0301-0082,

<https://doi.org/10.1016/j.pneurobio.2022.102386>.

L. Swiech, M. Heidenreich, A. Banerjee, N. Habib, Y. Li, J. Trombetta, M. Sur, F. Zhang

In vivo interrogation of gene function in the mammalian brain using CRISPR-Cas9

Nat. Biotechnol., 33 (2015), pp. 102-106, [10.1038/nbt.3055](#)