

Nov 20, 2020

Nuclei preparation from human lung using grinding for snRNA-seq

DOI

dx.doi.org/10.17504/protocols.io.bprdmm26

Center for Epigenomics UCSD¹

¹UCSD

LungMap2 Consortium

Tech. support email: lungmap2dcc@gmail.com

[Click here to message tech. support](#)



Sebastian Preissl

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bprdmm26

Protocol Citation: Center for Epigenomics UCSD 2020. Nuclei preparation from human lung using grinding for snRNA-seq. [protocols.io](https://dx.doi.org/10.17504/protocols.io.bprdmm26) <https://dx.doi.org/10.17504/protocols.io.bprdmm26>

Manuscript citation:

Wang, A. et al. Single cell multiomic profiling of human lung reveals cell type-specific and age-dynamic control of SARS-CoV2 host genes. Elife 9, doi:10.7554/eLife.62522 (2020).

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 16, 2020

Last Modified: November 20, 2020

Protocol Integer ID: 44549

Keywords: lung, nuclei preparation, single nucleus RNA-seq



Abstract

Protocol describing nuclei isolation from frozen tissue including lung using grinding for droplet-based single nucleus RNA-seq.

Attachments



Grinding_snRNA-

seq_0...

195KB

