

Jan 12, 2023

Version 2

Nucleofection (Amaxa) and electroporation (Biorad) of hPSCs V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwkkewvmk/v2>

Hanqin Li¹, Steven Poser², Dirk Hockemeyer¹, Frank Soldner¹, Oriol Busquets²

¹University of California, Berkeley; ²Albert Einstein College of Medicine



Hanqin Li

UC Berkeley

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwkkewvmk/v2>

External link: <https://doi.org/10.7554/eLife.79208>

Collection Citation: Hanqin Li, Steven Poser, Dirk Hockemeyer, Frank Soldner, Oriol Busquets 2023. Nucleofection (Amaxa) and electroporation (Biorad) of hPSCs. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwkkewvmk/v2> Version created by [Devin E Snyder](#)

Manuscript citation:

Hanqin Li, Oriol Busquets, Yogendra Verma, Khaja Mohieddin Syed, Nitzan Kutnowski, Gabriella R Pangilinan, Luke A Gilbert, Helen S Bateup, Donald C Rio, Dirk Hockemeyer, Frank Soldner (2022) Highly efficient generation of isogenic pluripotent stem cell models using prime editing eLife 11:e79208

<https://doi.org/10.7554/eLife.79208>

License: This is an open access collection distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this collection and it's working

Created: January 12, 2023

Last Modified: December 13, 2024

Collection Integer ID: 75230

Keywords: ASAPCRN, cultured hpscs for nucleofection, nucleofection nucleofection of hpscs electroporation, cultured hpscs for electroporation, plasmids for nucleofection, free hpscs for nucleofection, hpscs electroporation, plasmids for electroporation, preparing plasmid, mrna for nucleofection, cultured hpsc, preparing mrna, delivery of plasmid, transfected cell, enrichment of transfected cell, human pluripotent stem cell, rnp for nucleofection, using nucleofection, free hpsc, nucleofection nucleofection, hpsc, term hpsc, nucleofection, ribonucleoprotein, plasmid, mrna, described procedure, electroporation, standard procedure, biorad, standard procedure for the delivery

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000486



Abstract

This collection describes the standard procedure for the delivery of plasmids, mRNA or ribonucleoprotein (RNP) into human pluripotent stem cells (hPSCs) using nucleofection or electroporation.

Collection Overview

Preparing MEF-cultured hPSCs for nucleofection

Preparing Feeder-free hPSCs for nucleofection

Preparing plasmids for nucleofection

Preparing mRNA for nucleofection

In vitro assembling of RNP for nucleofection

Nucleofection of hPSCs

Electroporation of hPSCs

A. Preparation of MEFs-cultured hPSCs for electroporation

B. Preparing plasmids for electroporation

C. Electroporation

FACS-based enrichment of transfected cells

A. Preparation of samples for FACS sorting

B. After FACS

General notes

1. Throughout these protocols, the term hPSC is used to collectively refer to both hiPSCs and hESC. All described procedures have been tested and work equally well for hiPSCs and hESCs.

Materials

	Item	Vendor	Catalog #
	DMEM/F12	Thermo Fisher	11320082
	DPBS w/o calcium and magnesium	Corning	MT21031CV
	Fetal Bovine Serum (FBS)	Corning	35-011-CV
	Knockout Serum Replacement	Thermo Fisher	10828-028
	L-Glutamine	Sigma	G8540
	Penicillin & Streptomycin (100X)	Thermo Fisher	15140163
	MEM Non-Essential Amino Acids (100X)	Thermo Fisher	11140050
	Heat Stable Recombinant Human FGF2	Thermo Fisher	PHG0360
	Collagenase type IV	Thermo Fisher	17104019
	2-Mercaptoethanol	Sigma	M3148
	mTeSR-plus	STEMCELL Technologies	100-0276
	StemFlex	Thermo Fisher	A3349401
	Vitronectin (VTN-N) Recombinant Human Protein, Truncated	Thermo Fisher	A14700
	Accutase	Thermo Fisher	SCR005
	Dispase	STEMCELL Technologies	NC9995391
	Y-27632	Chemdea	CD0141
	Cas9, purified protein, 40uM	Macrolab, QB3 UC Berkeley	
	Synthetic pegRNAs	IDT or Synthego	
	Synthetic sgRNAs	Synthego	
	P3 primary Cell 4D X kit S	Lonza	V4XP-3032



	Item	Vendor	Catalog #
	Countess™ Cell Counting Chamber Slides	Thermo Fisher	C10228
	pCMV-PE2	Addgene	132775
	4D-Nucleofector TM Core + X Unit	Lonza	AAF-1002B, AAF-1002X
	5 ml polystyrene round-bottom tube with cell-strainer cap	Corning	352235
	Cell-strainer (70 µm)	Fisher	07201431
	Gene Pulser Xcell Eukaryotic System	Bio-Rad	1652661
	Gene Pulser Electroporation Cuvettes, 0.4 cm gap	Bio-Rad	1652081
	Exact N Amp Blood PCR Kit	Sigma	XNAB2-1TK



Files

 SEARCH

Protocol

NAME

Preparing MEF-cultured hPSCs for nucleofection

VERSION 2

CREATED BY



Hanqin Li
UC Berkeley

OPEN →

Protocol

NAME

Preparing feeder-free hPSCs for nucleofection

VERSION 1

CREATED BY



Hanqin Li
UC Berkeley

OPEN →

Protocol

NAME

Preparing plasmids for nucleofection of hPSCs

VERSION 1

CREATED BY



Hanqin Li
UC Berkeley

OPEN →

Protocol

NAME

Preparing mRNA for nucleofection of hPSCs

VERSION 1



CREATED BY



Hanqin Li
UC Berkeley

OPEN →

Protocol

NAME

In vitro assembling of RNP for nucleofection of hPSCs

VERSION 1

CREATED BY



Hanqin Li
UC Berkeley

OPEN →

Protocol

NAME

Nucleofection of hPSCs

VERSION 1

CREATED BY



Hanqin Li
UC Berkeley

OPEN →

Protocol

NAME

Electroporation of hPSCs

VERSION 1

CREATED BY



Hanqin Li
UC Berkeley

OPEN →

Protocol

NAME

FACS-based enrichment of transfected hPSCs

VERSION 1

CREATED BY



Hanqin Li
UC Berkeley

OPEN →

