



Jul 31, 2024

OHSU SenNet Senescence-Associated Beta-Galactosidase (SA b-Gal) Staining of Carboxymethyl Cellulose (CMC) Embedded Skin Tissue Sections on Polyethylene Naphthalate (PEN) Coated Slides

DOI

dx.doi.org/10.17504/protocols.io.81wgbzq7ygpk/v1

Julia R Unsworth¹, Xianmin Luo², Sheila A. Stewart^{2,3,4}, Megan K. Ruhland^{1,5}

¹Cell, Developmental and Cancer Biology Department, Oregon Health & Science University;

²Department of Cell Biology and Physiology, Washington University School of Medicine, St. Louis, Missouri;

³Siteman Cancer Center, Washington University School of Medicine, St. Louis, Missouri Department of Medicine, Washington University School of Medicine, St. Louis, Missouri;

⁴ICCE Institute, Washington University School of Medicine, St. Louis, Missouri;

⁵Dermatology Department, School of Medicine, Oregon Health & Science University

Cellular Senescence Net...



Julia R Unsworth

Oregon Health & Science University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS





DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbzq7ygpk/v1>

Protocol Citation: Julia R Unsworth, Xianmin Luo, Sheila A. Stewart, Megan K. Ruhland 2024. OHSU SenNet Senescence-Associated Beta-Galactosidase (SA b-Gal) Staining of Carboxymethyl Cellulose (CMC) Embedded Skin Tissue Sections on Polyethylene Naphthalate (PEN) Coated Slides. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.81wgbzq7ygpk/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 30, 2024

Last Modified: July 31, 2024

Protocol Integer ID: 104333

Keywords: galactosidase staining, gal staining of mouse skin tissue section, staining of carboxymethyl cellulose, mouse skin tissue section, senescent cell, skin tissue sections on polyethylene naphthalate, coated slides senescence, skin tissue section, gal staining, galactosidase, ohsu sennet senescence, staining technique, carboxymethyl cellulose, staining, use of laser capture microdissection, downstream mass spectrometry analysis

Funders Acknowledgements:

Spatially-resolved proteome mapping of senescent cells and their tissue microenvironment at single-cell resolution

Grant ID: UG3CA275697

Disclaimer

The **protocols.io** team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

Abstract

Senescence-associated beta-galactosidase staining is a common staining technique to detect senescent cells. This protocol describes SA b-Gal staining of mouse skin tissue sections from CMC embedded blocks onto PEN slides which enable the combined use of laser capture microdissection (LCM) and UV laser cutting for downstream mass spectrometry analysis.

Troubleshooting

TISSUE SECTIONING AND EMBEDDING PROCEDURE

- 1 Collecting the Skin
 - 1.1 From the wet lab take petri dishes for each mouse that you are collecting tissue from and place in a bucket of ice
 - 1.2 Once in mouse room, follow protocol for euthanizing the mice
 - 1.3 Once the mice are euthanized, use an electric razor to shave the backs of each mouse
 - 1.4 Apply a generous amount of Nair onto the back of each mouse and let sit for ~ 1 min
 - 1.5 Using a paper towel, carefully wipe the Nair off the mouse's back and spray with ethanol to clean area
 - 1.6 Using scissors and a pair of tweezers, cut a rectangle shape of the dorsal skin (from upper shoulders down to lower hips) and place epithelial side up onto the appropriately labeled petri dish
 - 1.7 Bring the tissue samples back to the wet lab and prepare an ice bucket of dry ice
 - 1.8 Using a clean razor blade, cut off the uneven edges of each tissue segment and then cut into 6-8 smaller rectangle pieces that are all roughly the same size
- 2 Embedding the Skin in CMC
 - 2.1 Label two Tissue-Tek Cryomolds per each mouse and pipette Carboxymethyl Cellulose (CMC) media into them until they are halfway full
 - 2.2 Using tweezers, carefully pick up the cut sections of tissue and place them vertically in the CMC media, making sure they are all oriented the same direction. Place 3-4 tissue pieces in each cryomold
 - 2.3 Add more CMC media once tissue pieces are in the cryomold.



Note: Try to have the tissue pieces not sink towards the bottom, and have them all relatively on the same plane

- 2.4 Once tissue sections are vertical, place the cryomolds onto dry ice to harden the CMC media. Add enough CMC media to where you do not see the tissue sections
- 2.5 Allow the cryomolds to harden on dry ice for at least 15 minutes before placing in -80°C freezer for long term storage

3 Tissue Sectioning on the Cryostat

- 3.1 Take CMC embedded tissue sections out of -80°C freezer, remove block from the cryomold and mount onto pedestal in the Cryostat chamber using a small amount of OCT media to adhere it. Allow to set for ~15-30min
Note: Cryostat should be set at -20°C
- 3.2 Cut tissue sections at 10 µm thickness in the cryostat at -20°C and place onto Polyethylene Naphthalate (PEN) slides
Note: Make sure that the tissue is placed onto the membrane bound side of the slid
- 3.3 Move slides into -80°C freezer until ready for staining procedure

STAINING PROCEDURE

4 Fixation of Slides

- 4.1 Take tissue slides from -80°C freezer and immediately place in fixation solution for 10 minutes at 4°C on shaker

Fixation Solution

0.2% glutaraldehyde
2% glycerol
in PBS

- 4.2 Wash slides in 2% glycerol in PBS for 2 minutes on shaker at room temperature

5 X-gal Staining

- 5.1 Make staining solution

	Actual Con.	Stocks	Amount needed for 5ml
	1mg/ml X-gal	40mg/ml	125µl
	150mM NaCl	5M	150µl
	2mM MgCl ₂	1M	10µl
	5mM K ₃ Fe(CN) ₆	100mM	250µl (protect from light)
	5mM K ₄ Fe(CN) ₆	100mM	250µl (protect from light)
	40mM NaPi pH 6.0	0.1M	2 ml
	H ₂ O		2.217 ml

* K₃Fe(CN)₆ and K₄Fe(CN)₆ should be protected from light

* Keep X-gal stock at -20°C

To make up 0.1M NaPi stock (remake every 2 months):

	Actual Con.	Stocks	Amount needed for 10ml
	Na ₂ HPO ₄	1M	125µl
	NaH ₂ PO ₄	1M	880µl
	H ₂ O		8.995 ml

- 5.2 Check the pH of the NaPi stock using a pH meter, it should be at 6.0
- 5.3 Check the pH of the final staining solution, for best results the pH should be at 6.0
- 5.4 Filter the staining solution using a .45µm filter to prevent crystals from forming. **Protect from light**
- 5.5 Use a Kimwipe to absorb excess PBS. Do not allow the tissue section to dry out



- 5.6 Using a PAP pen, draw a square around the tissue and place ~200ml of staining solution onto tissue and incubate for 6-24 hours at 37°C in a humidified chamber/slide holder
Note: The timing depends on the negative (non-senescent) control tissue, you will need to adjust the timing based on what is appropriate for your sample
- 5.7 To stop staining reaction, fix again for 10 minutes in fixation solution at room temperature on shaker
- 5.8 Wash slides in 2% glycerol in PBS for 5 minutes on shaker at room temperature
- 6 Nuclear counter stain: Nuclear fast red is a counterstain that targets nucleic acids and helps identify the nucleus of cells
- 6.1 Place slides in Nuclear Fast Red for 10 minutes at room temperature
- 6.2 Rinse slides with tap water 3 times
- 6.3 Place slides in tap water for 3 minutes on shaker
- 6.4 Let slides air-dry at room temperature
Note: For laser ablation, cover slip and mounting media were not used.