

Oct 24, 2025

On-bead digestion under denaturing conditions

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg32kypv25/v1>



Cristina CARDENAL PERALTA¹

¹University of Edinburgh



Cristina CARDENAL PERALTA

University of Edinburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg32kypv25/v1>

Protocol Citation: Cristina CARDENAL PERALTA 2025. On-bead digestion under denaturing conditions. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg32kypv25/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 24, 2024

Last Modified: October 24, 2025

Protocol Integer ID: 102292

Keywords: beads, denaturing conditions, urea, trypsin, bead digestion, mass spectrometry, mass spectrometry analysis, efficient protein digestion, accuracy of mass spectrometry analysis, peptide recovery, enhancing peptide recovery, sample preparation process

Abstract

On-bead digestion is used for mass spectrometry to streamline the sample preparation process, reducing sample loss and contamination risks. It allows for efficient protein digestion directly on the solid support, enhancing peptide recovery and improving the sensitivity and accuracy of mass spectrometry analysis by minimizing handling steps and ensuring more consistent results.

Guidelines

Perform all steps at  Room temperature .

Materials

Buffers:

Denaturation buffer: 8M urea in 50 mM ABC. (For membrane proteins use 7M urea/2M thiourea).

Digestion buffer: 50 mM ABC (ammonium bicarbonate) in water, pH 8.

Reagents:

DTT: 1 M in 50 mM ABC.

IAA: 1M in 50 mM ABC

Enzymes:

LysC: 1 µg/ µL in 0.1% TFA

Trypsin: 1 µg/ µL in 0.1% TFA

Safety warnings

 In this procedure all steps are done at  Room temperature to reduce unwanted denaturization of amino acid side-chains by denaturing agents. **Never heat your sample.**



Sample reception

5m

- 1 When the samples arrive isolate beads from the purification buffer using a magnetic tube rack. Allow 5 minutes of incubation and remove purification buffer.
- 2 Resuspend the beads in three times the original bead volume using denaturation buffer

5m

Reduction/Alkylation

1h

- 3 Add 1 μ L of [M] 1 Molarity (M) DTT per 10 μ L of sample volume and incubate for 30 min at Room temperature . Shake at 1200 rpm.
- 4 Add 1 μ L of [M] 1 Molarity (M) IAA per 10 μ L of sample volume and incubate for 30 min at Room temperature . Shake at 1200 rpm in the dark.

30m



30m



Digestion

4h

- 5 Add 1 μ g of LysC solution per 50 μ g of protein and incubate for at least 3-4 hours at Room temperature .
- 6 Dilute sample with 5x digestion buffer.
- 7 Add 4 μ g of trypsin per 50 μ g of protein and incubate O/N at 37 °C .
The trypsin and the LysC digestion steps are interchangeable, however, trypsin does not work in high urea concentrations, and dilution of the sample is needed before adding trypsin.

4h



Sample acidification

- 8 Separate beads from supernatant by using a magnetic tube rack (allow 5 min for all the beads to stick to the side of the tube). **Collect supernatant into a new Eppendorf tube.**
- 9 Stop the digestion by acidifying the sample to pH < 2.5 with 10 % TFA, to a final % of 1%.



- 10 Proceed to load the sample onto equilibrated C18 StageTips (samples can stay in the fridge until the tips are equilibrated).