

Jul 02, 2025

Open Field Behavior in Mice

 Forked from [Open Field Behavior in Mice](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl88rb9l2w/v1>

Yue Sun¹, Richard Roth¹, Fuu-Jiun Hwang¹, Jun B. Ding¹, Shenyu Zhai²

¹Stanford University; ²Northwestern University



Shenyu Zhai

Northwestern University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl88rb9l2w/v1>

Protocol Citation: Yue Sun, Richard Roth, Fuu-Jiun Hwang, Jun B. Ding, Shenyu Zhai 2025. Open Field Behavior in Mice. protocols.io <https://dx.doi.org/10.17504/protocols.io.5jyl88rb9l2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: July 02, 2025

Last Modified: July 02, 2025

Protocol Integer ID: 221428

Keywords: ASAPCRN, open field, mouse, open field, open field behavior in mice, open field test, like responses in rodent, spontaneous locomotor activity, open field behavior, mice, rodent, anxiety, exploratory behavior, field, field test, test

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-020551

Abstract

The open field test is widely used to assess spontaneous locomotor activity, exploratory behavior, and anxiety-like responses in rodents. In this fork, the test duration has been modified and the software used described.

Attachments



Ding_Lab_Open field ...

15KB

Materials

Open field arena: 50 × 50 × 50 cm chambers

Video recording: A camera mounted above the chambers, adjusted to ensure full coverage of the arena.

ANY-maze software (<https://www.any-maze.com/>)

70% ethanol spray and paper towel

Safety warnings

! Wear appropriate PPE as required by your institution.

Ethics statement

Prior ethics approval (e.g. IACUC) should be obtained before performing these experiments. Approval was obtained by the Stanford University IACUC before any procedures were performed.

Before start

Customize the settings on the ANY-maze software, including the type of animals used, parameters to be tracked, and calibration of the camera feed.

Procedure

- 1 Habituate the mice for 45-60 min in the room where the testing will be performed.
- 2 Open the ANY-maze software and start a new experiment.
- 3 Place the mouse gently in the center of the open field arena.
- 4 Start recording.
- 5 The duration of each test session is 5 minutes, but may be changed depending on the experimental design.
- 6 After allotted time, stop the recording and return the mouse to its home cage.
- 7 **Optional observation:** Record the number of fecal boli and urine spots as additional indicators of anxiety.
- 8 Thoroughly clean the chambers with 70% ethanol between trials and after testing to eliminate odor cues from previous animals.