



Dec 05, 2025

Open field test

 [PLOS One](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvmqbqo5g3p/v1>

Mario Alberto Bautista-Carro¹

¹Cinvestav



Mario Alberto Bautista-Carro

Cinvestav

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvmqbqo5g3p/v1>

External link: <https://doi.org/10.1371/journal.pone.0339028>

Protocol Citation: Mario Alberto Bautista-Carro 2025. Open field test. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxmvmqbqo5g3p/v1>

Manuscript citation:

Bautista-Carro MA, Sánchez-Teoyotl P, Juárez-Serrano D, Iannitti T, Díaz A, Flores G, Morales-Medina JC (2026) Olfactory bulbectomy induces neurobiological alterations in the prefrontal cortex and hyperlocomotion in male rats. PLOS One 21(1). doi: [10.1371/journal.pone.0339028](https://doi.org/10.1371/journal.pone.0339028)



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 08, 2025

Last Modified: December 05, 2025

Protocol Integer ID: 229260

Keywords: Locomotor activity, Exploration, Novel environment, Anxiety-like behavior, exploration in the open field test, open field test, open field test the open field test, test, field, open enclosure, laboratory, rodent, scientific research method, rat

Abstract

The open field test is a scientific research method that measures the general locomotion and willingness to explore of a rodent by placing it in a large, open enclosure. This protocol uses a black square arena without a top cover, in a field divided into 25 squares of equal size, with 300 LUX lighting. During the rat's exploration in the open field test, the researcher leaves the laboratory. To control for possible effects due to circadian cycles, all tests are performed at the same time of day. All tests are recorded on video with a central view of the field at a height of 2 m above the ground.

Guidelines

All tests should be conducted at approximately the same time of day. During the rat's exploration of the open field, the researcher should leave the laboratory.

Materials

Lux meter

Box with square base (85 × 85 cm) in black

Lamp

Canon video camera (model Vixia HF R70) with a 2.8-89.6 mm lens

Camera tripod

70% alcohol

Before start

Verify that there is uniform 300 LUX lighting within the open field and clean it with 70% alcohol on the base and walls. Transfer the rats to the behavior laboratory and wait 1 hour for habituation before beginning the open field test.

- 1 After 1 h of habituation in the laboratory, place the rat in the center of the field and allow it to explore for 5 min.
- 2 Remove the rat from the field and place it in a new box separate from its group.
- 3 After each trial, clean the base and walls of the open field test with 70% alcohol to eliminate olfactory cues.
- 4 Record locomotor activity as well as the frequency of rearing and grooming behaviors in the rats, using a blind observer.

Protocol references

Bautista-Carro, M. A., Flores, G., Iannitti, T., Galindo-Paredes, G., 26 Morales-Medina, J. C. (2021). Curcuma longa administration significantly reduces acute and persistent inflammatory pain measures in male and female rats. *Arch. Vet. Sci. Med.*, 4(02). doi:10.26502/avsm.022

Morales-Medina, J. C., Dumont, Y., Benoit, C. E., Bastianetto, S., Flores, G., Fournier, A., 26 Quirion, R. (2012). Role of neuropeptide Y Y1 and Y2 receptors on behavioral despair in a rat model of depression with co-morbid anxiety. *Neuropharmacology*, 62(1), 200208. <https://doi.org/10.1016/j.neuropharm.2011.06.030>