

Jul 29, 2024

Version 1

Opentrons Pipeline: DNA Extraction with the Mag-Bind Blood & Tissue DNA HDQ 96 Kit Tissue Protocol V.1

DOI

dx.doi.org/10.17504/protocols.io.bp2l62rqrgqe/v1

Thalía Silvestre¹

¹UNMSM



Thalía Silvestre

UNMSM

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l62rqrgqe/v1>

Protocol Citation: Thalía Silvestre 2024. Opentrons Pipeline: DNA Extraction with the Mag-Bind Blood & Tissue DNA HDQ 96 Kit Tissue Protocol. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l62rqrgqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol



Created: July 27, 2024

Last Modified: July 29, 2024

Protocol Integer ID: 104196

Keywords: Opentrons, OT-2, Mag-Bind Blood & Tissue DNA HDQ 96 Kit Tissue Protocol, tissue dna hdq 96 kit tissue protocol, tissue dna hdq, dna extraction, kit tissue protocol, pipeline, automated pipeline, automated pipetting robot, pipetting robot, well plate of dna, bind blood, dna, extraction, blood, dna extraction with the mag

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment.

Abstract

This protocol is an automated pipeline to extract a full 96-well plate of DNA from tissue lysates.

This protocol was developed and optimized for the following:

Platform: Opentrons OT-2 automated pipetting robot

Kit: Mag-Bind® Blood & Tissue DNA HDQ Prefilled 96 Kit

Recommended number of samples: 96

Two ellutions: Firts ellution(60μL) , Second ellution (30μL)

Troubleshooting



Samples Lysis

- 1 Prepare 96 samples lysis following the Mag-Bind[®] Blood & Tissue DNA HDQ 96 Kit Tissue Protocol:
 - 1.1 Mince up to 10 mg tissue and transfer to a 1.5 mL tubes (one for each sample). (Cut each 10 mg sample with sterilized cutters as small as possible)
 - 1.2 Add 250 μ L TL Buffer to each sample.
 - 1.3 Add 20 μ L Proteinase K Solution to each sample. Vortex to mix.
 - 1.4 Incubate at 55°C overnight in a shaking water bath.
 - 1.5 Centrifuge at maximum speed ($\geq 4,000 \times g$) for 5 minutes to pellet undigested tissue debris and transfer 200 μ L of the supernatant to the VWR 96-Well Deep Well Plate (1000 μ L) without disturbing the undigested pellet

List of materials to start to use OT-2

- 2 It is necessary to have this equipment to use this pipeline.
Opentrons Equipment List:

	A	B
	OT-2	
	OT-2 8-Channel Pipette P300	
	OT-2 Magnetic Module GEN2	Slot 7



3 Tips & Labware:

A	B
Four Opentrons 200µL Filter Tips	Slot: 5,6,8,9
Two Nest 1-Well Reservoirs (195 mL)	Slot: 10 and 11
VWR 96-Well Deep Well Plate (1000 µL)	Slot 7 on the Magnetic module
Two Nest Well plate (100 µL)	Slot 2 and 3
Two Nest 12 Well Reservoirs (15 mL)	Slot 1 and 4

Star OT-2 run

- 4 The pipeline start on the step number 7 of the Mag-Bind® Blood & Tissue DNA HDQ 96 Kit Tissue Protocol

Load the pipeline on the Opentrons app:



Extracción.py 11KB

The VWR 96-Well Deep Well Plate (1000 µL) labware can be download from:

dx.doi.org/10.17504/protocols.io.dm6gp39njvzp/v1

4.1 Order of the reagents on the Two Nest 12 Well Reservoirs (15 ml)

Slot 1:



	A	B	C	D	E	F	G	H	I	J	K	L
	Well 1	Well 2	Well 3	Well 4	Well 5	Well 6	Well 7	Well 8	Well 9	Well 10	Well 11	Well 12
	11.5 mL of AL buffer	11.5 mL of AL buffer	Mag - Bind mixed with HBQ buffer (20:340) in a total of 12 mL	Mag - Bind mixed with HBQ buffer (20:340) in a total of 12 mL	Mag - Bind mixed with HBQ buffer (20:340) in a total of 12 mL	15 ml of VHB buffer	15 ml of VHB buffer	15 ml of VHB buffer	15 ml of VHB buffer			

Slot 4:

	A	B	C	D	E	F	G	H	I	J	K	L
	Well 1	Well 2	Well 3	Well 4	Well 5	Well 6	Well 7	Well 8	Well 9	Well 10	Well 11	Well 12
	15 ml of VHB buffer	15 ml of VHB buffer	15 ml of VHB buffer	15 ml of VHB buffer	15 ml of SPM buffer	15 ml of SPM buffer	15 ml of SPM buffer	15 ml of SPM buffer				Ellution Buffer