

May 30, 2025

## Optimized heat shock transformation protocol for Escherichia coli TOP10, JM109, and BL21 strains

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm6wbv3p/v1>

Reclone Latin America Hub<sup>1</sup>

<sup>1</sup>Reclone (The Reagent Collaboration Network)

Reclone.org (The Reagent Collaboration Network)

Tech. support email: [protocols@reclone.org](mailto:protocols@reclone.org)

[Click here to message tech. support](#)



Karla Vanessa Gaona

Instituto de Bioquímica y Biotecnología, Facultad de Ciencia...

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.yxmvm6wbv3p/v1>

**Protocol Citation:** Reclone Latin America Hub 2025. Optimized heat shock transformation protocol for Escherichia coli TOP10, JM109, and BL21 strains. [protocols.io https://dx.doi.org/10.17504/protocols.io.yxmvm6wbv3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvm6wbv3p/v1)



**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 30, 2025

**Last Modified:** May 30, 2025

**Protocol Integer ID:** 219160

**Keywords:** optimized heat shock transformation protocol, escherichia coli top10, competent escherichia coli, heat shock approach, jm109, heat, escherichia coli

**Funders Acknowledgements:**

Chan Zuckerberg Initiative

## Abstract

This protocol describes the standard method for transforming chemically competent *Escherichia coli* using a heat shock approach.

## Materials

- Chemically competent *Escherichia coli* cells (e.g., TOP10, JM109, BL21)
- Plasmid DNA (10 pg – 100 ng)
- SOC medium
- Ice
- Sterile 1.5 mL microcentrifuge tubes
- Water bath or heat block at 42 °C
- Incubator or water bath at 37 °C
- Sterile LB agar plates with appropriate antibiotic
- Vortex mixer



## Protocol

- 1 Thaw 50  $\mu$ L of chemically competent *Escherichia coli* cells (tested with TOP10, JM109, and BL21 strains) on ice for 5–10 minutes.
- 2 Add 1  $\mu$ L of plasmid DNA (between 10 pg and 100 ng) and gently mix by tapping or pipetting.
- 3 Incubate the tube on ice for 30 minutes.
- 4 Heat shock the cells by placing the tube at 42 °C for 45 seconds in a water bath or incubator.
- 5 Immediately return the tube to ice and incubate for 2 minutes.
- 6 Add 400  $\mu$ L of SOC medium to the tube.
- 7 Incubate the cells at 37 °C for 1 hour with shaking or static incubation in a water bath or incubator (to allow expression of antibiotic resistance genes).
- 8 Plate 100  $\mu$ L on an agar plate containing the corresponding antibiotic for selection, depending on the plasmid used.

## Results

- 9 The image below shows bacterial colonies grown on an LB agar plate supplemented with ampicillin, resulting from the transformation of *E. coli* BL21 cells.



**Figure 1:** *Bacterial colonies obtained after transformation of Escherichia coli BL21 cells.*