



Feb 27, 2025

Optimized Mung Bean Nuclease Nanoseq Library Preparation Protocol for Whole Genome Sequencing

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj97wwlk5/v1>

Hsu Chao¹, Kavya Kottapalli¹, Qi Jiang¹, Muchun Niu², Chenghang Chuck Zong², Harsha Doddapaneni¹

¹Human Genome Sequencing Center, Baylor College of Medicine;

²Department of Molecular and Human Genetics, Baylor College of Medicine



hchao Chao

BCM

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj97wwlk5/v1>

Protocol Citation: Hsu Chao, Kavya Kottapalli, Qi Jiang, Muchun Niu, Chenghang Chuck Zong, Harsha Doddapaneni 2025. Optimized Mung Bean Nuclease Nanoseq Library Preparation Protocol for Whole Genome Sequencing. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvj97wwlk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 07, 2025

Last Modified: February 27, 2025

Protocol Integer ID: 123163

Keywords: Duplex Sequencing, NanoSeq, Mung bean nuclease (MBN) , optimized mung bean nuclease nanoseq library preparation protocol, whole genome sequencing duplex sequencing, hpych4v enzyme for genomic dna fragmentation, using dideoxy nucleotide, mechanical fragmentation of genomic dna, genomic dna fragmentation, stranded dna, genomic dna, more genome coverage, mung bean nuclease, sequencing accuracy, sequencing method, dideoxy nucleotide, sequencing protocol, sequencing, genome, parent dna strand, dna, hpych4v enzyme, combining botseq, nanoseq method, identifying true mutation, nanoseq, generation sequencing

Funders Acknowledgements:

SMAHT Network (NIH)

Grant ID: 5UM1DA058229; UG3NS132132

Abstract

Duplex Sequencing (DS) is an ultra-accurate next-generation sequencing method that independently sequences both parent DNA strands, identifying true mutations shared between them and eliminating errors from artifacts or misincorporation. There are multiple duplex sequencing protocols, Nanoseq being one such method (PMID: 33911282). NanoSeq enhances duplex sequencing accuracy by combining BotSeqS (PMID: 27528664) with a restrictive end-repair method using dideoxy nucleotides, overcoming limitations of shallow sequencing and targeted panels for high-sensitivity applications. Here we describe step-by-step, the Nanoseq method using the mung bean nuclease, an enzyme that preferentially cleaves single-stranded DNA while preserving double-stranded DNA. Combining use of this enzyme following mechanical fragmentation of genomic DNA ensures nearly 2x more genome coverage compared to using the HpyCH4V enzyme for genomic DNA fragmentation. This method has been optimized and tested across >15 human tissues, including heart, muscle, colon, blood, cord blood, liver, and brain etc. to generate reproducible results.

Guidelines

Additional Notes

The accuracy of the qPCR-1 is critical quantification for optimizing the complexity of the sequencing library given the desired sequencing coverage, maximizing duplex coverage.

High adapter dimmer in step 7 library PCR could results in inefficient amplification.



Materials

Reagents

- TE Buffer **Thermo Fisher Scientific** Catalog # 12090015
- Mung Bean Nuclease **New England Biolabs** Catalog # M0250S
- Mung bean nuclease buffer **New England Biolabs** Catalog # B0250S
- UltraPure™ SDS Solution, 10% **Thermo Fisher Scientific** Catalog # 15553027
- T4 Polynucleotide Kinase **New England Biolabs** Catalog # M0201S
- NEBuffer™ 4 **New England Biolabs** Catalog # B7004S
- Klenow Fragment (3'→5' exo-) **New England Biolabs** Catalog # M0212L
- Dideoxynucleoside Triphosphate Set **Millipore Sigma** Catalog # 3732738001
- dATP Solution (100mM) **New England Biolabs** Catalog # N0440S
- T4 DNA Ligase **New England Biolabs** Catalog # M0202L
- Adenosine 5'-Triphosphate (ATP, 10mM) **New England Biolabs** Catalog # P0756S
- CS adapter (15 µM) **Integrated DNA Technologies (IDT)** Catalog # 1080799
- xGen™ UDI 10nt Primer Plates 1-4 **Integrated DNA Technologies (IDT)** Catalog # 10008052
- NEBNext® Ultra™ II Q5® Master Mix **New England Biolabs** Catalog # M0544L
- AMPure XP **Beckman Coulter** Catalog # A63882
- DEPC-Treated Water **Thermo Fisher Scientific** Catalog # AM9906
- 10 mM Tris-HCl Elution buffer **Invitrogen** Catalog# A33566
- DNA 7500 kit **Agilent Technologies** Catalog# 5067-1506
- High Sensitivity DNA Kit **Agilent Technologies** Catalog# 5067-4626
- KAPA Library Quant Complete kit (ABI Prism®) **Roche** Catalog # 07960204001
- Nuclease free water, **Thermo Fisher Scientific**, Catalog# AM9937
- 100 µM NanoqPCR primer 1: 5'-ACACTCTTTCCCTACACGAC-3' **IDT HPLC grade**
- 100 µM NanoqPCR primer 2: 5'-GTGACTGGAGTTCAGACGTG-3' **IDT HPLC grade**
- Ethanol absolute (200 Proof) **Avantor** Catalog # 89125-170

Equipments

- Covaris E220 Focused-ultrasonicator system **Covaris, Catalog # 500239**
- QuantStudio™ 6 Flex Real-Time PCR System,
384-well **Thermo Fisher Catalog # 4485691**
- PTC-200 Peltier Thermo Cycler **MJ Research Catalog # 8252-30-0001**

⊗ 1 x TE Buffer **Thermo Fisher Scientific** Catalog #12090015

⊗ Mung Bean Nuclease - 1,500 units **New England Biolabs** Catalog #M0250S

⊗ UltraPure® SDS Solution, 10% **Thermo Fisher** Catalog #15553027

⊗ T4 Polynucleotide Kinase - 500 units **New England Biolabs** Catalog #M0201S

⊗ NEBuffer 4 - 5.0 ml **New England Biolabs** Catalog #B7004S

⊗ Klenow Fragment (3'-5' exo-) - 1,000 units **New England Biolabs** Catalog #M0212L

⊗ Dideoxynucleoside Triphosphate Set **Merck MilliporeSigma (Sigma-Aldrich) Catalog #03732738001**

⊗ dATP Solution - 25 umol **New England Biolabs Catalog #N0440S**

⊗ T4 DNA Ligase - 100,000 units **New England Biolabs Catalog #M0202L**

⊗ Adenosine-5 Triphosphate (ATP) - 1 ml **New England Biolabs Catalog #P0756S**

⊗ xGen™ NGS Adapters & Indexing Primers for Illumina® sequencing **Integrated DNA Technologies, Inc. (IDT) Catalog #1080799**

⊗ xGen™ NGS Adapters & Indexing Primers for Illumina® sequencing **Integrated DNA Technologies, Inc. (IDT) Catalog #10008052**

⊗ NEBNext Ultra II Q5 Master Mix - 250 rxns **New England Biolabs Catalog #M0544L**

⊗ AMPure XP **Bechman Coulter Catalog #A63882**

⊗ DEPC-Treated Water **Thermo Fisher Catalog #AM9906**

⊗ Elution Buffer for Dynabeads®; mRNA Purification Kits **Thermo Fisher Catalog #A33566**

⊗ Bioanalyzer DNA Analysis DNA 7500 Kit **Agilent Technologies Catalog #5067-1506**

⊗ Agilent High Sensitivity DNA Kit **Agilent Technologies Catalog #5067-4626**

⊗ KAPA Library Quantification Kits **Roche Catalog #07960204001**

⊗ Nuclease-Free Water (not DEPC-Treated) **Thermo Fisher Scientific Catalog #AM9937**

⊗ Ethanol absolute, KOPTEC, meets analytical specification of BP, Ph. Eur., USP (200 Proof) **Avantor Sciences Catalog #89125-176**

Equipment

E220 Focused-ultrasonicator

NAME

Covaris

BRAND

500239

SKU

<https://www.covaris.com/e220-focused-ultrasonicator-500239>^{LINK}

Equipment

QuantStudio™ 6 Flex Real-Time PCR System, 384-well, laptop	NAME
Applied Biosystems™	BRAND
4485691	SKU
https://www.thermofisher.com/order/catalog/product/4485691	LINK

Equipment

PTC-200 Thermal Cycler	NAME
MJ Research	BRAND
8252-30-0001	SKU
https://www.gmi-inc.com/product/mj-research-ptc-200-thermal-cycler/?srsItid=AfmBOopNoU-2KhYScE8×1qQsrdebGEXKkaNSnnebjUbsH9pBJbJhxn0J	LI NK



Library Preparation: DNA Shearing

5m

- 1 Prepare  50-100 ng DNA sample into  50 μ L TE buffer in PCR plate.
- 2 Gently pipette up and down to mix. 
- 3 Transfer entire sample volume into microTUBE plate (Covaris, 520078), spin plate at  600 x g, 00:05:00 . 
- 4 Shear DNA to about 450-bp on Covaris E220 Focused-ultrasonicator system (Covaris, Inc. Woburn, MA) with the following setting:

	A	B
	Water level	6
	Power	140
	Duty cycle	10%
	Cycle Per Burst	200
	Treatment time	40 seconds
	Temperature	7°C+/-2°C

Library Preparation: DNA sample purification

13m

- 5 Spin the microTUBE plate and transfer sheared DNA samples into PCR strip tube. 
- 6 Purify sheared DNA by adding  90 μ L (1.8X) AMPure XP beads (Beckman Coulter Catalog # A63882) to tube and pipette up and down 8 times followed by  00:08:00 incubation at  Room temperature . 

7 Put sample tube on magnet stand for  00:05:00 .

5m

8 Remove and discard the supernatant.

9 Wash beads twice by adding  150 μ L 80% ethanol.



10 Dry beads for 3-5 min at  Room temperature .

11 Elute DNA into  22 μ L DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906).

12 DNA QC on Bioanalyzer2100: optional

*

12.1 Dilute sheared DNA 10X using DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906).

12.2 Load  1 μ L diluted DNA on Agilent HS DNA kit and QC on Bioanalyzer 2100.

Mung Bean Nuclease Treatment

45m

13 Dilute Mung Bean Nuclease to  0.25 μ L .

	A	B
	Mung Bean Nuclease (New England Biolabs Catalog # M0250S)	1 μ L
	Mung bean nuclease buffer (New England Biolabs Catalog # B0250S)	4 μ L
	DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906)	35 μ L

13.1 Mix well by pipetting up and down 10 times.





14 Prepare 0.3% SDS Solution.

A	B
DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906)	32 μ L
UltraPure™ SDS Solution, 10% (Thermo Fisher Scientific Catalog # 15553027)	1 μ L

15 Mung bean nuclease reaction.

A	B
Sheared DNA	20 μ L
Mung bean nuclease buffer (New England Biolabs Catalog # B0250S)	2.8 μ L
DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906)	5.2 μ L
diluted Mung bean nuclease (0.25 U/ μ L)	2.0 μ L

15.1 Mix well by pipetting up and down 10 times.



15.2 Incubate reaction at 30 °C 00:30:00 on thermocycler.

30m

15.3 At the end of incubation, add 1 μ L 0.3% SDS immediately and pipetting up and down to mix.

16 Sample Purification.

- 16.1
- Add 60 μ L (2X) AMPure XP beads (Beckman Coulter Catalog # A63882) to sample.
 - Pipetting up and down 8 times to mix.



16.2 Incubate at Room temperature for 00:10:00 , then put on magnet stand for 00:05:00 .

15m



- 16.3
- Wash beads with 150 μL with 80% Ethanol.
 - Dry beads at Room temperature for 3-5 min.
- 16.4 Elute DNA with 11.5 μL DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906).
- 17 **Sample QC:** Optional.
- 17.1 Dilute MBN treated DNA sample 10X and load 1 μL on High Sensitivity DNA Kit (Agilent Technologies Catalog# 5067-4626).

T4 PNK Treatment

30m

- 18 Prepare T4PNK reaction master mix:

A	B
NEBuffer™ 4 (New England Biolabs Catalog # B7004S)	1.5 μL
Adenosine 5'-Triphosphate (ATP) (New England Biolabs Catalog # P0756S)	1.5 μL
T4 Polynucleotide Kinase (New England Biolabs Catalog # M0201S)	1.0 μL
DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906)	1.0 μL

- 19 T4PNK reaction

A	B
MBN-treated DNA	10 μL
T4PNK master mix	5 μL

- 19.1 Gently pipetting up and down to mix.

19.2 Incubate at  37 °C for  00:30:00 on a thermocycler.

30m



A-Tailing

30m

20 Prepare  1 millimolar (mM) equimolar dATP/ddBTPs.

A	B
dATP Solution (New England Biolabs Catalog # N0440S)	2.5 µL
ddTTP Solution Dideoxynucleoside Triphosphate Set (Millipore Sigma Catalog # 3732738001)	25 µL
ddCTP Solution Dideoxynucleoside Triphosphate Set (Millipore Sigma Catalog # 3732738001)	25 µL
ddGTP Solution Dideoxynucleoside Triphosphate Set (Millipore Sigma Catalog # 3732738001)	25 µL
DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906)	172.5 µL

20.1 Mix by gentle vortex.



21 Prepare A-Tailing reaction mix:

A	B
NEBuffer™ 4 (New England Biolabs Catalog # B7004S)	0.5 µL
Klenow Fragment (3'→5' exo-) (New England Biolabs Catalog # M0212L)	1.5 µL
1 mM equimolar dATP/ddBTPs	1.5 µL
DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906)	1.5 µL

21.1 Mix by gently pipetting up and down 6X.



22 A-Tailing reaction



A	B
DNA sample from previous step.	15 µL
A-Tailing reaction mix	5.0

22.1 Mix gently by pipetting up and down.



22.2 Incubate reaction at 37 °C for 00:30:00 on a thermocycler.

30m



Adapter Ligation

20m

23 Prepare Ligation master mix.

A	B
NEBuffer™ 4 (New England Biolabs Catalog # B7004S)	1.74 µL
Adenosine 5'-Triphosphate (10mM ATP) (New England Biolabs Catalog # P0756S)	4.0 µL
CS adapter (Integrated DNA Technologies (IDT) Catalog # 1080799)	0.4 µL
DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906)	9.76 µL

23.1 Mix gently by pipetting up and down 6 times.



24 Adapter ligation reaction

A	B
DNA sample from A-Tailing step	20 µL
Ligation master mix.	15.9 µL
T4 DNA Ligase (New England Biolabs Catalog # M0202L)	1.5 µL

24.1 Mix gently by pipetting up and down 6 times.



24.2 Incubate reaction at  20 °C for  00:20:00 followed by  4 °C hold on a thermocycler.

20m



24.3 Purify sample with 1X ( 37.4 µL) AMPure XP beads (Beckman Coulter Catalog # A63882) and eluted into  25 µL H₂O.

24.4 **Optional:** Dilute  1 µL adapter-ligated DNA 5X and load  1 µL library on High Sensitivity DNA Kit (Agilent Technologies Catalog# 5067-4626).



Note

If high adapter dimer was seen on QC, perform a second AMPure XP beads purification.

Post Adapter-Ligation qPCR-1 Quantification

25 Prepare qPCR master mix.

A	B
Primer Mix (KAPA Library Quant Kit, Roche Catalog # 07960204001)	1 mL
KAPA SYBR qPCR Master Mix (KAPA Library Quant Kit, Roche Catalog # 07960204001)	5 mL
100 µM Nanoseq PCR primer1: 5'-ACACTCTTTCCCTACACGAC-3' (IDT, HPLC grade)	20 µL
100 µM Nanoseq qPCR primer2: 5'-GTGACTGGAGTTCAGACGTG-3' (IDT, HPLC grade)	20 µL

26 Dilute adapter-ligated DNA 20x by adding  2 µL adapter-ligated DNA to  38 µL DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906).



27 qPCR reaction setup.



A	B
20X diluted library	1 µL
NFW Nuclease free water (Thermo Fisher Scientific, Catalog# AM9906)	3 µL
qPCR master mix	6 µL

27.1 Mix by pipetting up and down 10X.



27.2 Set reactions as triplicates in 384 well plates.

28 qPCR Standards.



A	B
Standard 1-6 (KAPA Library Quant Kit, Roche Catalog # 07960204001)	2 µL
NFW Nuclease free water (Thermo Fisher Scientific, Catalog# AM9906)	2 µL
qPCR master mix	6 µL

28.1 Set reactions as triplicates in 384 well plates.

28.2 Mix by pipetting up and down 10X.



29 qPCR amplification.



Run qPCR reactions on QuantStudio™ 6 Flex Real-Time PCR System (Applied Biosystems™) using the following PCR program:

A	B	C	D
Step	Temp	Time	Cycle #



	A	B	C	D
	1	95 °C	5 min	1
	2	95°C	30 s	35
		60°C	45 s	
	3	95°C	15 s	1
	4	65°C	1 min	1
	5	95°C	15 s	1

30 Determining the sample concentration.

The concentration (nM (fmol/μl)) is determined as follows:

mean of sample concentration × dilution factor (20) × 452/573/1,000.

Note

where 452 is the size of the standard in bp; 573 is the proxy for the average library size in bp.

Library PCR

- 31 Based on qPCR quantification result, dilute **1M 0.1 Mass Percent** adapter-ligated DNA in **23.5 μL** using DEPC-Treated Water (Thermo Fisher Scientific Catalog # AM9906).

- 32 Set up library amplification mix.

A	B
adapter-ligated DNA (0.1 fmol)	23.5 μL
IDT Index primer, xGen™ UDI 10nt Primers (IDT Catalog # 10008052)	1.5 μL



A	B
NEBNext® Ultra™ II Q5® Master Mix (New England Biolabs Catalog # M0544L)	25 µL

Oligonucleotides sequences

P7 CAAGCAGAAGACGGCATAACGAGAT[i7-NNNNNNNNNN]GTCTCGTGGGCTCGG

P5 AATGATACGGCGACCACCGAGATCTACAC[i5-NNNNNNNNNN]TCGTCGGCAGCGTC

33 Library amplification program.

A	B	C	D
Step	Temp	Time	Cycle #
1	98 °C	30 s	1
2	98°C	10 s	16*
	65°C	75 s	
3	65°C	5 min	1
4	4°C	Hold	1

*The number of PCR cycles is dependent on the input: 0.1 fmol, 16 cycles; 0.3 fmol, 14 cycles; 0.6 fmol, 13 cycles; 5 fmol, 10 cycles.

34 Library Purification.

34.1 Purify amplified library twice with 0.75X AMPure XP beads (Beckman Coulter Catalog # A63882).

34.2 Elute library into  20 µL  10 millimolar (mM) Tris-HCl Elution buffer (Invitrogen Catalog# A33566).

35 Library Quantification and normalization.

35.1 QC  1 μ L library on DNA 7500 kit (Agilent Technologies Catalog# 5067-1506).

35.2 Dilute library to  5 nanomolar (nM) in  10 millimolar (mM) Tris-HCl Elution buffer (Invitrogen Catalog# A33566) for pooling.

Final Library Quantification by qPCR-2

36 Library dilution.

36.1 Dilute  5 nanomolar (nM) library 3000X by sequential 100 X and 30X dilutions using Nuclease free water (Thermo Fisher Scientific, Catalog# AM9937).

37 Prepare qPCR-2 master mix.



A	B
Primer Mix (KAPA Library Quant Kit, Roche Catalog # 07960204001)	1 mL
KAPA SYBR qPCR Master Mix (KAPA Library Quant Kit, Roche Catalog # 07960204001)	5 mL

38 qPCR-2 reaction setup.



A	B
3000x diluted library or Standard 1-6 (KAPA Library Quant Kit) Roche Catalog # 07960204001	4 μ L
qPCR master mix	6 μ L

38.1 Set reactions as triplicates in 384 well plates.

38.2 Mix by pipetting up and down 10X.





- 39 qPCR-2 amplification.
Carry out qPCR reactions on QuantStudio™ 6 Flex Real-Time PCR System (Applied Biosystems™) using the following PCR program:

	A	B	C	D
	Step	Temp	Time	Cycle #
	1	95 °C	5 min	1
	2	95°C	30 s	35
		60°C	45 s	
	3	95°C	15 s	1
	4	65°C	1 min	1
	5	95°C	15 s	1

Library Pooling

- 40 Pool libraries at equimolar concentrations based on qPCR-2 results and sequence libraries on Illumina platform for 30X coverage.



Protocol references

Federico Abascal, Luke M. R. Harvey, Emily Mitchell, Andrew R. J. Lawson, Stefanie V. Lensing, Peter Ellis, Andrew J. C. Russell, Raul E. Alcantara, Adrian Baez-Ortega, Yichen Wang, Eugene Jing Kwa, Henry Lee-Six, Alex Cagan, Tim H. H. Coorens, Michael Spencer Chapman, Sigurgeir Olafsson, Steven Leonard, David Jones, Heather E. Machado, Megan Davies, Nina F. Øbro, Krishnaa T. Mahubani, Kieren Allinson, Moritz Gerstung, Kourosh Saeb-Parsy, David G. Kent, Elisa Laurenti, Michael R. Stratton, Raheleh Rahbari, Peter J. Campbell, Robert J. Osborne & Iñigo Martincorena Somatic mutation landscapes at single-molecule resolution *Nature*. 2021 May;593(7859):405-410, PMID: 33911282

Margaret L Hoang, Isaac Kinde, Cristian Tomasetti, K Wyatt McMahon, Thomas A Rosenquist, Arthur P Grollman, Kenneth W Kinzler, Bert Vogelstein, Nickolas Papadopoulos; Genome-wide quantification of rare somatic mutations in normal human tissues using massively parallel sequencing. *Proc Natl Acad Sci USA*. 2016 Aug 30;113(35):9846-51, PMID: 27528664