



Aug 24, 2022

Version 2

Optimized protocol for transcriptome analysis of mouse brain endothelial cells V.2


PLOS One

Peer-reviewed method

DOI

[dx.doi.org/10.17504/protocols.io.8epv59or6g1b/v2](https://doi.org/10.17504/protocols.io.8epv59or6g1b/v2)Won-Jong OH¹, Namsuk Kim¹, Mi-Hee Jun¹, Jin-Young Jeong¹¹Korea Brain Research InstitutePLOS ONE Lab Protocols
Tech. support email: plosone@plos.org

Neurovascularlab KBRI



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN ACCESS

DOI: <https://dx.doi.org/10.17504/protocols.io.8epv59or6g1b/v2>External link: <https://doi.org/10.1371/journal.pone.0275036>

Protocol Citation: Won-Jong OH, Namsuk Kim, Mi-Hee Jun, Jin-Young Jeong 2022. Optimized protocol for translome analysis of mouse brain endothelial cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv59or6g1b/v2> Version created by **Neurovascularlab KBRI**

Manuscript citation:

Kim N, Jun MH, Jeong JY, Oh WJ (2022) Optimized protocol for translome analysis of mouse brain endothelial cells. PLOS ONE 17(9): e0275036. <https://doi.org/10.1371/journal.pone.0275036>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 23, 2022

Last Modified: August 24, 2022

Protocol Integer ID: 69048

Keywords: Brain endothelial cells, RiboTag, Translatome, RNA sequencing, pure mrnas from bec, many transcriptome analyses for bec, mouse brain endothelial cells brain endothelial cell, specificity of the mrna pool, hemagglutinin epitope under cre recombinase activation, molecular analysis of bec, mrna, becs at the molecular level, mrna level, optimized protocol for translome analysis, mrna pool, translated protein, diverse neurodegenerative disease, many transcriptome analysis, molecular pathogenesis of these neurological disease, cre recombinase activation, low yield of mrna, pure mrna, bulk gene expression analysis, gene expression, protein, impairment of bec, parenchyma cells from toxin, quality cdna library for rna, translome analysis, rna, using ribotag mice, molecular pathogenesis, ribotag mice, brain homeostasis, hemagglutinin epitope, ribosomal component

Abstract

Brain endothelial cells (BECs) are important conduits that deliver oxygen and nutrients, protect parenchyma cells from toxins, and drain wastes to maintain brain homeostasis. Impairment of BECs has been implicated in diverse neurodegenerative diseases, including Alzheimer's disease and Parkinson's disease. Therefore, molecular analysis of BECs is important for understanding the molecular pathogenesis of these neurological diseases. Even though many transcriptome analyses for BECs have been developed, mRNA levels do not necessarily correlate with the levels of actively translated proteins. Translatome analysis using RiboTag mice, in which Rpl22, a ribosomal component, is tagged by the hemagglutinin epitope under Cre recombinase activation, could serve as an excellent tool that overcomes these caveats. However, implementation of this technique is limited by high noise-to-signal ratios as well as the low yield of mRNAs from BECs, which limits bulk gene expression analysis. In this study, we established a protocol to isolate highly pure mRNAs from BECs in the cortex of eight- to twelve-week-old male *Tie2-Cre; Rpl22^{HA/HA}* mice by using a cell strainer to trap blood vessels prior to immunoprecipitation. According to the results of RT-PCR, the specificity of the mRNA pools isolated by our protocol was much higher than that of the pools isolated by the standard protocol. We were also able to generate a high-quality cDNA library for RNA-seq with the small amount of mRNA isolated with our protocol. Thus, this optimized method will be useful for future studies of BECs at the molecular level.



Guidelines

1. An RNase-free environment is essential. Use barrier pipet tips to avoid RNase contamination. Wipe down the surface of an experimental table and all equipment including surgical tools, pipets, etc., with RNase Zap.
2. Homogenization buffer and high-salt buffer should be freshly prepared.
3. Washes should be done in cold conditions.
4. Tissue samples should be processed fresh directly to RNA yield.
5. The average amount of BEC mRNA \Rightarrow whole cortex (8-12 weeks): 7.3 ng, visual cortex (8-12 weeks): 1.05 ng, visual cortex (2 weeks): 0.25 ng per mouse.

Materials

A	B	C
REAGENT or RESOURCE	SOURCE	IDNETIFIER
Antibodies		
Mouse anti-HA	Millipore	Cat# 05-904 RRID: AB_417380
Chemicals, peptides, and kits		
TRIzol™ Reagent	Thermo	Cat# 15596026
Cycloheximide	Sigma-Aldrich	Cat# 1810
Magnesium chloride	Sigma-Aldrich	Cat# M8266
Potassium chloride	Sigma-Aldrich	Cat# P9333
DNase1	Invitrogen	Cat# 18068015
Pierce™ Protein A/G Magnetic Beads	Thermo	Cat# 88803
Chloroform	Sigma	Cat# C2432
20X TE Buffer (pH 7.5)	Promega	Cat# A2651
Ethyl alcohol, Pure	Sigma	Cat# E7023
Glycogen, Molecular Biology Grade	Roche	Cat# 10901393001
RNasin® Ribonuclease Inhibitor	Promega	Cat# N2115
Pierce™ Protein A/G Magnetic Beads	Thermo Fisher Scientific	Cat# 88803
Halt™ Protease and Phosphatase Inhibitor Cocktail	Thermo Fisher Scientific	Cat# 78444
NEBNext® Single Cell/Low Input RNA Library Prep Kit for Illumina®	NEB	Cat# E6420L
NEBNext® Multiplex Oligos for Illumina®	NEB	Cat# E7600S
High Sensitivity D5000 Screen tape	Agilent	Cat# 5067-5592
High Sensitivity D1000 Screen tape	Agilent	Cat# 5067-5584
High Sensitivity RNA Screentape	Agilent	Cat# 5067-5579



A	B	C
High Sensitivity D5000 Screen tape Reagent	Agilent	Cat# 5067-5593
High Sensitivity D1000 Screen tape Reagent	Agilent	Cat# 5067-5585
High Sensitivity RNA Screen tape Reagent	Agilent	Cat# 5067-5580
High Sensitivity D5000 Screen tape ladder	Agilent	Cat# 5067-5594
High Sensitivity D1000 Screen tape ladder	Agilent	Cat# 5067-5587
High Sensitivity RNA Screen tape ladder	Agilent	Cat# 5067-5581
Glass homogenizer	WHEATON	Cat# 357542
Disposable scalpel	Bard-Parker	Cat# 371611
Experimental models: Organisms/strains		
Mouse: Tie2-Cre	The Jackson Laboratory	Stock# 008863
Mouse: Ai9	The Jackson Laboratory	Stock# 007909
Mouse:RiboTag mice (Rpl22HA/HA)	The Jackson Laboratory	Stock# 011029

Troubleshooting

Safety warnings

⚠️ TRIzol is a highly corrosive and toxic chemical that can cause burns on contact with the skin as well as systemic poisoning.

Chloroform can cause a person to become unconscious and even be fatal at high doses.









Before start

1. An RNase-*free* environment is essential. Use barrier pipet tips to avoid RNase contamination. Wipe down the surface of an experimental table and all equipment including surgical tools, pipets, etc., with RNase Zap.
2. Homogenization buffer and high-salt buffer should be freshly prepared.












Vessel isolation

2h

- 1 The whole mouse cortex of a *Tie2-Cre; Rpl22^{HA/HA}* mouse is isolated in the chilled DMEM. Then, tissues are dissociated by using a glass homogenizer (WHEATON, 357542) in  10 mL of chilled DMEM. 10m
- 2  1000 x g, 4°C, 00:10:00 10m
- 3 After discarding the supernatants, the pellets are resuspended in  15 mL of 20 % BSA-DMEM to avoid myelin contamination. 5m
- 4  2500 x g, 4°C, 00:15:00 15m
- 5  This process is repeated three times. 40m
- 6 After discarding the supernatants, the pellets are resuspended in  5 mL of chilled PBS. 5m
- 7 PBS containing blood vessels is passed through a 40-micrometer cell strainer. 5m

Immunoprecipitation

1d 1h

- 8 The strainer mesh containing vessels is then cut with a disposable scalpel (Bard-Parker, 371611) and transferred into a microcentrifuge tube for lysis in  600 μ L of homogenization buffer containing  1 % (v/v) NP-40,  100 millimolar (mM) KCl,  50 millimolar (mM) Tris ( 7.4),  12 millimolar (mM) $MgCl_2$, cycloheximide ( 100 mg/mL), heparin ( 1 mg/mL), Halt Protease and Phosphatase Inhibitor Cocktail (Thermo Fisher Scientific, 78444), RNA inhibitor (5 units/ml, Promega, N2615), and  1 millimolar (mM) DTT). 5m



- 9 The lysates are incubated On ice for 00:05:00 . 5m
- 10 12000 x g, 4°C, 00:10:00 10m
- 11 After being transferred to a new 1.5 ml microcentrifuge tube, the supernatants are incubated with a mouse monoclonal antibody against the HA epitope tag (1:200, Millipore, 05-904) for 04:00:00 at 4 °C with rotation by using a multimixer (NanoEnTek, 4519). 4h
- 12 Protein A/G magnetic beads (Thermo Fisher Scientific, 88803) equilibrated in homogenization buffer for 30 min are added to the antibody-lysate solution and incubated Overnight at 4 °C with gentle rotation. 16h
- 13 The next day, after a brief spin-dwon, the magnetic beads are washed five times with 1 mL high salt buffer (1% NP-40, 300 millimolar (mM) KCL, 50 millimolar (mM) Tris (7.4), 12 millimolar (mM) MgCl₂, cycloheximide (100 mg/mL), and 0.5 millimolar (mM) DTT). 30m

mRNA isolation

5h 55m

- 14 After the last wash, all supernatants are removed and 1 mL of TRIzol reagent (Invitrogen, 15596026) is added to the bead-antibody-tissue homogenate, followed by 200 µL of chloroform (Sigma-Aldrich, C2432). 10m

Safety information

TRIzol is a highly corrosive and toxic chemical that can cause burns on contact with the skin as well as systemic poisoning.

Chloroform can cause a person to become unconscious and even be fatal at high doses.



- 15 12000 x g, 4°C, 00:10:00 10m
- 16 The upper aqueous layer (approximately 600 µL) is transferred into a new 15 ml conical tube, and 60 µL of 4 Molarity (M) LiCl, 120 µL of 20 X TE (0.2 Molarity (M) Tris-HCl, 20 millimolar (mM) EDTA, 7.5 , Promega, A2651), 1.8 mL of 100% ethyl alcohol (Sigma, E7023), and 3 mL of glycogen (Roche, 10901393001) are added for RNA precipitation. 15m
- 17 The mRNA mixture is incubated Overnight at -20 °C 16h
- 18 The following day, samples are centrifuged 12000 x g, 4°C, 00:10:00 . After the supernatants are discarded, 1 mL of 75% ethyl alcohol is added to the pellets for washing. 10m
- 19 After centrifugation 7500 x g, 4°C, 00:05:00 and subsequent supernatant removal, the samples are air-dried for 00:05:00 at Room temperature . Do not overdry the beads. 10m
- 20 The dried pellets are then resuspended in 16 µL of RNase-free water. 5m
- 21 2 µL of DNase I and 2 µL of 10X DNase I Reaction Buffer (Invitrogen, 18068-015) are added to the reaction mixture, which is then incubated for 00:15:00 at Room temperature . 15m
- 22 DNaseI is inactivated by adding 25 millimolar (mM) of EDTA and heating at 65 °C for 00:10:00 . 10m
- 23 For RNA precipitation, 2.2 µL of 4 Molarity (M) of LiCl, 4.8 µL of 20 X TE (0.2 Molarity (M) Tris-HCl, 20 millimolar (mM) EDTA, 7.5 (Promega, A2651), 66 µL of 100% ethyl alcohol (Sigma-Aldrich, E7023), and 1 µL of 16h



glycogen (Roche, 10901393001) are added to the RNA mixture, followed by

Overnight incubation at -20 °C .

24 The next day, the RNA mixture is centrifuged 12000 x g, 4°C, 00:10:00 10m

25 After removing the supernatants, 1 mL of 75% ethyl alcohol is added to the pellets for washing. 5m

26 After centrifugation 7500 x g, 4°C, 00:05:00 , the supernatants are discarded. 5m

27 The pellets are then air-dried and finally resuspended in 10 µL of RNase-free water. 5m

Generation of cDNA library

2h 25m

28 The amount of isolated mRNA is measured by using High Sensitivity RNA ScreenTape (Agilent, 5067-5579), High Sensitivity RNA ScreenTape Reagent(Agilent, 5067-5580), and a High Sensitivity RNA ScreenTape ladder(Agilent, 5067-5581) from the Agilent 4200 TapeStation System according to the manufacturer's instructions.



29 One nanogram of mRNA obtained from RiboTag immunoprecipitation is reverse-transcribed into cDNA using the NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina (NEB, E6420L) according to the manufacturer's protocol.







30 One nanogram of mRNA is added to the mixture containing 1 µL of NEBNext Single Cell RT (Reverse Transcription) Primer Mix. The 9 µL of the final volume is achieved by adding nuclease-free water.

31 The mixture is incubated at 70 °C for 00:05:00 with the heated lid set to 105 °C for annealing and then held at 4 °C . 5m




32 The RT mixture is prepared in a separate tube as follows On ice ; 5 µL of NEBNext Single Cell RT buffer, 1 µL of NEBNext Template Switching Oligo, 2 µL of NEBNext Single Cell RT Enzyme Mix, 3 µL of nuclease-free water. It is important to vortex the NEBNext Single Cell RT buffer prior to use for optimal performance.





33 The RT mixture ( 11 μL) is combined with the annealed sample ( 9 μL). Mix well by pipetting up and down at least 10 times.

34 The reaction is incubated in a thermocycler with the following steps: the heated lid is set to  105 °C , followed by  01:30:00 at  42 °C and  00:10:00 at  70 °C , and then held at  4 °C



1h 40m

35 The cDNA amplification mix is prepared as follows:  50 μL of NEBNext Single Cell cDNA PCR Master Mix,  2 μL of NEBNext Single Cell cDNA PCR Primer, and  28 μL of nuclease-free water.

36  80 μL of cDNA amplification mix are added to  20 μL of the sample with pipetting.






















37 The reaction is performed in a thermocycler with the following PCR cycling conditions.

	Cycle step	Temperature	Time	Cycles
	Initial Denaturation	98 °C	45 sec	1
	Denaturation	98 °C	10 sec	32
	Annealing	62 °C	15 sec	
	Extension	72 °C	3 min	
	Final Extension	72 °C	5 min	1
	Hold	4 °C		

38 For the next step, the NEBNext Bead Reconstitution Buffer and the SPRI (Solid Phase Reversible Immobilization) beads should be warmed to  Room temperature for at least  00:30:00 before use.

30m



- 39  60 μL SPRI beads are added to the PCR. (mix well by pipetting up and down at least 10 times).
- 40 The samples are incubated for at least  00:05:00 at  Room temperature . 5m
- 41 The samples are placed on the magnetic stand (Promega, Z5342) to separate the beads from the supernatant.
- 42 After  00:05:00 , the supernatant is removed. then,  200 μL of 80% freshly prepared ethanol is added for washing. The samples are Incubated at  Room temperature for  00:00:30 , and then the supernatant is carefully removed and discarded. 5m 30s
- 43  . This process is repeated twice. The samples are air-dried for  00:05:00 at  Room temperature . Do not overdry the beads. 5m
- 44  50 μL of 0.1X TE (diluted from 1X TE buffer) is added to the samples to elute the cDNA from the beads. 2m
The samples are mixed well and incubated for at least  00:02:00 at  Room temperature .
- 45 Next,  45 μL of NEBNext Bead Reconstitution Buffer is added to the cDNA-Bead mixture. Mix well by pipetting up and down at least 10 times and incubate for at least  00:05:00 at  Room temperature . 5m
- 46 The samples are placed on a magnetic stand to separate the beads.
- 47 After  00:05:00 , the supernatant is carefully removed. 5m
- 48 Then,  200 μL of 80% freshly prepared ethanol is added to the tube to wash the beads. After  00:00:30 of incubation at  Room temperature ,  . This process is repeated twice. 30s




- 49 The beads containing cDNA are air-dried for 00:05:00 at Room temperature . 5m
Do not overdry the beads.
- 50 cDNA is eluted from the beads by adding 33 μL of 1X TE. Mix well by pipetting up and down at least 10 times. The sample is incubated for at least 00:02:00 at Room temperature . 2m
- 51 The sample is placed on the magnetic stand. After 00:05:00 of incubation at Room temperature , 30 μL of the solution is transferred to a new tube. 5m
- 52 The cDNA quality and quantity can be assessed by using High Sensitivity D5000 ScreenTape (Agilent, 5067-5592), High Sensitivity D5000 ScreenTape Reagent (Agilent, 5067-5593), and a High Sensitivity D5000 ScreenTape ladder (Agilent, 5067-5594) in the Agilent 4200 TapeStation System.
- 53 40 ng of cDNA is used for Illumina NGS (Next Generation Sequencing) library preparation.
- 54 40 ng of cDNA in 1X TE is mixed with 7 μL of NEBNext Ultra II FS Reaction Buffer and 2 μL of NEBNext Ultra II FS Enzyme Mix in a PCR tube. The final volume of the mixture is brought to 35 μL , and the sample is vortexed for 00:00:05 . 5s
- 55 In a thermocycler, with the heated lid set to 75 $^{\circ}\text{C}$, the following program is performed: 00:25:00 at 37 $^{\circ}\text{C}$ and 00:30:00 at 65 $^{\circ}\text{C}$. 55m
- 56 While the PCR is running, prepare the solution for the next step. NEBNext Adaptor for Illumina is diluted by 25-fold in the NEBNext Adaptor Dilution Buffer.
- 57 The following components should be added directly to the above sample (35 μL). The adaptor should be added separately to each sample (DO NOT premix with ligation master mix and enhancer).





Component	Volume
FS Reaction Mixture	35 µl
NEBNext Ultra II Ligation Master Mix	30 µl
NEBNext Ligation Enhancer	1 µl
NEBNext Adaptor for Illumina (dilluted 1:25)	2.5 µl

- 58 The samples are mixed well by using pipetting the entire volume up and down at least 10 times. The ligation mixture is incubated at 20 °C for 00:15:00 in a thermocycler without the heated lid. 15m
- 59 3 µL of USER Enzyme, a mixture of uracil DNA glycosylase (UDG) and the DNA glycosylase-lyase endonuclease VIII, is added to the ligation mixture and mixed well. Incubate at 37 °C for 00:15:00 . 15m
- 60 For the next step, the NEBNext Bead Reconstitution Buffer and the SPRI beads should be warm to Room temperature for at least 00:30:00 before use. 30m
- 61 57 µL of SPRI beads are added to the PCR reaction. The sample is incubated for at least 00:05:00 at Room temperature . 5m
- 62 The sample is placed on a magnetic stand.
- 63 After 00:05:00 incubation, the supernatant is removed. Then, 200 µL of 80% freshly prepared ethanol is added to the tube. After incubation at Room temperature for 00:00:30 , the supernatant is removed. . This process is repeated twice. 5m 30s
- 64 The beads containing cDNA are air-dried for 00:05:00 at Room temperature . Do not overdry the beads. 5m



65  17 μL of 0.1X TE is added to resuspend the beads. The cDNA-bead mixture is incubated for at least  00:02:00 at  Room temperature .

2m

66 The sample is placed on a magnetic stand. After  00:05:00 ,  15 μL of the cleared solution is transferred to a new PCR tube.



5m

67 The following components are combined into a new PCR tube.

Component	Volume
Adaptor Ligated DNA Fragments	15 μl
NEBNext Ultra II Q5 Master Mix	25 μl
Index Primer / i7	5 μl
Index Primer / i5	5 μl













68 Labelling with dual barcodes is performed by using the following PCR cycling conditions.

Cycle step	Temperature	Time	Cycles
Initial Denaturation	98 °C	30 sec	1
Denaturation	98 °C	10 sec	8
Annealing	65°C	75 sec	
Final Extension	65 °C	5 min	
Hold	4 °C		

69 For the next step, the NEBNext Bead Reconstitution Buffer and the SPRI beads should be warmed to  Room temperature for at least  00:30:00 before use.

30m



- 70 The PCR mixture is resuspended in  45 μL of SPRI beads. The sample is incubated for at least  00:05:00 at  Room temperature . 5m
- 71 The cDNA-bead mixture is placed on a magnetic stand to separate the beads from the supernatant.
- 72 After  00:05:00 , the supernatant is removed and discarded. 5m
- 73  200 μL of 80% freshly prepared ethanol are added to the tube in the magnetic stand.  . This process is repeated twice.
- 74 The beads containing cDNA are air-dried on a magnetic stand for  00:05:00 at  Room temperature . 5m
- 75 The cDNA library is eluted by adding  33 μL of 0.1X TE. Mix well by pipetting up and down 10 times.
- 76 The sample is placed on a magnetic stand. After  00:05:00 ,  30 μL of the sample containing the cDNA library is transferred to a new tube. Libraries can be stored at  -20 $^{\circ}\text{C}$. 5m
- 77 Before NGS, the quality of the final cDNA libraries is checked by using High Sensitivity D1000 ScreenTape (Agilent, 5067-5584), High Sensitivity D1000 ScreenTape Reagent(Agilent, 5067-5585), and a High Sensitivity D1000 ScreenTape ladder (Agilent, 5067-5587) in the Agilent 4200 TapeStation System.