

Dec 23, 2025

Opto-Open Field

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge164yv47/v1>

Asa Mackenzie¹

¹Lund University



Åsa Wallén-Mackenzie

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge164yv47/v1>

Protocol Citation: Asa Mackenzie 2025. Opto-Open Field. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261ge164yv47/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 23, 2025

Last Modified: April 21, 2026

Protocol Integer ID: 235678

Keywords: ASAPCRN, open field opto, open field from serra et al, open field, opto, serra et al

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020600

Abstract

Opto-Open Field from Serra et al 2023

- 1 Mice were individually placed in neutral cages for 3 min in order to recover after connecting the optic cables.
- 2 Mice were subsequently placed in the central zone of the open field arena and allowed to freely explore it for 5 min before starting the test.
- 3 The open-field chamber consisted in a 50 cm, squared, transparent, plastic arena with a white floor.
- 4 The Opto-Open Field test consisted of a 20-min session divided into 5-min alternating time intervals (OFF-ON-OFF-ON).
- 5 Each 5-min interval was then divided into two 2.5-min intervals for analysis.
- 6 The patch cable was connected to a rotary joint, which was attached on the other end to a laser that was controlled by an Arduino Uno card (Arduino).
- 7 During the light-ON trials, blue light was delivered according to the behavioral stimulation protocol (20Hz, 5ms, 5mW).
- 8 Self-grooming and escape behaviors were manually recorded by an experimenter blind to the experimental groups using the EthoVision XT 14.0 tracking software (Noldus Information Technology, The Netherlands).