



May 23, 2024

Orexin A Analysis in Plasma

DOI

dx.doi.org/10.17504/protocols.io.kxygx365zg8j/v1

daniel.dautan daniel^{1,2}, Per Svenningsson^{1,2}, Wojciech Paslawski^{1,2}

¹Department of Clinical Neuroscience, Karolinska Institutet, 171 76 Stockholm, Sweden;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA



Jacquelyn Haytayan

Weill Cornell Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx365zg8j/v1>

Protocol Citation: daniel.dautan daniel, Per Svenningsson, Wojciech Paslawski 2024. Orexin A Analysis in Plasma. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx365zg8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: February 16, 2024

Last Modified: September 24, 2024

Protocol Integer ID: 95308

Keywords: ASAPCRN, Serum, orexin, using novus biologicals elisa kit, novus biologicals elisa kit, plasma mouse plasma analysis, mouse plasma, elisa kit

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: 020608

Abstract
















Mouse plasma analysis of Orexin A using Novus Biologicals ELISA kit (NBP2-80231).

Materials

ELISA Orexin A kit- Novus Biologicals, NBP2-80231, MO, US

Troubleshooting



- 1 Centrifuge samples for  00:15:00 at 1000×g at  4 °C . 15m
- 2 Add  50 µL of standards' working solutions and samples to 96-well plate in duplicates.
- 3 Add  50 µL of biotinylated detection antibody working solution to each well. Cover with the provided plate sealer and incubate for  00:45:00 at  37 °C . 45m
- 4 Decant the solution from each well. Add  350 µL of wash buffer added to each well and allowed to soak for  00:01:00 . 1m
- 5 Decant the solution from each well and pat dry against clean absorbent paper. Repeat this wash step (steps 4-5) 3 times.
- 6 Add  100 µL of HRP conjugate working solution to each well. Cover plate with the sealer and incubate for  00:30:00 at  37 °C . 30m
- 7 The wash step was performed as described above. Add  90 µL of substrate reagent to each well. Incubate for  00:15:00 at  37 °C protected from light. 15m
- 8 Add  50 µL of stop solution to each well. Measure the optical density using a micro-plate reader with absorbance set to 450 nm.