

Sep 20, 2023

Organoid Electroporation using CRISPR RNP method

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvor74xv4o/v1>

mns^{1,2}

¹CRUK Cambridge Institute; ²Cambridge Stem Cell Institute



Nefeli Skoufou

ICR

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvor74xv4o/v1>

Protocol Citation: mns 2023. Organoid Electroporation using CRISPR RNP method. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5qpvor74xv4o/v1>

Manuscript citation:
[10.1242/dmm.050279](https://doi.org/10.1242/dmm.050279)

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 06, 2023

Last Modified: September 20, 2023

Protocol Integer ID: 76464

Keywords: using crispr rnp method organoid electroporation, crispr rnp method organoid electroporation, efficient genome editing, crispr, organoid electroporation, using ribonucleoprotein, genome, rnp, electroporation

Abstract

Organoid electroporation using ribonucleoprotein (RNP) CRISPR based approach for highly efficient genome editing.

Materials

Table 1: WENRAFI media composition. Modified from Fujii et al., 2018. Cell Stem Cell.

	A	B	C	D
	Optimised organoid media (replacement of p38i)- WENRAFI	Stock concentration	Volume	Final concentration
	ADF+++	pure	13.08	
	Wnt3a conditioned medium	pure	25 ml	50%
	R-spo conditioned medium	pure	10 ml	20%
	Primocin (Invitrogen #ant-pm-1)	50mg/ml (500x)	100µl	500 µg/mL
	B-27® Supplement (Invitrogen #17504-044)	50x	1000 µl	1x
	Nicotinamide (Sigma #N0636, in water)	1 M (100x)	500 µl	10 mM
	N-Acetylcysteine (Sigma # A9165, in water)	500 mM (400x)	125 µl	1.25 mM
	A3801 (Tocris #2939, in DMSO,)	5 mM (10,000x)	5 µl	500 nM
	mEGF (Invitrogen Biosource #PMG8043)	100 ng/µl (2,000x)	25 µl	50 ng/mL
	mNoggin (Peprotech #250-38)	100 ng/µl (1,000x)	50 µl	100 ng/mL
	IGF-1 (Biolegend, 590904)	100 ng/µl	50 µl	100 ng/mL
	FGF-2 (Peprotech, #100-18B)	100 ng/µl	25 µl	50 ng/mL

	A	B	C	D
	Total		50 mL	

Table 2: ENAFI media composition

	A	B	C	D	E
	ENAFI (EGF, Noggin, ADF, FGF2, IGF1)	Stock concentration	Final concentration	ENAFI+ Y+Chir (48h before)	ENAFI+ Y+Chir+DMSO (24h before and elec day)
	ADF+++	pure		24010	23697.5
	Primocin (Invivogen #ant-pm-1)	50mg/ml (500x)	500 µg/mL	50	50
	B-27® Supplement (Invitrogen #17504-044)	50x	1x	500	500
	Nicotinamide (Sigma #N0636, in water)	1 M (100x)	10 mM	250	250
	N-Acetylcysteine (Sigma # A9165, in water)	500 mM (400x)	1.25 mM	62.5	62.5
	A3801 (Tocris #2939, in DMSO,)	5 mM (10,000x)	500 nM	2.5	2.5
	mEGF (Invitrogen Biosource #PMG8043)	100 ng/ µl (2,000x)	50 ng/mL	12.5	12.5
	mNoggin (Peprotech #250-38)	100 ng/ µl (1,000x)	100 ng/mL	25	25
	IGF-1 (Biolegend, 590904)	100 ng/ µl	100 ng/mL	25	25



	A	B	C	D	E
	FGF-2 (Peprotech, #100-18B)	100 ng/ μl	50 ng/mL	12.5	12.5
	Y-27632	10 mM	10 uM	25	25
	CHIR99021	10 mM	5 uM	25	25
	DMSO		1.25%		312.5
	Total			25mL	25mL



Organoid Expansion (Day -5)

3h

- 1 Expand organoids as previously described. Aim for 10-20 wells of organoids for sufficient cell numbers, depending on the numbers of conditions you want to test. Feed organoids with WENRAFI media (Table 1).

Media Preparation (Day -2)

10m

- 2 48 h before electroporation, replace the medium with 250 μ l of ENAFI medium supplemented with 5 μ M CHIR99021 and 10 μ M Y-27632 (Table 2).

Media Preparation (Day -1)

10m

- 3 24 h before electroporation, replace the medium with 250 μ l of ENAFI medium supplemented with 5 μ M CHIR99021, 10 μ M Y-27632 and 1.25% (vol/vol) DMSO (Table 2).

Single Cell Dissociation (Day 0)

1h

- 4 Remove the medium from the organoids and add  500 μ L of TrypLE Express supplemented with 10 μ M Y-27632 to each well. Scrape the Matrigel off the bottom of the wells with a 1,000- μ l pipette. Split the organoids in 2-4 15 mL Falcons to have smaller volume for the dissociation process.
- 5 Place the tubes in a  37 $^{\circ}$ C water bath for  00:30:00 . Pipette vigorously every 5 min, 10 times with 10 mL pipette and 10 times with a 1,000- μ l pipette with broken tip.
- 6 Thaw Cas9 and guide  On ice .
- 7 Add basal medium up to 10 ml and centrifuge for 4 minutes at 500g. Combine separate Falcon tubes at this stage to have a bigger pellet.
- 8 Aspirate and discard the supernatant. If pellet is loose, do a second centrifugation step in an ependorf with 500-1000 μ l of media left.

30m



- 9 Aspirate and discard supernatant. Add  500 μL of Opti-MEM media and pipette well to mix.
- 10 Count number of cells with a haemocytometer (take 10 μL). Determine number of conditions (100,000 cells per condition). You will need to include negative control, no Cas9.

Making RNP complex

30m

- 11 Mix  1 μL of Cas9 and  1 μL of guide (1:3.33 ratio), you will need to add 2 μL per condition.

20m

Standard concentration: 5ug True Cut Cas9 v2 (Invitrogen, A36499- 500 μg at 5 $\mu\text{g}/\mu\text{L}$) and 100pmol synthetic guide (Synthego- custom made, supplied 3 nmol lyophilised reconstituted with 30 μL water for 100pmol/ μL).

Make complex and leave  00:20:00 at  Room temperature .

- 12 Spin and pellet correct number of cells before washing with  300 μL of PBS. Centrifuge at 500g for 4 minutes.
- 13 While washing with PBS make P3 suppl buffer (20 μL /reaction) (Lonza, V4XP-3032).

Buffer P3: 16.4 μL and Supplement 1: 3.6 μL for a total of 20 μL per condition (recommended to make at least 10% excess for pipetting error). Supplemented with 10 μM Y-27632, leave at  Room temperature .
- 14 Completely remove and discard the supernatant. Resuspend in  20 μL of P3 buffer supplemented with 10 μM Y-27632 per condition.
- 15 Add  2 μL of RNP complex per condition.
- 16 Mix well and load  20 μL into electroporation chamber (16-well nucleovette strips).



Electroporation

30m

17 Leave ⌚ 00:10:00 at 🌡 Room temperature before electroporation.

10m

18 Perform electroporation on Lonza Amaxa 4D Nucleofector with program DS-138.

19 Incubate at 🌡 37 °C for ⌚ 00:10:00 .

10m

Seeding cells

20m

20 Add 🧴 80 µL of warm ENAFI media+ Y+ Chir+DMSO to each chamber. Remove 100 µL into separate ependorfs. Wash each chamber with another 🧴 100 µL of media to ensure you have taken all cells.

21 Centrifuge at 500g for 4 minutes.

22 Remove and discard the supernatant and suspend the pellet with 20-25 µL Matrigel per well. Set up 2 wells per condition.

23 Place the plate in a 🌡 37 °C incubator for ⌚ 00:10:00 to solidify the Matrigel.

10m

24 Once matrigel has solidified, add 250 µL of ENAFI medium supplemented with 5 µM CHIR99021, 10 µM Y-27632 and 1.25% (vol/vol) DMSO (Table 2) to each well.

Media change (Day +1)

10m

25 Next Day: change media back to WENRAFI (Wnt and Rspo conditioned, Table 1) supplemented with 10 µM of Y-27632.

DNA extraction and screening (Day +7)

5h

26 7 days after electroporation, extract DNA from half or a third of the well using PicoPure DNA extraction kit (Invitrogen, KIT0103). Perform 🌡 65 °C lysis step for ⌚ 03:00:00

3h



hours.

- 27 Perform PCR using primers that span the guide (500-800 bp) and submit for Sanger sequencing in both directions.
- 28 Analyse sanger trace using ICE Synthego. You will need to upload a control trace (No Cas9) for each edited trace.
<https://ice.synthego.com/#/>