



Oct 23, 2021

Version 3

## P1 Kidney Cold-Active Protease Single Cell Dissociation V.3

DOI

[dx.doi.org/10.17504/protocols.io.bzeap3ae](https://dx.doi.org/10.17504/protocols.io.bzeap3ae)

Andrew Potter<sup>1</sup>, Steve Potter<sup>1</sup>

<sup>1</sup>CCHMC



Andrew Potter

CCHMC

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bzeap3ae>

**Protocol Citation:** Andrew Potter, Steve Potter 2021. P1 Kidney Cold-Active Protease Single Cell Dissociation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bzeap3ae> Version created by **Andrew Potter**

**Manuscript citation:**

Psychrophilic proteases dramatically reduce single-cell RNA-seq artifacts: a molecular atlas of kidney development

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 22, 2021

**Last Modified:** October 23, 2021

**Protocol Integer ID:** 54434

**Keywords:** CAP, cold-active protease, bacillus licheniformis, single cell dissociation, kidney, single cell suspension from p1 mouse kidney, active protease single cell dissociation method, p1 mouse kidney, p1 kidney cold, single cell suspension, cell suspension, reduced artifact gene expresion change, cell, artifact gene expression change

## Abstract

Method used to derive single cell suspension from P1 mouse kidneys on ice, generating a cell suspension with greatly reduced artifact gene expression changes and suitable for downstream analysis using 10x Chromium or DropSeq.

## Guidelines

### Storage Conditions of Reagents

A	B
Reagent	Storage Condition
DPBS (ThermoFisher, 14190144)	4°C
1 M CaCl <sub>2</sub>	room temp.
BSA (Sigma, A8806)	4°C
Protease from <i>Bacillus Licheniformis</i> (Sigma, P5380)	Store 100 µL aliquots (100 mg/mL) in DPBS at -80°C
DNAse 1 (StemCell, 07469)	Store 10 µL aliquots (250 U/10 µL) in DPBS at -80°C

### Required Equipment

Equipment	Supplier	Catalog no.
gentleMACS dissociator	Miltenyi	130-093-235

### The protocol workflow is as follows:

- A. Isolate Kidney
- B. Initial digestion: triturate on ice
- C. Perform gentleMACS
- D. Continue triturating on ice
- F. Preparing cells for Chromium/DropSeq

## BEFORE STARTING

**Prepare *Bacillus Licheniformis* enzyme mix just prior to starting dissociation:**

	A	B	C
	Volume (μl)	Reagent	Final concentration
	890	DPBS	1X
	5	1 M CaCl <sub>2</sub>	5 mM
	5	DNase 1 (250 U/10 μL)	125 U / mL
	100	<i>B. Lich</i> (100 mg/mL)	10 mg/mL

**+25 mg tissue / 1 mL enzyme mix**

To prepare 0.01% BSA/PBS:

**Make stock of 10% BSA in DPBS and store at -20 °C. To make PBS/BSA 0.01% aliquot 50 mL of DPBS in 50 mL conical and pipet in 50 μL of 10% BSA stock.**

**Prepare 10% FBS/PBS with heat-inactivated FBS.**

Troubleshooting




1 Extract & isolate P1 kidneys in ice-cold PBS.


2 Mince kidneys on top of petri dish, on ice, using razor blade.

3 Weigh out 25 mg of tissue for each tube of B. Lich. enzyme mix (2 tubes total).

 25 mg

4 Incubate tissue + enzyme on ice for 7 minutes while triturating 15 strokes using 1 mL pipet every 2 minutes set to 700  $\mu$ L - first with tip cut off.

 00:07:00

 00:02:00


5 After 7 minutes, take the digest mix (combine the two tubes) and pipet into Miltenyi C-tube (placed on ice); take C-tube to gentleMACS placed in 4° cold room. Run program brain\_03 two times.


 4 °C

6 After MACS, briefly quick spin the MACS tube (to 300 G) at 4 °C to ensure contents are in the bottom of the tube.


 4 °C

7 Re-suspend and visualize cells using scope by taking small aliquot and using a slide; continue digesting cells in C-tube on ice for 8 additional minutes while triturating every 2 min 15 strokes using a 1 mL pipet.

 00:08:00

 00:02:00

8 Add 3 mL ice-cold 10% FBS/PBS to digest mix in C-tube to inhibit the protease.

 3 mL ice-cold 10% FBS/PBS

9 Transfer digest mix to a 15 mL conical. Spin 300 G for 5 minutes at 4 °C; discard supernatant; re-suspend cell pellet in 2 mL ice-cold PBS/BSA.

5m



4 °C

00:05:00 300 g spin

2 mL re-suspend in PBS/BSA

- 10 Filter re-suspended cells using 30  $\mu$ M filter on sterile 15 mL conical on ice - rinse filter with 8 mL ice-cold PBS/BSA.

8 mL rinse filter with PBS/BSA

- 11 Spin 15 mL conical tube containing filtered cells 300 G for 5 minutes at 4 °C; discard supernatant and re-suspend pellet in 10 mL ice-cold PBS/BSA.

5m

4 °C

00:05:00 300 g spin

10 mL PBS/BSA

- 12 Repeat rinse/spin in ice-cold PBS/BSA.
- 13 Remove supernatant and re-suspend in 1-2 mL ice-cold PBS/BSA.
- 14 Examine using hemocytometer and adjust concentration to 100 cells/ $\mu$ L for DropSeq or 1,000 cells/ $\mu$ L for 10X Chromium.