

Aug 08, 2025

(P34) immunostaining of free floating sections

DOI

<https://dx.doi.org/10.17504/protocols.io.261gek32og47/v1>

Yunming Wu¹

¹Stanford University



Yunming Wu

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261gek32og47/v1>

Protocol Citation: Yunming Wu 2025. (P34) immunostaining of free floating sections. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261gek32og47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working. But there is always room for improvement.

Created: August 08, 2025



Last Modified: August 13, 2025

Protocol Integer ID: 224354

Keywords: performing routine immunostaining, routine immunostaining, floating section, p34, protocol, optional section

Disclaimer

The method used here is widely used in many labs and is not invented in this protocol. This protocol is a brief and concise version of the the commonly used immunostaining protocol for brain slices for use as a routine procedure.

Abstract

This is a protocol for performing routine immunostainings of free floating sections. It includes 2 optional sections for enhancing weak signals.



Materials

- 1X PBS

Dilute 10X PBS 1:10 with ddH₂O.

Filter with 0.24um Vacuum filter.

This solution can be stored under room temperature for 1 year.

- PBST

Add 1 mL TritonX100, 100mL 10X PBS, to 500mL ddH₂O.

Stir mix with a stir bar until tritonX is fully dissolved.

Add water to 1L.

- 1M Citric acid

dissolve 192.1g citric acid in 1L ddH₂O.

- 1M Sodium citrate

dissolve 25.1g sodium dihydrogen citrate in 1L.

- 10X Citrate buffer

In a beaker, add 82 mL 1M sodium citrate

add 18 mL 1M Citric acid

add ddH₂O to 1L.

- 1X Citrate buffer

dilute 10X citrate buffer to 1L with ddH₂O.

- **ABC amplification kit:** VECTASTAIN™ Elite ABC-HRP Kit, Peroxidase (Standard), PK-6100

- **TSA amplification kit:** Alexa Fluor™ 488 Tyramide Reagent, Catalog number: B40953

- HRP Stop buffer:

Make 5M sodium azide by dissolving 3g sodium azide (Sigma S2002) into 10mL ddH₂O.

Make 0.2M Trolox by dissolving 1g Trolox (Sigma 238813) into 20mL H₂O.

Make 1M sodium ascorbate by dissolving 1g sodium ascorbate (Sigma 11140) into 5mL ddH₂O.

Store the stock solution at -20C.

At the day of experiment:

In 10mL PBST, add 10uL 1M sodium ascorbate, 100uL 0.2M Trolox, 10uL 5M sodium azide.



Transcardial perfusion surgery

16h

1



Note

All experiment involving the use of live animals must be reviewed and approved by the IACUC, and comply with NIH and institutional guideline.

2 Weight the animal.

3 Anesthetize animals 0.1 mL 2.5% Avertin per 10 g bodyweight i.p. using a sterile tuberculin syringe (26g needle).

4 Perfuse with 10 mL ice cold 1X PBS.

5 Perfuse with 10 mL ice cold 4% PFA.

6 Postfix the sample in 4% PFA 4 °C Overnight with gentle agitation.

16h

7 Store the samples at 4 °C in 1X PBS supplied with 1 µL of 4% PFA. The samples can be stored for upto 1 week.



Preparing free floating sections with a vibratome

1h 48m

8 Melt 4% low melting point agarose with a microwave briefly and keep it at 37 °C afterwards.

2m

9 Pour a thin layer of melt agarose to the peel-a-way mold.

1m

10 Keep the mold in 4 °C until it turns solid.

5m



- 11 Make a mark on one side of the brain.

5m

Note

This can be done by cutting a small piece off on one side, depending on what brain region will be used.

- 12 Embed the brain sample in the mold with the melt agarose. The A-P axis should be vertical to one edge of the mold. The brain is placed upside down so that bragma and lambda are in contact of the surface.

5m

- 13 Section into $\pm 50 \mu\text{m}$ with a vibratome.

1h

- 14 Place the sections into 24 well plate (for serial sectioning), or 6 well plate (selected sections without serial sectioning).

30m

Antigen retrival (Optional)

16m

- 15

Note

This is an optional step to enhance antibody binding. It is recommended for antibodies with weak signals or no signal.

- 16 Transfer the brain sections into a beaker with enough citrate buffer to cover all sections. e.g., 20 sections in a 500 mL beaker filled with 100 mL citrate buffer.

5m

- 17 Microwave 00:01:00 , or until boiling (depends on the power).

1m

- 18 Wait until the citrate buffer cools down to Room temperature .

10m

Immunostaining

16h 38m

- 19 Transfer the brain sections back to the multi-well plate if the antigen retrival step is performed.

5m



- 20 Dilute Primary antibody in PBST 1:1000. 1m
- 21 Incubate the sections with the diluted primary antibodies at room temperature Overnight with gentle agitation. Each well of 24 well plate uses 100uL. 16h
- 22 Wash PBST 5min X3. 15m
- 23 Dilute Secondary antibody in PBST 1:1000. 1m
- 24 Incubate the sections with the diluted secondary antibodies at room temperature Overnight with gentle agitation. Each well of 24 well plate uses 100uL. 1m
- 25 Wash PBST 5min X3. 15m
- 26 Mount in Ymount and image.

Signal amplification through biotin-avidin-hydrogen peroxidase

2h 23m

- 27
- ### Note
- This step can further enhance weak signals. The secondary antibody must be conjugated with biotin.
- 28 Deactivate endogenous peroxidase activity with 3 % (v/v) H₂O₂ in PBST 01:00:00 with gentle agitation. 1h
- 29 Wash with PBST 00:05:00 . 5m



- 30 Dilute reagent A and reagent B 1:50 in PBST. 1m
- 31 Mix and RT  00:30:00 . 30m
- 32 Add A B mix to sample,  00:30:00 (Optional overnight for complete penetration). 30m
- 33 Wash with PBST  00:05:00 . repeat 3 times. 5m
- 34 Incubate with 1X fluorescently conjugated TSA reagent with 0.0015% H2O2 in PBST in dark,  00:02:00 . 2m
- 35 Wash with HRP stop buffer  00:05:00 X3. 5m
- 36 Wash with PBST  00:05:00 . 5m
- 37 Mount in Ymount and image.