

Sep 03, 2024

Version 1

## Pan-microbial metagenomics protocol V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g717w1gwz/v1>

Adela Alcolea-Medina<sup>1,2</sup>, Luke Blagdon Snell<sup>1,3</sup>, Chris Alder<sup>4,3</sup>, Rahul Batra<sup>3</sup>

<sup>1</sup>King's College London; <sup>2</sup>Synnovis LLP; <sup>3</sup>Guy's & St. Thomas' NHS Foundation Trust;

<sup>4</sup>Oxford Nanopore Technologies



Adela Alcolea-Medina

Infectious research

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.q26g717w1gwz/v1>

**Protocol Citation:** Adela Alcolea-Medina, Luke Blagdon Snell, Chris Alder, Rahul Batra 2024. Pan-microbial metagenomics protocol. protocols.io <https://dx.doi.org/10.17504/protocols.io.q26g717w1gwz/v1>

#### Manuscript citation:

Unified metagenomic method for rapid detection of microorganisms in clinical samples. Commun Med (Lond). 2024 Jul 7;4(1):135. doi: 10.1038/s43856-024-00554-3. PMID: 38972920

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** September 03, 2024

**Last Modified:** September 03, 2024

**Protocol Integer ID:** 106879

**Keywords:** metagenomics, respiratory metagenomics, clinical metagenomics, human depletion

**Funders Acknowledgements:**

**Medical Research Council**

Grant ID: MR/W025140/1

**Medical Research Council**

Grant ID: MC\_PC\_19041

## Disclaimer

This workflow is not accredited.

This mechanical human DNA depletion method is patent pending (PCT/GB2023/051417)

## Abstract

DOI: 10.1038/s43856-024-00554-3



## Materials

### Panmetagenomics protocol

| Reagent                                                                                                                                                                                                                                                                                                                                    | Company                     | Lot | Expiry date/received date | Extraction day |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|-----|---------------------------|----------------|
| Lysis Matrix D 2mL Tube                                                                                                                                                                                                                                                                                                                    | MP Biomedicals              |     |                           |                |
| HL-SAN enzyme                                                                                                                                                                                                                                                                                                                              | ArcticZymes Technolo        |     |                           |                |
| LunaScript® RT SuperMix Kit (E3010)                                                                                                                                                                                                                                                                                                        | New England Biolabs         |     |                           |                |
| Sequenase Version 2.0 DNA Polymerase                                                                                                                                                                                                                                                                                                       | ThermoFisher Scientific     |     |                           |                |
| LongAmp® Taq 2X Master Mix                                                                                                                                                                                                                                                                                                                 | New England Biolabs         |     |                           |                |
| Rapid barcoding Kit SQK-RPB004                                                                                                                                                                                                                                                                                                             | Oxford Nanopore Thechnology |     |                           |                |
| AMPure XP Beads                                                                                                                                                                                                                                                                                                                            | BeckmanCoulte r             |     |                           |                |
| Qubit™ 1X dsDNA High Sensitivity (HS) (Q33230)                                                                                                                                                                                                                                                                                             | ThermoFisher Scientific     |     |                           |                |
| NATtrol™ Respiratory Panel 2.1 (RP2.1) Controls                                                                                                                                                                                                                                                                                            | Zeptometrix                 |     |                           |                |
| Viasure respiratory panel III (PCR for PC)                                                                                                                                                                                                                                                                                                 | Pro-labs                    |     |                           |                |
| Tobacco mosaic virus (PC-0107)                                                                                                                                                                                                                                                                                                             | DSMZ                        |     |                           |                |
| TMV primers:<br>TMV_fwd_aps<br>-Catalog number: 10336022<br>GGATATGTCTAAGTCTGTTGC<br>10629186 Nucleotides (25 nmole)<br>11732013 Desalted574910 W7149 (E12)<br>TMV_fwd_aps<br>574910 W7149 (F01)<br>3 TMV_rvr_aps<br>TMV_rvr_aps<br>-Catalog number: 10336022<br>CAGACAACTCGGGTGCG<br>10629186 Nucleotides (25 nmole)<br>11732013 Desalted | Invitrogen                  |     |                           |                |
| Fast SYBR™ Green Master Mix                                                                                                                                                                                                                                                                                                                | Applied Biosystem           |     |                           |                |

|                                                 |                            |  |  |  |
|-------------------------------------------------|----------------------------|--|--|--|
| Microbial DNA-Free Water<br>(ID: 338132) Qiagen | ThermoFisher<br>Scientific |  |  |  |
| RESPIRATORY SWAB<br>MATRIX NEGATIVE<br>CONTROL  | Vircell                    |  |  |  |
| Ethanol                                         |                            |  |  |  |
| Elution buffer                                  | ONT                        |  |  |  |

## Preparation of quality controls:

- 1 **Positive control:** NATtrol™ Respiratory Panel 2.1 (RP2.1) Controls (Zeptometrix)
- 2 Mix **300µL of control 1** and **300µL of control 2** in an Eppendorf. Label with date and LOT number.
- 3 Vortex for 1 minute and Centrifuge at **1,200xg for 5min** immediately prior to each use.
- 4 Aliquot the volume specified by the manual or automated extraction method and process in parallel with equivalent clinical sample volumes.
- 5 Store the mixed positive controls at 4°C, it can be used for two extractions.
- 6 Run a target PCR (**Viasure respiratory panel III**) anytime a new box is opened and record it on the trend analysis spreadsheet.

## Negative control: Respiratory Swab Matrix Negative Control (Vircell)

- 7 Add **500µL** of Microbial DNA-Free Water to a vial and mix until completely reconstituted as per manufacturer's instructions. Vortex for **30 sec** to dissolve and homogenise completely. Spike with **1µL** of TMV.

## Internal control (IC): Tobacco Mosaic Virus (TMV) PC-0107 (DSMZ)

- 8 Reconstitute a new vial of TMV every Monday with **3mL** of Microbial DNA-Free Water (Qiagen). Vortex for 1 min.
- 9 Aliquot **500µL** of IC into six Eppendorf tubes, labelling each tube with the preparation date and lot number. Store at 4°C.
- 10 Use one aliquot per day to spike your samples and negative control, vortex for 1 min before spiking



- 11 For every new vial reconstituted, TMV target PCR should be performed and Ct values should be recorded for trend analysis.
  - 200ul of neat was extracted directly from the vial (no depletion performed).
  - In our hands the Ct value was 5-6
- 12 Discard any remaining volume from used aliquot each day and any remaining aliquots by Friday.
- 13 **Spiking samples with the internal control (IC):**  
  
Spike samples with **1µL of TMV** in the following samples:
  - Broncho-alveolar lavage (**BAL**)
  - Non-directed Broncho-alveolar lavage (**NBL**)
  - Endotracheal tube (**ETT**) aspirates
  - Sputum (**SPT**)  
Spike with **0.5µL of TMV** the following samples:
  - Pleural fluids (**PF**)
  - Nasal throat swabs (**NTS**)

## Sample preparation:

- 14 -Respiratory samples validated for this protocol are: **BAL, NBL, ETT, NTS, sputum and PF**.
  - Samples processed should be fresh up to four days old storage in the fridge, older than that it could not be guarantee the RNA virus detection.
  - Samples are processed in a different way depending on the sample type.

## Sputum and ETT:

- 15 Mucoïd samples such as **sputum or ETT** should be mucolysed before starting the Human DNA depletion process as follows:  
  
Add approximately **2 mL** sputum into a numbered, sterile, plastic universal container or if 2 mL already in original container then do the following:
- 16 Add equal volume (2mL) of mucolytic agent Mucolyse to the sample container and Vortex in the safety cabinet for **30 sec**.



- 17 Leave until liquefaction at room temperature (RT).
- 18 Gently agitate for a further **15 sec** before to proceed to the human DNA depletion
- 19 If after mucolysing the sample it is still mucoid, pipette **500µL** of the sample and mix with **500µL** of PBS, vortex well and proceed with the human DNA depletion.

## Human DNA depletion:

- 20 The human DNA depletion is three steps: **centrifugation, bead-beating and endonuclease treatment**.  
However, BAL, NBL and ETT aspirates are processed in a different way depending on the appearance of the sample. This is due to the different composition of the sample.
- 21 **For turbid and mucoid samples:** the human DNA depletion is performed as usual, centrifugation followed by bead-beating and endonuclease treatment.
- 22 **Transparent samples with some mucoid debris:** the centrifugation step should be skipped and only bead-beating and endonuclease treatment should be performed.
- 23 **For clear samples with watery appearance:** the human DNA depletion should be avoided and the sample should just be extracted.

## 1) Centrifugation:

- 24 Vortex the sample
- 25 Centrifuge speed **1,200xg for 10 min.**
- 26 Aliquot **500µL** of the supernatant into an Eppendorf tube  
**If centrifugation is being removed then vortex sample well and add 500ul of whole sample.**
- 27 Spike **1µL** of TMV to the sample (as indicated above) and **0.5µL** for processing a PF and NTS. Make sure to vortex the TMV aliquot before use.

- 28 Vortex the sample with the TMV
- 29 Proceed with bead beating host-depletion step.

## 2) Bead beating and depletion of free nucleic acid, followed by extraction:

- 30 Transfer sample from Eppendorf (500ul) to the Matrix **Lysing D tube (MP biomedical<sup>TM</sup>)** after vortexing for 1 min.
- 31 Place the matrix tube into the **TissueLyser LT (Qiagen)**. Bead-beat at **50 osc for 3min**.  
Note: sometimes the sample will become foamy after bead-beating, this is ok.
- 32 Transfer **200µL** of the bead-beaten sample to a new Eppendorf tube, placing the tip at the bottom of the tube to take-up the liquid instead of the foam.
- 33 Add **10µL of HL-SAN** enzyme (ArticZymes Technology).
- 34 Mix for **10' at 37°C at 1000rpm** in the Eppendorf ThermoMixer (Eppendorf<sup>TM</sup>).
- 35 NOTE: The HL-SAN is stopped by the Magnapure extraction when the instrument adds proteinase k. However, not all the extraction instruments/kits add proteinase k. If the extraction instrument/kit does not have this step, add the proteinase K manually as followed:  
  
-Add **10µL of proteinase K** to 200ul of the samples previously treated with HL-SAN.  
-Vortex for 10 sec.  
-Incubate at **56°C for 10 min at 1000rpm**.
- 36 Transfer **200µL to Magnapure 24 (Roche<sup>TM</sup>)** for extraction or proceed using alternative extraction methods validated in your own laboratory
- 37 Use total nucleic acid extraction kit 1.1 with pre-set BAL sample parameters and 50µL elution volume.  
  
Fast pathogen 200 1.1 is used for processing <8 samples



Pathogen 200 3.2 for  $\geq 8$  samples. Settings for other extraction robots will need to be established in user laboratories.

## 1. Microbial nucleic acid processing and PCR amplification

### 38 RT-dsDNA:

39 Transfer **16 $\mu$ L** of the 50 $\mu$ L extract (from Magnapure) to a PCR tube and add **4 $\mu$ L of LunaScript® RT SuperMix**.

40 Mix by flicking the tube and briefly spinning down in a small bench top Eppendorf centrifuge.

41 Place the tube in the thermocycler and incubate at **25°C for 2 min, 55°C for 10 min, 95°C for 1 min**.

42 Prepare the master mix for the double strand synthesis using the Sequenase Version 2.0 DNA Polymerase using the following volumes:

| dsDNA synthesis           | Volume | x ____ samples |
|---------------------------|--------|----------------|
| 5x Sequenase Buffer       | 2ul    |                |
| Microbial DNA-Free Water  | 7.7ul  |                |
| Sequenase Dilution Buffer | 0.9ul  |                |
| Enzyme                    | 0.6ul  |                |
| Total volume              |        |                |

Table 1: Reagents required for DsDNA master mix.

43 Add **11.2 $\mu$ L** of the master mix to the **20 $\mu$ L** of product post RT step. Mix by flicking the tube and briefly spin down.

44 Place the tube in the thermocycler. Incubation at **37°C for 8 min**.



- 45 Prepare stock of **70% (700 $\mu$ L absolute ethanol+ 300 $\mu$ L of Microbial DNA-free water for 1mL)** ethanol in Microbial DNA-Free Water.
- 46 Add **45 $\mu$ L** of resuspended AMPure XP Beads (AXP) into a new Eppendorf tube for each sample.
- 47 Add all of the product post dsDNA synthesis (31.2 $\mu$ L) into the Eppendorf tube containing the beads, mix well by flicking the tube and incubate at RT for **5 min**.
- 48 Briefly spin down the samples and pellet the beads on a magnet until the eluate is clear and colourless. Keep the tubes on the magnet and discard the supernatant.
- 49 Keep the tubes on the magnet and wash the beads with **200 $\mu$ L** of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.
- 50 Repeat the previous step.
- 51 Briefly spin down and place the tubes back on the magnet. Pipette off any residual ethanol with the 10p pipette. Allow to dry for **30 sec**, but do not dry the pellet to the point of cracking.
- 52 Remove the tubes from the magnetic rack and resuspend the pellet in **12 $\mu$ L** of Microbial DNA free water. Spin down and incubate for **2 min at RT**.
- 53 Pellet the beads on a magnet until the eluate is clear and colourless.
- 54 Remove and retain **10 $\mu$ L** of eluate for each sample into clean 1.5 mL Eppendorf DNA LoBind tubes, individually.
- 55 **NOTE:** the sample can be store at -80°C after the RT or dsDNA step

## 2. Library preparation Rapid PCR Barcoding Kit (SQK-RPB004) Oxford Nanopore:

- 56 Add **3 $\mu$ L** of the previous product (dsDNA) to a new PCR tube.

- 57 Add **1µL of FRM** (from ONT **Rapid PCR barcoding kit SQK -RPB004**). Mix by flicking the tube and briefly spin down.
- 58 Place the tube in the thermocycler and incubate at **30°C for 1 min and then 80°C for 1 min.**
- 59 Prepare the master mix for the PCR, adding **20µL of water** and **25µL of LongAmp Taq 2X Master Mix** per sample and **1µL** of the barcode primers to **4µL** of the product obtained after the FRM step.
- 60 Mix by gently flicking the tube and briefly spin down.
- 61 Place in the thermocycler. Incubate @ **95° 3', (95° 15'', 56° 15'', 65° 4') x30, 65° 4'**

## Post PCR:

- 62 Perform Qubit™ 1X dsDNA High Sensitivity (HS) on the PCR products as the manufacturer's instructions.

63

| Sample number | Barcode | Qbit concentration (ng/ul) | Volume (ul) to be added to the pool |
|---------------|---------|----------------------------|-------------------------------------|
|               |         |                            |                                     |
|               |         |                            |                                     |
|               |         |                            |                                     |
|               |         |                            |                                     |
|               |         |                            |                                     |
|               |         |                            |                                     |
|               |         |                            |                                     |

**Table 2:** Data collection for qubit concentrations and volume (µL) per sample being added to the pool.

- 64 **Pool samples as follows:**

Based on the qubit readings, prioritise samples with a qubit reading below 1 ng/µL. For these samples, normalise the qubit concentration of your PCR products to 100 ng/µL.

n=sample

e.g.  $100/(n)0.6 = 166.67$  add the total volume of PCR product.

- 65 Clean up 1:1 by adding the same volume of resuspended AMPure XP Beads (AXP) to the same volume of the pooled samples. for example, **50ul of pool = 50ul of beads**
- 66 Mix well by flicking the tube, incubate for **5 min** in the Hula mixer (rotor mixer) at RT.
- 67 Briefly spin down the samples and pellet the beads on a magnet until the eluate is clear and colourless. Keep the tubes on the magnet and discard the supernatant.
- 68 Keep the tubes on the magnet and wash the beads with **200µL** of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.
- Repeat the previous step.
- 69 Briefly spin down and place the tubes back on the magnet. Pipette off any residual ethanol with the 10p pipette. Allow to dry for 30 sec, but do not dry the pellet to the point of cracking.
- 70 Remove the tubes from the magnetic rack and resuspend the pellet in **15µL** Microbial DNA-Free Trish buffer. Spin down and incubate for 2 min at RT.
- 71 Pellet the beads on a magnet until the eluate is clear and colourless.
- 72 Remove and retain **10µL** of eluate into clean 1.5 mL Eppendorf DNA LoBind tubes, individually. Make sure no beads are remaining.
- 73 -Add RAP **1µL** to **10µL** of template, mix gently by flicking the tube, spin down and incubate at **RT for 5 min**.  
- Add **25.5µL** of LB  
- Add **34µL** of SQB  
- Add **4.5µL** of water  
- Mix gently pipetting up and down

### Flowcell loading:

- 74 Thaw FLB and FLT at RT if not thawed from before.
- 75 Add **30µL** of FLT (vortex and spin before by pipetting) to a FLB (vortex before by pipetting). Mix by vortex.



- 76 Open the GridION device lid and slide the flow cell under the clip. Press down firmly on the flow cell to ensure correct thermal and electrical contact.
- 77 Slide the flow cell priming port cover clockwise to open the priming port.
- 78 After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles.
- 79 Set a P1000 pipette to 200µL and insert the tip into the priming port.
- 80 Turn the wheel until the dial shows **220-230µL** , to draw back **20-30µL** , or until you can see a small volume of buffer entering the pipette tip  
  
**Note:** Visually check that there is continuous buffer from the priming port across the sensor array.
- 81 Load **800µL** of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for five min.
- 82 Gently lift the SpotON sample port cover to make the SpotON sample port accessible.
- 83 Load **200µL** of the priming mix into the flow cell priming port (not the SpotON sample port), avoiding the introduction of air bubbles.
- 84 Mix the prepared library gently by pipetting up and down just prior to loading.
- 85 Add **75µL** of the prepared library to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next.
- 86 Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION or GridION device lid.

### GridION settings:

- 87 Select the position of your flow cell on MinKNOW.



- 88 Type your **experiment name** and **sample ID**.
- 89 Continue to kit selection and select **Rapid PCR barcoding kit SQK -RPB004**.
- 90 Continue to run options and change to **24hr**.
- 91 Continue to analysis and edit options and select **barcode both ends** and **mid read barcode filtering**.
- 92 Continue to output and unselect fast 5 and edit FastQ to **100 reads per file**.

### Bioinformatic pipeline and reporting:

- 93 Curated bioinformatic pipeline will generate automatic reports at 30 min (for organism identification) and 2 hr (for AMR determinants). 16-24 hr reports can be generated manually for SNP typing and genomic alignment.
- 94 This is the link for the bioinformatic pipeline analysis [GitHub - GSTT-CIDR/RespiratoryCmg](#) .  
**NOTE:** A new version will be available soon.

### qPCR for the TMV

- 95 The RT step is performed using the LunaScript master mix explained previously in the step 39.
- 96 Master-mix calculations as followed:



| Reagent                  | Volumen | X__samples |
|--------------------------|---------|------------|
| Fast Syber-green         | 10ul    |            |
| Primer<br>Forward (10uM) | 1ul     |            |
| Primer<br>Reverse (10uM) | 1ul     |            |
| DNAse-free water         | 3ul     |            |
| Template (cDNA)          | 5ul     |            |
| Total                    |         |            |

**Table 3:** Reagents required for TMV qPCR

1094/PHYTO-06-19-0201-FI. Epub 2019 Nov 18. PMID: 31502520.

**Primer sequence for TMV:**

**-TMV\_fwd\_aps:**

Catalog number: 10336022

GGATATGTCTAAGTCTGTTGC

10629186 Nucleotides (25 nmole) 11732013 Desalted.

**-TMV\_rvr\_aps:**

Catalog number: 10336022

CAGACAACTCGGGTGCG

10629186 Nucleotides (25 nmole) 11732013 Desalted.

- 97 1 cycle at 95°C for 20 s, 40 cycles at 95°C for 1 s, 40 cycles at 60°C for 20 s, 1 cycle at 95°C for 15 s, 1 cycle at 60°C for 1 min, and 1 cycle at 95°C for 15 s
- 98 The validation of this method is available on [Unified metagenomic method for rapid detection of microorganisms in clinical samples | Communications Medicine \(nature.com\)](#).