



Mar 20, 2024

## Passaging and Seeding Mouse Embryonic Fibroblasts (MEFs) for Experiments

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr3jdbvmk/v1>

madalynn.erb<sup>1</sup>

<sup>1</sup>Van Andel Institute

Team Lee



Jane Balster

ASAP - Team Lee

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.6qpvr3jdbvmk/v1>

**Protocol Citation:** madalynn.erb 2024. Passaging and Seeding Mouse Embryonic Fibroblasts (MEFs) for Experiments. protocols.io <https://dx.doi.org/10.17504/protocols.io.6qpvr3jdbvmk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 09, 2024

**Last Modified:** May 31, 2024

**Protocol Integer ID:** 95134

**Keywords:** ASAPCRN, seeding mouse embryonic fibroblast, mouse embryonic fibroblast, embryonic fibroblast, seeding mouse, mef

## Abstract

This protocol details the passaging and seeding mouse embryonic fibroblasts (MEFs) for experiments.

## Attachments



[1003-2590.docx](#)

18KB



## Passaging and Seeding Mouse Embryonic Fibroblasts (MEFs) for Experiments

15m

- 1 Passage cells when confluency exceeds 60%.
- 2 Remove media.
- 3 Wash cells with  5 mL prewarmed PBS without Calcium and magnesium. 
- 4 Aspirate and discard PBS.
- 5 Add  3 mL -  5 mL  Room temperature TrypLE Express (enough to ensure complete coverage of cells). 
- 6 Incubate cells at  37 °C until they have detached. 
- 6.1 Check at  00:03:00 -  00:05:00 intervals under an inverted microscope, gently tap flask to dislodge cells if needed.  8m
- 7 Add  5 mL -  10 mL (or at least 1:1 ratio of media : TrypLE) of pre-warmed complete media to flask. 
- 7.1 Pipette up and down to dislodge cells as needed. 
- 8 Transfer cell suspension to conical tube. 
- 9 Spin down at  1300 rpm -  1800 rpm (  300 x g ) for 3 mins -  00:05:00 .  5m



9.1 During this time label one microcentrifuge tube for each cell line being used and add  20  $\mu\text{L}$  of Trypan blue.



10 Check for cell pellet then aspirate and discard supernatant.

10.1 Avoid disturbing the cell pellet.

11 Resuspend cells in appropriate volume of prewarmed complete media.

12 Take  20  $\mu\text{L}$  of resuspended cell mix (gently pipetting up and down 10 times) and add to  20  $\mu\text{L}$  of trypan blue.



12.1 Pipette up and down 10 times to mix cells with trypan blue.



13 Count the cells using a Countess 3 (Thermofisher Scientific).

13.1 Clean the glass hemocytometer slide and add  10  $\mu\text{L}$  of cell / trypan blue mix to chamber A and B.



13.2 Avoid bubbles and dirt.

14 Perform cells count ensuring working in same units.

#### Note

e.g.  $\times 10^5$  cells per mL or  $\times 10^6$  cells per mL.

A=\_\_\_\_\_ B=\_\_\_\_\_ Average A & B =\_\_\_\_\_

15 Once counts are performed mix cell suspension again as cells will have settled to the bottom of the tube by the time counts are done.

Add \_\_\_\_\_  $\mu\text{L}/\text{mL}$  of to \_\_\_\_\_  $\mu\text{L}/\text{mL}$  of media



16 Seed cells in appropriate volumes.



- 17 Record the following:
  - Cell type, strain and genotype\*
  - Number of wells
  - Density per well\*
  - Passage\*
  - Time and Date of Seeding\*
  - Temperature(if not at standard room temp)
- 18 If left over cell are needed, place in an appropriate volume of media and transfer to a fresh flask for future experiments
- 19 Incubate plate and flasks  Overnight (O/N) in  37 °C incubator.

5m

