

Dec 15, 2023

Patient PBMC flow cytometry

DOI

dx.doi.org/10.17504/protocols.io.n92ldm688l5b/v1

rebeccawallings¹

¹University of Florida



rebeccawallings

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldm688l5b/v1>

Protocol Citation: rebeccawallings 2023. Patient PBMC flow cytometry. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldm688l5b/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 14, 2023

Last Modified: December 15, 2023

Protocol Integer ID: 92335

Keywords: cytometry

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020527

Abstract

Patient PBMC flow cytometry

Troubleshooting

- 1 1×10^6 PBMCs were taken for flow cytometry and transferred to a v-bottom 96-well plate (Sigma, CLS3896-48EA) and centrifuged at $300\times g$ for 5 minutes at 4°C .

Cells were resuspended in 50 μL of PBS containing diluted fluorophore-conjugated antibodies and incubated in the dark at 4°C for 20 minutes. Cells were centrifuged at $300\times g$ for 5 minutes at 4°C and washed in PBS $\times 2$. Cells were fixed in 50 μL of 1% paraformaldehyde (PFA) at 4°C in the dark for 30 minutes.

Cells were centrifuged at $300\times g$ for 5 minutes and resuspended in 200 μL FACS buffer (PBS, 0.5 mM EDTA, 0.1% sodium azide).

Cells were taken for flow cytometry on a MACS Quant Analyzer (Miltenyi). A minimum of 100,000 events were captured per sample and data were analyzed using FlowJo version 10.6.2 software (BD Biosciences).

When validating flow cytometry panels and antibodies, fluorescence minus one controls (FMOs) were used to set gates and isotype controls were used to ensure antibody-specific binding.