


Mar 27, 2023

Version 2

Pavlovian Conditioned Approach V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v9yjsx4g3e/v2>



lukehann¹, cpking¹, Paul Meyer¹

¹University at Buffalo



lukehann Hannan

SUNY The University at Buffalo

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9yjsx4g3e/v2>

Protocol Citation: lukehann , cpking , Paul Meyer 2023. Pavlovian Conditioned Approach. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.x54v9yjsx4g3e/v2>Version created by **lukehann Hannan**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 27, 2023

Last Modified: March 27, 2023

Protocol Integer ID: 79501

Keywords: incentive motivational property, association with reward, motivational property, incentive, reward, food cue, incentive value, pavlovian conditioned approach, cue, associated stimuli, behavior, construction of the pavca apparatus, pavca apparatus, association, tracker, low tendency, pavca, more info on pavca

Abstract

Through their association with rewards, reward associated stimuli “cues” can acquire incentive motivational properties (aka incentive value), and thus the ability to attract, reinforce, and motivate behavior (Robinson and Berridge 1993; Robinson et al. 2013). However, there are substantial individual differences in the ability of cues to acquire this incentive value. Through the use of a Pavlovian Conditioned Approach (PavCA) paradigm, we have identified subgroups of individuals that differentially attribute incentive value to a food cue. Below, we describe the construction of the PavCA apparatus using modular hardware purchased from Med-Associates, and the procedures used to identify subgroups with high, medium, or low tendencies to attribute incentive value to a food cue (“sign-trackers”, “intermediates”, and “goal-trackers”, respectively.)

OpenBehavior Link (contains info about chamber construction):

<https://edspace.american.edu/openbehavior/project/pavca/>

GitHub Link (contains all necessary programs):

<https://github.com/MeyerLabScience/Pavlovian-Conditioned-Approach>

More info on PavCA:

[Meyer et al. \(2012\)](#)

[King et al. \(2021\)](#)

Guidelines

All programs used during this procedure can be found on the linked GitHub folder. The custom-built chamber's description can be found in the linked OpenBehavior article, which contains a PDF with all specifications.

All conditioning chambers are tested daily to ensure that the equipment is functional.

Materials

All programs used during this procedure can be found on the linked GitHub folder. The custom-built chamber's description can be found in the linked OpenBehavior article, which contains a PDF with all specifications.



Before start

Prior to beginning any pre-training or experimentation, allow any subjects at least 14 days to habituate to their home cages. Within this period, we give all subjects ad-lib food and water. Additionally, provide a reverse day/night cycle during this habituation to ensure maximum efficiency during testing.



Reward Habituation

2d

- 1 Give each subject ~25 banana flavored pellets/day in their home cage for the 2 days preceding the start of the Pre-training.

Pre-Training

1d

- 2 Pre-training consists of a single, ~17.5 minute session (time varies between banana pellets on a VI-30s schedule).
- 3 Prior to testing, download the .txt file named "PavCA Pre-training.txt" from GitHub and create a program within MED-PC V. Copy the text in this file into the program. Save and compile this program. 
- 4 Go to the computer that operates the conditioning chambers, and ensure PavCA Pre-Training is compiled in Trans. Then, open Med-PC-V (or Med-PC IV for older testing setups).
 - 4.1 Select "Sessions" in the upper left corner, and then "Load Box" (Or press Ctrl + O). A "Load Box" window broken down into 9 columns should appear.
 - 4.2 For all boxes that will contain subjects to pre-train, go to the "Procedure" column, click on the unit that aligns with the desired "Box" # (Chamber #), and select the program containing the aforementioned "PavCA Pretraining.txt" text.
 - 4.3 Select "Ok" at the bottom right of the "Load Box" window.
- 5 Remove subject from their home cage and place in transportation cart.
- 6 Transport the rats to the testing room.
- 7 Place subjects into their assigned testing chambers.
 - 7.1 Ensure both the operant chamber's door and the sound attenuating chamber's door are fully closed.



8 At this point the subjects are prepared to begin, and the program is ready to run. Select the traffic light icon at the top of the MED-PC V window. A pop-up window named "Send Signals to Boxes" will appear.

8.1 Select all boxes that contain subjects, select "Start Command" or "Start" and press "Issue" at the bottom right hand of the "Send Signals to Boxes" window.

9 The subjects are now undergoing the pre-training procedure.
Description: Within this session the subjects are afforded 5 minutes in the chamber to familiarize themselves, during which the red house-light is off. After this, the light turns on signaling the start to the pre-training period. Throughout this period, the rats are delivered 25 banana-flavored pellets on a 1-30s variable interval (VI) into the food cup. This ensures that they are familiar with the location at which the reward is delivered, along with the general auditory experience they should expect while in the chamber. There are no cues (levers) provided during this session.

17m 30s

10 After the session is done, ensure that within the "Box Status" section the "Progress" column lists all of the boxes/chambers as "Closed".

11 Remove all subjects from their testing chambers and place into a transportation cart.

12 Transport subjects back to their home cage room.

13 Place subjects back into their assigned home cages.

14 Before running additional rats in the same chambers, clean all walls and flooring with a 70% ethanol solution. Ensure any droppings are removed as the floor bars are often too close for them to fall through to the catch-tray.



Conditioning Sessions

5d

15 The day following pre-training, the conditioning sessions begin. This occurs for 5 consecutive days.

16 Prior to collecting any subjects, download the .txt file named "PavCA Conditioning.txt" from GitHub and create a program within MED-PC V. Copy the text in this file into the program. Save and compile this program.



17 Go to the computer that operates the chambers.



- 17.1 Select "Sessions" in the upper left corner, and then "Load Box" (Or press Ctrl + O). A "Load Box" window broken down into 9 columns should appear.
- 17.2 For all boxes containing subjects to pre-train, go to the "Procedure" column, click on the unit that aligns with the desired "Box" # (Chamber #), and select the program containing the aforementioned "PavCA Conditioning.txt" text.
- 17.3 Select "Ok" at the bottom right of the "Load Box" window
- 17.4 Return to the computer, and direct your attention to the "Box Status" section of the MED-PC hub. Confirm that all boxes/chambers are listed as "Closed" under the "Progress" column.
- 18 Remove subjects from home cages and place them in a transportation cart.
- 19 Transport the rats to the testing room.
- 20 Place subjects in their assigned testing chambers.
- 20.1 Ensure both the operant chamber's door and the sound attenuating chamber's door are fully closed.
- 21 At this point the subjects are prepared to begin, and the program is ready to run. Select the traffic light icon at the top of the MED-PC V window. A pop-up window named "Send Signals to Boxes" will appear.
- 21.1 Select all boxes that contain subjects, select "Start Command" and press "Issue" at the bottom right hand of the "Send Signals to Boxes" window.
- 22 The subjects are now undergoing the conditioning procedure.
Description: After 1 minute in the chamber, the red house light turns on signaling the start of the procedure. From here, 25 banana-flavored pellets were delivered into the food cup on a 30-150s variable interval (VI). For 8 seconds prior to the pellet being delivered a lever is extended and a white LED light located above the lever is illuminated. After these 8 seconds, the lever retracts, the light turns off, and a pellet is distributed into the adjacent food cup. During this session, the number of lever deflections as well as food cup entries (produced when the photo-beam spanning the cup is broken by rats' heads)

33m 30s

are collected. The delivery of food pellets occurs independent of the subjects in-chamber behavior.

23 After ensuring all subjects are finished with their session (session time is no longer increasing and the "Box Status" lists them as closed), remove all subjects from their testing chambers and place into a transportation cart.

24 Transport rats to their home cage room.

25 Place subjects back into their assigned home cage.

Data Analysis

26 Make sure to download the MPC2XL program sold by Med Associates. This will be used to export the data collected during conditioning to Microsoft Excel for analysis.

26.1 This program requires a specific, unique row profile for each experiment. Download the "PavCA Med Associates Row Profile" .txt file.

27 Make sure to download the MPC2XL program sold by Med Associates. This will be used to export the data collected during conditioning to Microsoft Excel for analysis.

28 Open the MPC2XL program and a blank Microsoft Excel sheet.

29 In the "Med-PC to Excel" window, under the "Profile" section, click the "Select" button.

30 Select the previously mentioned "PavCA Med Associates Row Profile" and open the file.

31 Under the "Transfer" section, check the boxes labeled "Column Labels" and "Data" and click the "Transfer!" button.

32 In the "Specify Data Files to Transfer" window that pops-up, navigate to the location where your data is stored and select all data files (no limit) you wish to perform the analysis on.



- 33 Select the cell that you wish to be the top left data containing cell in your excel sheet. DO NOT DO THIS WITH AN EXCEL SHEET ALREADY CONTAINING DATA AS THE TRANSFER IS IRREVERSIBLE AND ANY DATA OVERWRITTEN CANNOT BE RECOVERED.
- 34 In the "Specify Data Files to Transfer" window, select "Done" and all data will be copied into your excel sheet as directed by the row profile.
- 35 Download the file "PavCA Microsoft Access Data Calculator.accdb" from github and open the file using Microsoft Access.
- 35.1 IGNORE THE DATA ALREADY IN THE FILE, THIS IS PLACEHOLDER DATA TO MAINTAIN THE CELLS CALCULATIONS.
- 36 In the left column under the "Unassigned Objects" section, select "PavCA Raw".
- 37 Go to the tab labeled "External Data" in the tool bar.
- 37.1 Select "New Data Source" in the top left corner, then select "From File" and finally, select "Excel".
- 37.2 From here, the program will open a window titled "Get External Data - Excel Spreadsheet" in which you can navigate to and open the excel sheet created in steps 28 - 34.
- 37.3 Make sure the option titled "Import the source data into a new table in the current database." is selected before clicking "OK".
- 38 A window titled "Import Spreadsheet Wizard" will open.
- 38.1 Make sure the box labelled "First Row Contains Column Headings" is checked, then select "Next >."
- 38.2 Select "Next >" on the following section.
- 38.3 Make sure the point labelled "Let access add primary key" is selected in blue, then click "Next >."



- 38.4 Name the file and select "Finish."
- 39 You should now see a second unassigned object underneath "PavCA Raw" with the title given in step 38.4.
- 40 Open both unassigned objects (you will see them in the tabs below the tool bar if they are both open).
- 41 Right click the "PavCA Raw" tab and select "Design View".
- 41.1 This will show each column in the sheet and the type of value contained.
- 42 Scroll to the bottom of the "PavCA Raw" data in design mode.
- 43 Using the black squares located to the left of the column, select all rows from "Lever Press Probability" to "PavCA Index" and copy. (All of these rows' data type should be listed as "Calculated")
- 44 Switch to the tab containing the desired collected data imported from Excel. Right-click on the tab and select design mode.
- 45 Scroll to the bottom of the data and paste the copied lines at the bottom of the data.
- 46 Right-click on the tab containing your data and select "Datasheet View"
- 46.1 Make sure to save the changes when prompted or else no values will be calculated.
- 47 Now, you have calculated a PavCA index which can be found at the end of the dataset.
- 48 From here, in the "Unassigned Objects" column, right-click the data set just created, and select export to export the data into the desired form (multiple options).

