



Feb 24, 2025

PBMCs isolation with CT Tubes (8mL Whole blood)

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l291qp1y/v1>

Sara Lucas Del Pozo^{1,2}, Kai-Yin Chau^{1,2}, Jane Macnaughtan^{1,2}, Giuseppe Uras^{1,2}, Roxana Mezabrovschi^{1,2}

¹University College London; ²Aligning Science Across Parkinson's



Sofia Koletsi

University College London, University of London

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l291qp1y/v1>

Protocol Citation: Sara Lucas Del Pozo, Kai-Yin Chau, Jane Macnaughtan, Giuseppe Uras, Roxana Mezabrovschi 2025. PBMCs isolation with CT Tubes (8mL Whole blood). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l291qp1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 06, 2025



Last Modified: February 24, 2025

Protocol Integer ID: 119991

Keywords: ASAPCRN, pbmcs isolation with ct tube, pbmcs isolation, isolation of peripheral blood mononuclear cell, cpt tube, peripheral blood mononuclear cell, pbmc, ct tube, whole blood, isolation

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000420

Abstract

This protocol describes the isolation of peripheral blood mononuclear cells with CPT tubes.

Guidelines

The protocol needs prior approval by the users' Institutional Review Board (IRB), Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee(s) as applicable.

Safety warnings

- ! The protocol needs prior approval by the users' Institutional Review Board (IRB), Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee(s) as applicable.

Before start


The protocol needs prior approval by the users' Institutional Review Board (IRB), Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee(s) as applicable.

Isolation of human peripheral blood mononuclear cells (PBMC) (adapted from Hugues et al 2020, DOI:10.21769/BioProtoc.3572):

35m

- 1 Perform venepuncture with the arm in a downward position. Collect venous blood into 8 ml BD Vacutainer Cell Preparation Tubes (CPT), keeping tubes upright following blood collection. Immediately prior to centrifugation, invert tubes 8 times to remix the blood sample.
- 2 Centrifuge at  1800 x g, Room temperature, 00:20:00 ( 21 °C) in a swing bucket centrifuge with acceleration and deceleration set at 9.
- 3 Collect PBMC layer with sterile glass pipette (whitish layer just beneath the plasma) and transfer to a 15 ml tube.

20m






Figure 1. Empty CPT (left), after blood draw (middle), and after centrifugation (right). Location of gel plug and sample layers after centrifugation are shown.

- 4 Add room temperature 1x DPBS to bring volume to 🧴 15 mL and invert tube 5 times, followed by centrifugation for 🌀 300 x g, Room temperature, 00:15:00 . 15m
- 5 Aspirate supernatant and resuspend pellet in 🧴 10 mL of warmed (🌡️ 37 °C) culture media/DPBS for cell counting.
- 6 **Take a sample for cell counting:** Dilute PBMC suspension 1:1 in Trypan blue, transfer mixture to a cell counting chamber slide, and count PBMCs using an automated cell counter.
- 7 Determine the required volume to obtain 1×10^6 cells.

Protocol PBMC isolation from Puleo et al; 2017 (DOI:10.21769/BioProtoc.2103)

40m

8 Sample preparation



Mix the blood by gentle inversion several times to ensure a homogeneous suspension (manufacturer recommends 8 times).



- 8.1 Centrifuge the tubes at  1800 x g, Room temperature, 00:20:00 (brake can be on). Be sure that the tubes can clear the rotor in swinging-bucket configurations; some tube locations can cause collision with the rotor and result in broken tubes.

20m

**Note**

- If shipping to another location, wrap in protective foam and ship via same-day or overnight courier to the processing lab at  Room temperature .
- Shipping sample at  4 °C (or on wet ice), can result in platelet activation and other unwanted physiological and phenotypic changes.

9 PBMC isolation

In a hood, gently invert the CPT tubes several times (10 times to resuspend PBMCs and plasma), then pipette the plasma and PBMCs into a sterile 50 ml conical tube.

- 9.1 Fill tube containing the plasma and PBMCs with PBS for a final volume of  50 mL in conical. Centrifuge at  250 x g, Room temperature, 00:10:00 with the brake on.

10m



- 9.2 Aspirate the supernatant carefully to not disturb cell pellet.

- A Pasteur pipette and vacuum can be used for aspirating as well as manually aspirating with serological pipettes and a pipette gun.

- 9.3 Gently re-suspend the cells in  10 mL of PBS and remove an aliquot for cell counting. Use a  10 µL sample for a microscope using trypan blue exclusion, or amount needed for automated cell counter.

- 9.4 After removing the aliquot for counting, fill the tube to final volume of  50 mL with PBS and centrifuge for  250 x g, 00:10:00 with the brake on.

10m



**Note**

After that step, the protocol recommends to discard the supernatant and resuspend the pellet in the appropriate volume of freezing media. Cells should be frozen in approximately 10^7 cells per ml of freezing media.

PBMC isolation with CPT (Chen et al. BMC Immunology (2020))

1d 0h 45m

<https://doi.org/10.1186/s12865-020-00345-0>:

- 10 Collect the whole blood into an 8 mL sodium heparinized CPT vacutainer and invert multiple times to ensure homogenization of the sodium heparin anticoagulant and blood.
- 11 Centrifuge the vacutainer at  1700 x g, 21°C, 00:20:00 , resulting in the separation of contents into layers: the upper layer containing plasma with a cloudy band of PBMCs, the middle layer containing a thick polyester resin, and the lower layer containing erythrocytes and granulocytes. 
- 12 After centrifugation, gently invert the CPT 10 times to resuspend PBMCs and plasma, then decant into a 15 mL conical tube pre-filled with  8 mL of Dulbecco's Phosphate-Buffered Saline (PBS, Thermo Scientific).
- 13 Mix the capped 15 mL conical tube by inversion and centrifuge at  300 x g, 4°C, 00:15:00 . 

- 14 Carefully aspirate the supernatant, then resuspend cell pellet in  10 mL of fresh PBS, and centrifuge for a subsequent  300 x g, 4°C, 00:10:00 . 

- 15 Following this centrifugation step, carefully aspirate the supernatant again without disturbing the cell pellet, and resuspend the pelleted PBMCs in  3 mL of fetal bovine serum (FBS, HyClone).
- 16 Aliquot five hundred microliter of cells into six 2 mL cryovials, each pre-filled with  500 μ L of freezing medium composed of FBS and 20% DMSO (Sigma Hybri-Max D2650) for a final DMSO concentration of 10%. 

- Place the filled cryovials in CoolCell freezing containers (Biocision) at -80°C for 24:00:00 prior to transfer to liquid nitrogen for long-term storage. Retain residual volume of PBMCs in FBS and use for counting and viability assessment.

Manufactures instructions (BD):

1d 0h 50m


- 17 Same first centrifugation as previous protocols (Mix the blood by gentle inversion 8-10 times to ensure a homogeneous suspension and centrifuge the tubes at

20m

1800 x g, Room temperature, 00:20:00 (brake can be on).

- 18 After centrifugation, mononuclear cells and platelets will be in a whitish layer just under the plasma layer. Two methods:

Layering of Formed Elements in the BD Vacutainer® CPT™ Tube



- 18.1 Aspirate approximately half of the plasma without disturbing the cell layer. Collect cell layer with a Pasteur Pipette and transfer to a 15 mL size conical centrifuge tube with cap. Collection of cells immediately following centrifugation will yield best results.

- 18.2 An alternative procedure for recovering the separated mononuclear cells is to resuspend the cells into the plasma by inverting the unopened BD Vacutainer® CPT™ Tube gently 5



to 10 times. This is the preferred method for storing or transporting the separated sample for up to  24:00:00 after centrifugation. To collect the cells, open the BD Vacutainer® CPT™ Tube and pipette the entire contents of the tube above the gel into a separate vessel.

19 Suggested Cell Washing Steps:



Add PBS to bring volume to  15 mL . Cap tube. Mix cells by inverting tube 5 times.

19.1 Centrifuge for  300 rcf, 00:15:00 . Aspirate as much supernatant as possible without disturbing cell pellet.

15m



19.2 Resuspend cell pellet by gently vortexing or tapping tube with index finger.

19.3 Add PBS to bring volume to  10 mL . Cap tube. Mix cells by inverting tube 5 times.



19.4 Centrifuge for  300 rcf, 00:10:00 . Aspirate as much supernatant as possible without disturbing cell pellet. Resuspend cell pellet in the desired medium for subsequent procedure.

10m



20 Average number of mononuclear cells (Lymphocytes & Monocytes) recovered per millilitre of whole blood for each method was:

5m

- BD Vacutainer® CPT™ 1.30×10^6 cells vs FICOLL™ Hypaque™ 1.40×10^6 cells.
- PRP isolated from this tubes: can be snap-frozen  00:05:00 in Liquid nitrogen and transferred to -80° .