

May 15, 2024

Version 2

PBS with 0.1% Sodium Azide V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vz8q28vx1/v2>

Allen Institute for Brain Science¹

¹Allen Institute

Allen Institute for Brain Science
Tech. support
[Click here to message tech. support](#)



Allen Institute

Allen Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vz8q28vx1/v2>

Protocol Citation: Allen Institute for Brain Science 2024. PBS with 0.1% Sodium Azide. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vz8q28vx1/v2> Version created by **Allen Institute**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 15, 2024

Last Modified: May 15, 2024

Protocol Integer ID: 99872

Keywords: RP0068, storage, processing, fixed tissue, preparing pb, pb, tissue,

Abstract

This protocol describes the method for preparing PBS with 0.1% Sodium Azide for storage and/or processing of fixed tissue.

Note: Research reported in this publication was supported by the National Institute Of Mental Health of the National Institutes of Health under Award Number U19MH114830. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Attachments



MSDS_Sodium_Azide.pdf

59KB



RP0068_PBS_with_0.1_..

26KB

