

May 13, 2020

PCR cleanup

DOI

dx.doi.org/10.17504/protocols.io.bgaqjsdw

Alex Zegarra¹

¹BYU



Alex Zegarra

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bgaqjsdw

Protocol Citation: Alex Zegarra 2020. PCR cleanup . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bgaqjsdw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 12, 2020

Last Modified: May 13, 2020

Protocol Integer ID: 36912



- 1 Column Equilibration
- 2
 - A. Insert DNA Mini Column into a 2 mL Collection Tube
- 3 Add 100 uL 3M NaOH to the DNA Mini Column
- 4 Centrifuge at maximum speed for 30-60 seconds
- 5 Discard the filtrate and reuse the collection tube
- 6 Rinse with 200 uL of Neutralization Buffer (P3) or 3M sodium acetate pH=4.5
- 7 2. Save 2 microliters of your PCR reaction mix in a separate tube for a pre-cleanup sample in your gel. Dilute the rest of your PCR reaction mixture 20 fold (22 uL of PCR reaction to 440 uL) in Mini-Prep PB buffer. Transfer the entire diluted PCR reaction into the DNA Mini Column
- 8 Centrifuge at maximum speed for 1 minute
- 9 Save the filtrate in a separate tube and reuse the collection tube
- 10 Add 500 uL PB Buffer
- 11 Centrifuge at maximum speed for 1 min
- 12 Save the filtrate in the tube from step 9 and reuse collection tube
- 13 Add 700 uL PE



NOTE: (ENSURE YOU ADDED THE 100% ETHANOL TO THE PE WHEN YOU MADE IT)

- 14 Centrifuge at maximum speed for 1 min
- 15 Discard the filtrate and reuse the collection tube

OPTIONAL: (REPEAT STEP 13-15 FOR A SECOND DNA WASH BUFFER (PE))

- 16 Centrifuge the empty DNA Mini Column for 2 minutes at maximum speed to dry the column matrix

note: (It's important to dry the DNA Mini column matrix before elution. Residual ethanol may interfere with elution and

- 17 Transfer the DNA Mini Column to a clean 1.5 mL microcentrifuge tube
- 18 Add 40 uL Elution Buffer or sterile deionized water directly to the center of the column membrane

Note: (Efficiency of eluting DNA from the DNA Mini Column is dependent on pH. If using sterile deionized water, makes sure that the pH is around 8.5)

- 19 Let sit at room temperature for 1 min
- 20 Centrifuge at maximum speed for 1 min.

Note: (This represents approx 70% of bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration) (Concentration decreases as you yield more-start with enough for your desired yield)

- 21 Store DNA at -20 C, or use immediately in a gel. If you see your PCR amplicon in the gel, throw away the wash you saved. If you know it amplified (because of the 2 microliters you saved) but don't see it in your elution, then it may be in the wash. You can add more HBC buffer to the wash and try 1 more time.