

Oct 25, 2020

## PCR (Error-prone PCR)

DOI

[dx.doi.org/10.17504/protocols.io.bnxfm6](https://dx.doi.org/10.17504/protocols.io.bnxfm6)

Zhujun Wei<sup>1</sup>

<sup>1</sup>2020 iGEM NEFU China

2020 iGEM NEFU China



Zhujun Wei

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bnxfm6](https://dx.doi.org/10.17504/protocols.io.bnxfm6)

**Protocol Citation:** Zhujun Wei 2020. PCR (Error-prone PCR). [protocols.io https://dx.doi.org/10.17504/protocols.io.bnxfm6](https://dx.doi.org/10.17504/protocols.io.bnxfm6)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** October 25, 2020

**Last Modified:** October 25, 2020

**Protocol Integer ID:** 43729

### Safety warnings

- Please wear gloves during experiments. Don't touch the lid after PCR program initiation.  
Synthesize primers in advance before starting.

- 1 Set up a small box with ice, put the tubes of DNA, 2 x Mut Random System, Mut Enhancer and ddH<sub>2</sub>O into it before loading them into the Bio-rad S1000TM Thermo Cycler.
- 2 Add the following reagents to a PCR tube (20 µL).

Ingredient	Volume
Template Plasmid(1-10 ng/µL)	0.4 µL
2 x Mut Random System	10 µL
Forward Primer (10 µM)	0.4 µL
Reverse Primer (10 µM)	0.4 µL
Mut Enhancer	0-20 µL
ddH <sub>2</sub> O	To 20 µL

- 3 Program the thermocycler as follows:

Temperature	Time
95°C	2 min
94°C	30 s
55°C	1 min
72°C	1 min/kilo base pairs
72°C	7 min
4°C	∞

Note: Higher initial template concentrations can lead to lower mutation rates. More amplification cycles may cause higher mutation rates.

- 4 Build a mutant library, expand the strains, and add the corresponding inducers, collect data using a Microplate reader operating procedure V.1, and use the GraphPad Prism for data processing and analysis.
- 5 Select the positive results for DNA sequencing to determine the mutated nucleotides in the DNA promoter or coding sequence.

