

Oct 15, 2019

PCR (Toehold) Fragment Protocol for a 25ul reaction

DOI

<https://dx.doi.org/10.17504/protocols.io.786hrze>

Janet Standeven¹

¹Lambert iGEM



Lambert Igem

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.786hrze>

Protocol Citation: Janet Standeven 2019. PCR (Toehold) Fragment Protocol for a 25ul reaction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.786hrze>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol in our group and it is a working PCR protocol.

Created: October 15, 2019



Last Modified: October 15, 2019

Protocol Integer ID: 28670

Keywords: fragment protocol, pcr, reaction, protocol

Materials

Q5 reaction buffer- 5ul

dNTPs- 1ul

Primer mix- 2.5ul

- (1.25ul per primer)

Q5 High Fidelity Polymerase- .25ul

Template DNA- *** varies 1000ng<

- Lambert iGEM: 1.5ul

H2O- to 25ul



- 1 Pipette all reaction agents into pcr tube
- 2 put pcr tube in thermal cycler and set the annealing temperature to 55 degrees celsius
- 3 Annealing temperatures vary depending on DNA
- 4 Set elongation time to 30 seconds
- 5 1 min/2kb for Q5 High Fidelity → may vary