

May 16, 2020

Version 1

# PCR Using Q5U Hot Start High-Fidelity DNA Polymerase (NEB #M0515): *General PCR, USER®Cloning, dUTP incorporation/Carryover prevention* V.1



In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.7schnaw](https://dx.doi.org/10.17504/protocols.io.7schnaw)

New England Biolabs<sup>1</sup>

<sup>1</sup>New England Biolabs

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: [info@neb.com](mailto:info@neb.com)



New England Biolabs

New England Biolabs



## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7schnaw>



**External link:** <https://www.neb.com/protocols/2019/07/02/pcr-using-q5u-hot-start-high-fidelity-dna-polymerase-neb-m0515>

**Protocol Citation:** New England Biolabs 2020. PCR Using Q5U Hot Start High-Fidelity DNA Polymerase (NEB #M0515): General PCR, USER®Cloning, dUTP incorporation/Carryover prevention. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.7schnaw>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 28, 2019

**Last Modified:** March 30, 2021

**Protocol Integer ID:** 28196

**Keywords:** Q5U, USER cloning, dUTP incorporation, novel thermostable dna polymerase, fidelity dna polymerase, using q5u hot start high, mutation in the uracil, carryover prevention q5u hot start high, q5u hot start high, pcr, q5u, uracil, general pcr, dna, containing uracil

## Abstract

Q5U Hot Start High-Fidelity DNA Polymerase is a modified version of Q5® High-Fidelity DNA Polymerase, a novel thermostable DNA polymerase that possesses 3' to 5' exonuclease activity, and is fused to a processivity-enhancing Sso7d domain. Q5U contains a mutation in the uracil-binding pocket that enables the ability to read and amplify templates containing uracil and inosine bases.

# Guidelines

## General Guidelines:

### 1. Template:

Use of high quality, purified DNA templates greatly enhances the success of PCR.

Recommended amounts of DNA template for a 50 µl reaction are as follows:

	DNA	AMOUNT
	DNA Genomic	1 ng – 1 µg
	Plasmid or Viral	1 pg – 1 ng

### 2. Primers:

Oligonucleotide primers are generally 20 – 40 nucleotides in length and ideally have a GC content of 40 – 60 %. Computer programs such as Primer3 can be used to design or analyze primers. The best results are typically seen when using each primer at a final concentration of 0.5 µM in the reaction.

### 3. USER DNA Engineering

Target DNA molecules and cloning vector are generated by PCR with 8 – 12 bases of homology between two fragments. PCR primers start with a 5' A and contain a single deoxyuracil residue (dU) flanking the 3' end of the homology region, and can be designed to accommodate multiple-fragment assembly, nucleotide substitutions, insertions and/or deletions. We recommend using the GeneDesign (<http://genedesign.thruhere.net/gd/>) software to design primers for USER junctions. The best results are typically seen when using each primer at a final concentration of 0.5 µM.

### 4. Mg<sup>++</sup> and additives:

Typically, the Mg<sup>++</sup> concentration for Q5U Hot Start High-Fidelity DNA Polymerase should be 2.0 mM. When used at a final concentration of 1X, the Q5U Reaction Buffer provides this optimal Mg<sup>++</sup> concentration. The addition of common PCR additives such as DMSO may improve amplification of certain difficult or long targets. In these cases, we recommend the addition of up to 2 % DMSO.

### 5. Deoxynucleotides:

The final concentration of dNTPs is typically 200 µM of each deoxynucleotide.

### 6. dUTP Incorporation/Carryover Prevention

Q5U Hot Start High-Fidelity DNA Polymerase is a dUTP-tolerant DNA polymerase that efficiently incorporates dUTP and amplifies uracil-containing substrates. To prevent carryover contamination, dUTP and Antartic Thermolabile UDG (NEB #M0372) can be added to the reaction. dTTP can be fully replaced by dUTP in the amplification of

certain targets. For best results, we recommend adding dUTP at a final concentration of 200  $\mu$ M. For UDG activation, a 10 minute, 25 °C incubation step should be added before the initial denaturation step. Typical cycling parameters can be used thereafter.

### **7. Q5U Hot Start High-Fidelity DNA Polymerase concentration:**

We generally recommend using Q5U Hot Start High-Fidelity DNA Polymerase at a final concentration of 20 units/ml (1.0 unit/50  $\mu$ l reaction). However, the optimal concentration of Q5U Hot Start High-Fidelity DNA Polymerase may vary from 10 – 40 units/ml (0.5 – 2.0 units/50  $\mu$ l reaction) depending on amplicon length and difficulty. It is rarely helpful to exceed 2.0 units/50  $\mu$ l reaction, especially for amplicons longer than 5 kb.

### **8. Buffers:**

The 5X Q5U Reaction Buffer provided with the enzyme is recommended as the first-choice buffer for robust, high-fidelity amplification. The 5X Q5U Reaction Buffer contains 2.0 mM  $Mg^{++}$  at a final (1X) concentration.

### **9. Denaturation:**

*Q5U Hot Start High-Fidelity DNA Polymerase does not require a separate activation step.*

An initial denaturation of 30 seconds at 98 °C is sufficient for most targets being amplified from pure DNA templates. Longer initial denaturation times can be used (up to 3 minutes) for templates that require it. During thermocycling, the denaturation step should be kept to a minimum. Typically, a 5 – 10 second denaturation at 98 °C is recommended for most templates.

### **10. Annealing:**

Optimal annealing temperatures for Q5U Hot Start High-Fidelity DNA Polymerase tend to be higher than for other PCR polymerases. The NEB  $T_m$  Calculator should be used to determine the annealing temperature when using this enzyme. A temperature gradient can also be used to optimize the annealing temperature for each primer pair.

For high  $T_m$  primer pairs, two-step cycling without a separate annealing step can be used (see note 11).

### **11. Extension:**

The recommended extension temperature is 72 °C. Extension times are generally 20 – 30 seconds per kb for complex, genomic samples. Extension time can be increased to 1 minute per kb for long, complex templates, if necessary.

A final extension of 5 minutes at 72 °C is recommended.

### **12. Cycle number:**

Generally, 30–35 cycles yield sufficient product. For genomic amplicons, 30 cycles are recommended.

### **13. 2-step PCR:**

When primers with annealing temperatures  $\geq$  72 °C are used, a 2-step thermocycling protocol (combining annealing and extension into one step) is possible.



#### 14. Amplification of long products:

When amplifying products > 6 kb, it is often helpful to increase the extension time to 1 minute /kb.

#### 15. PCR product:

The PCR products generated using Q5U Hot Start High-Fidelity DNA Polymerase have blunt ends. If cloning is the next step, then blunt-end cloning is recommended. If T/A-cloning is preferred, the DNA should be purified prior to A-addition, as Q5U Hot Start High-Fidelity DNA Polymerase will degrade any overhangs generated.

Addition of an untemplated -dA can be done with *Taq* DNA Polymerase (NEB [#M0267](#)) or Klenow  $\text{exo}^-$  (NEB [#M0212](#)).

## Materials

### MATERIALS

 Q5U® Hot Start High-Fidelity DNA Polymerase [New England Biolabs Catalog #M0515](#)

## Troubleshooting

## Safety warnings

 Please see SDS (Safety Data Sheet) for hazards and safety warnings.

## Before start

Please note that protocols with *Q5U Hot Start High-Fidelity DNA Polymerase* may differ from protocols with other polymerases. Conditions recommended below should be used for optimal performance.

### Reaction Setup:

*Q5U Hot Start High-Fidelity DNA Polymerase is inhibited at room temperature, allowing flexible reaction setup (room temperature or ice).*

**All components should be mixed prior to use.**

## General PCR, USER®Cloning, dUTP incorporation/Carryover prevention

1 Set up the reaction using the following table:

Component	25 µl Reaction	50 µl Reaction	Final Concentration
5X Q5U Reaction Buffer	5 µl	10 µl	1X
10 mM dNTPs	0.5 µl	1 µl	200 µM
10 µM Forward Primer	1.25 µl	2.5 µl	0.5 µM
10 µM Reverse Primer	1.25 µl	2.5 µl	0.5 µM
Template DNA	variable	variable	< 1,000 ng
Q5U Hot Start High-Fidelity DNA Polymerase	0.25 µl	0.5 µl	0.02 U/µl
Nuclease-Free Water	to 25 µl	to 50 µl	

Gently mix the reaction. Collect all liquid to the bottom of the tube by a quick spin if necessary. Overlay the sample with mineral oil if using a PCR machine without a heated lid.

2 Transfer PCR tubes to a PCR machine and begin thermocycling.

### Note

*Q5U Hot Start High-Fidelity DNA Polymerase does not require a separate activation step.*

STEP	TEMP	TIME
Initial Denaturation	98 °C	30 seconds



	30 Cycles	98 °C	5 – 10 seconds
		*55 – 72 °C	20 seconds
		72 °C	20 – 30 seconds/kb
	Final Extension	72 °C	5 minutes
	Hold	4 – 10 °C	

**Thermocycling Conditions for a Routine PCR**

\*Use of the [NEB T<sub>m</sub> Calculator](#) is highly recommended.