

Aug 11, 2025

PEACH: Purifying Elongation complex with Associated factors from Chromatin

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqjpj6ovk5/v1>

 PEACH: Purifying Elongation complex with Associated factors from Chromatin

Yuki Aoi¹, Ali Shilatifard²

¹Simpson Querrey Institute for Epigenetics, Department of Medicine/Biochemistry and Molecular Genetics, Feinberg School of Medicine, Northwestern University, Chicago, IL 60611, USA;

²Simpson Querrey Institute for Epigenetics, Department of Biochemistry and Molecular Genetics, Feinberg School of Medicine, Northwestern University, Chicago, IL 60611, USA



Yuki Aoi

Northwestern University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqjpj6ovk5/v1>

Protocol Citation: Yuki Aoi, Ali Shilatifard 2025. PEACH: Purifying Elongation complex with Associated factors from Chromatin. protocols.io <https://dx.doi.org/10.17504/protocols.io.36wgqjpj6ovk5/v1>

Manuscript citation:

Citation

Aoi Y, Iravani L, Mroczek IC, Gold S, Howard BC, Shilatifard A (2025)
. SPT5 regulates RNA polymerase II stability via Cullin 3-ARMC5 recognition.

<https://doi.org/10.1126/sciadv.adt5885>

LINK

Citation

Aoi Y, Shah AP, Ganesan S, Soliman SHA, Cho BK, Goo YA, Kelleher NL, Shilatifard A (2022)
. SPT6 functions in transcriptional pause/release via PAF1C recruitment.

<https://doi.org/10.1016/j.molcel.2022.06.037>

LINK

Citation

Aoi Y, Takahashi YH, Shah AP, Iwanaszko M, Rendleman EJ, Khan NH, Cho BK, Goo YA, Ganesan S, Kelleher NL, Shilatifard A (2021). SPT5 stabilization of promoter-proximal RNA polymerase II.

<https://doi.org/10.1016/j.molcel.2021.08.006>

LINK

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 08, 2025

Last Modified: August 11, 2025

Protocol Integer ID: 224290



Keywords: nuclei, mammalian, cell culture, immunoprecipitation, IP, chromatin, RNA Polymerase II, Pol II, proteomics, purification, transcription, elongation, purification, benzonase, CTD, PEACH, SPT5, SPT6, NELF, ARMC5, RPB1, biochemistry, antibody, magnetic bead, sucrose cushion, RNAPII, mass spectrometry, IP-MS, PAF1, transcription elongation factor, vivo functions of transcription elongation factor, chromatin solubilization via enzymatic digestion, purifying elongation complex with associated factor, factors from native chromatin, rna polymerase ii elongation complex, following chromatin solubilization, purifying elongation complex, native chromatin, chromatin, proteomic analysis, quality nuclei from mammalian cell culture, rna polymerase ii, nucleic acid, interacting protein, enzymatic digestion, purification, protocol for purification, mammalian cell culture, nuclei

Funders Acknowledgements:

National Cancer Institute

Grant ID: R35-CA197569

Abstract

Here we provide a protocol for purification of RNA Polymerase II elongation complex with associated factors from native chromatin. To isolate high-quality nuclei from mammalian cell culture, adherent cells are lysed directly on a cell culture plate and nuclei are separated from cytoplasmic and other fractions by centrifugation and sucrose cushion. Nuclei can be stored in a freezer. Following chromatin solubilization via enzymatic digestion of nucleic acids, RNA Polymerase II and interacting proteins are enriched using antibody-conjugated magnetic beads. An elution fraction can be used for western blot and proteomic analysis. This method is particularly useful when combined with rapid depletion systems (e.g. auxin-inducible degron/AID, dTAG, PROTAC) to study in vivo functions of transcription elongation factors, including SPT5, SPT6, PAF1, and NELF.



Guidelines

The brief workflow of this protocol:

Day1: Nuclei isolation (2h)

Day2: Purification (~6h)

- Chromatin digestion (2.5h)
- Immunoprecipitation (1.5h)
- Wash and elution (1.5h)

Cell culture




Rapid depletion



Nuclei


Chromatin digestion



Immunoprecipitation of RNA Pol II



Elution





Mass spec

PEACH workflow: A schematic of the workflow of PEACH experiment. Cells (e.g. SPT5-AID) are treated with auxin for rapid depletion of AID-tagged proteins (this rapid depletion step using AID, dTAG, PROTAC, or other systems is optional, but highly recommended.). Nuclei are isolated and can be stored at -80°C .

A chromatin fraction that contains RNA Polymerase II (Pol II) elongation complex is generally insoluble. To solubilize chromatin, benzonase nuclease is used to digest nucleic acid (DNA and RNA). The solubilized chromatin fraction is then used for immunoprecipitation to purify Pol II elongation complex with associated elongation factors. The elution fraction can be used for the subsequent mass spectrometric analysis.

Experimental design:

Cells:

The use of a rapid depletion system is highly recommended: Auxin-inducible degron (AID), dTAG, PROTAC, or other systems can be used to rapidly and specifically degrade a protein of interest within a few hours inside the cells. A minimum of two conditions (e.g. with and without depletion) with four biological replicates are generally required to obtain high-quality and reproducible proteomic data.

Proteomics:

The strategy of proteomic analysis can vary depending on proteomics cores. In our cases at IDeA National Resource for Quantitative Proteomics, briefly, enriched proteins are reduced, alkylated, and digested using filter-aided sample preparation. Peptides are then analyzed on Orbitrap Exploris with data-independent acquisition (DIA). Please see our recent [paper](#) for the detailed protocol regarding our proteomic analysis.



Materials

- ⊗ Benzonase® Nuclease **Merck Millipore (EMD Millipore) Catalog #E1014-25KU**
- ⊗ Protein LoBind tubes **Eppendorf Catalog #022431081**
- ⊗ Dynabeads Protein G for Immunoprecipitation **Thermo Scientific Catalog #10003D**
- ⊗ Ammonium bicarbonate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A6141**
- ⊗ Rpb1 CTD (4H8) Mouse mAb **Cell Signaling Technology Catalog #2629S**
- ⊗ AccuGENE® 1M TRIS-HCl Buffer, pH 7.4, 1 L **Lonza Catalog #51237**
- ⊗ 0.5M EDTA **Invitrogen Catalog #AM9260G**
- ⊗ Potassium Chloride **Fisher Scientific Catalog #P217**
- ⊗ Magnesium Chloride Hexahydrate **Fisher Scientific Catalog #BP214**
- ⊗ D-Sucrose **Fisher Scientific Catalog #BP220-1**
- ⊗ Glycerol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G5516**
- ⊗ Triton™ X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-100ML**
- ⊗ Pierce™ Protease Inhibitor Tablets, EDTA-free **Thermo Fisher Catalog #A32965**
- ⊗ Pierce™ Protease Inhibitor Mini Tablets, EDTA-free **Thermo Fisher Catalog #A32955**
- ⊗ Pierce™ Phosphatase Inhibitor Mini Tablets **Thermo Fisher Catalog #A32957**
- ⊗ Roche DTT 14-Dithiothreitol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #10197777001**
- ⊗ Phosphate buffered saline powder, pH 7.4, for preparing 1 L solutions **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P3813**
- ⊗ Cell scraper, 2-position blade, size: L **Sarstedt Catalog #83.3952**
- ⊗ 5M Sodium Chloride **Invitrogen Catalog #AM9759**
- ⊗ IGEPAL® CA-630 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I8896**
- ⊗ Sodium dodecyl sulfate dust-free pellets **Merck MilliporeSigma (Sigma-Aldrich) Catalog #75746**
- ⊗ Proteinase K **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3115879001**
- ⊗ Gel Loading Dye Purple (6X) - 4.0 ml **New England Biolabs Catalog #B7024S**
- ⊗ Invitrogen™ DynaMag™-2 Magnet **Thermo Fisher Scientific Catalog #12321D**
- ⊗ 2-Mercaptoethanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M3148**

6x SDS loading buffer

2x SDS loading buffer (reducing)



Safety warnings

- ! Sodium dodecyl sulfate (SDS) and Protease/Phosphatase inhibitors are toxic and they may cause skin/eye damage. Care should be taken when handling these reagents.

Before start

Cell culture: Prepare cells in a 15-cm plate at 80-90% confluency.

Buffers:

Lysis buffer (Filter and store at 4°C)

	Stock conc	Volume	Final conc
Tris-HCl pH7.4	1 M	10 mL	10 mM
KCl	2 M	5 mL	10 mM
MgCl ₂	1 M	1.5 mL	1.5 mM
Sucrose	(powder)	120 g	12% (w/v)
Glycerol	100%	100 mL	10% (v/v)
Triton X-100	10%	20 mL	0.2% (v/v)
H ₂ O		Up to 1 L	
<i>Add fresh:</i>			
Protease inhibitor EDTA-free		1 tablet/50 mL	
Phosphatase inhibitor mini		1 tablet/10-50 mL	
DTT	1 M	25 µl/50 ml	0.5 mM

Sucrose Cushion (Filter and store at Room Temperature)

	Stock conc	Volume	Final conc
Tris-HCl pH7.4	1 M	10 ml	10 mM
KCl	2 M	5 ml	10 mM
MgCl ₂	1 M	1.5 ml	1.5 mM
Sucrose	(powder)	300 g	30% (w/v)
H ₂ O		Up to 1 L	
<i>Add fresh:</i>			
DTT	1 M	25 µl/50 ml	0.5 mM

**Freeze buffer** (Filter and store at Room Temperature)

	Stock conc	Volume	Final conc
Tris-HCl pH7.4	1 M	500 µl	10 mM
KCl	2 M	250 µl	10 mM
MgCl ₂	1 M	75 µl	1.5 mM
Glycerol	100%	20 ml	40% (v/v)
H ₂ O		Up to 50 mL	
<i>Add fresh:</i>			
DTT	1 M	1 µl/2 ml	0.5 mM

Chromatin digestion buffer (Filter and store at 4°C)

	Stock conc	Volume	Final conc
Tris-HCl pH7.4	1 M	20 mL	20 mM
NaCl	5 M	30 mL	150 mM
MgCl ₂	1 M	1.5 mL	1.5 mM
Glycerol	100%	100 mL	10% (v/v)
IGEPAL	10 %	5 ml	0.05% (v/v)
H ₂ O		Up to 1 L	
<i>Add fresh:</i>			
Protease inhibitor EDTA-free mini		1 tablet/10 mL	
Phosphatase inhibitor mini		1 tablet/10 mL	
DTT	1 M	5 µl/10 mL	0.5 mM

Mild Elution buffer (Filter and store at Room Temperature)

	Stock conc	Volume	Final conc
SDS	10%	100 µl	0.1% (w/v)
Tris-HCl pH 8.0	1 M	500 µl	50 mM
EDTA	0.5 M	20 µl	1 mM
H ₂ O		Up to 10 mL	

For Day1:

- Keep PBS in a cold room.
 - Cool a centrifuge with swing buckets for 50-ml tubes to 4°C.
 - Prepare buffers ~2h before starting nuclei isolation: Dissolve inhibitor tablets completely on nutator in cold room.
- Add DTT to the buffers. Keep working solutions on ice.



For Day2:

- Cool a centrifuge with a fixed-angle rotor for 1.5-ml tubes to 4°C.
- Preparation of inhibitor-containing chromatin digestion buffer:
 - a. Per 10 ml Chromatin digestion buffer,
 - b. Add 1 Protease inhibitor mini tablet EDTA-free (Thermo, Cat# A32955, 4°C)
 - c. Add 1 Phosphatase inhibitor mini tablet EDTA-free (Thermo, Cat# A32957, 4°C)
 - d. Incubate on nutator at 4°C for around 30 min until tablets get completely dissolved.
 - e. Thaw an aliquot of 1 M DTT (stock at -20°C) and add 5 µl 1 M DTT to Chromatin digestion buffer + inhibitors.Keep on ice.
- Preparation of Benzonase-containing chromatin digestion buffer:
 - a. Per 1 ml Chromatin digestion buffer + inhibitors + DTT,
 - b. Add 4 µl 250 U/µl Benzonase (Sigma, Cat# E1014-25kU, -20°C). Adjust the volume as needed.
 - c. Invert 5 times to mix. Keep on ice.

Note

Once cells are lysed, samples need to be chilled  On ice or at  4 °C to minimize protein degradation during sample preparation.



Day1: Nuclei isolation (2h)

1 Move cell culture (15-cm plates) to a cold room.

2 Pour off media, then add  10 mL ice-cold PBS .

Note

Complete this step for all plates before starting the next step.

3 Wash with  10 mL ice-cold PBS again.

4 After the 2nd wash, tilt plates, then carefully discard the supernatant by aspirating with 10-ml pipette.

5 Add  7 mL Lysis buffer per plate.

Note

Complete this step for all plates before starting the next step.

6 Scrape and transfer cells to a 50-ml tube.

7 Add  2 mL lysis buffer to the plate. Collect and pool lysates. Total lysate =  9 mL

8  00:05:00 Incubation  On ice

9  800 x g, 4°C, 00:05:00 , Swinging bucket

10 Discard the supernatant by aspirating.

**Note**

Do not disturb the pellets.

- 11 Add  5 mL lysis buffer . Carefully resuspend the pellet with a wide-bore P1000 tip.

- 12 Layer nuclei suspension over  25 mL Sucrose Cushion in a 50-ml tube.

Note

Pipette slowly to add cushion under the nuclei suspension.

- 13  800 x g, 4°C, 00:20:00 , swinging bucket .

Note

If a large number of nuclei keep floating after the centrifuge, repeat this step.

- 14 Discard the supernatant.

Note

White debris may be floating in the cushion. Make sure to aspirate them.

- 15 Add  500 µL freeze buffer to the nuclei pellet. Carefully resuspend the pellet with a wide-bore P1000 tip. Transfer nuclei to a 1.5-ml Protein LoBind tube.

- 16 Snap freeze in liquid nitrogen.

Note

PAUSE POINT: Store the nuclei at  -80 °C .



Day2: Chromatin digestion (2.5h)

17 See the "Before start" section above. Prepare the buffers and cool the centrifuge.

18 Thaw frozen nuclei on ice (~20 million nuclei/sample).  On ice

19  500 x g, 4°C, 00:05:00

20 Discard the supernatant using pipette gently.

21 Resuspend in  1 mL ice-cold Benzonase-containing Chromatin digestion buffer .

Note

Use a wide-bore tip to pipet gently.

22 Incubate at  4 °C for  02:00:00 with rotation.

Note

During this incubation step, start beads-Ab preparation (see the section below). Also, heat and stir 6x SDS loading buffer at 80°C to dissolve SDS precipitates completely.

23 Add  2 µL 0.5 M EDTA . Invert to mix.

Note

EDTA concentration can be increased up to 5 mM to inhibit metalloproteases. Note that >1 mM EDTA will also inhibit benzonase.

24  20000 x g, 4°C, 00:05:00 .



- 25 Collect the supernatant in a new 1.5-ml Protein LoBind tube and keep on ice. On ice
Also, keep the insoluble pellet on ice.

Note

Skip step 26 and go to step 27 when opting out the optional recovery from pellet.

- 26 (Optional recovery from pellet)
Resuspend the pellet in
 100 μ L high-salt Chromatin digestion buffer + inhibitors + DTT that contains 500 mM NaCl and 3 mM EDTA.

- 26.1 (Optional recovery from pellet)
Incubate at 4 °C for 00:30:00 with rotation

- 26.2 (Optional recovery from pellet)
Mix with 900 μ L Chromatin digestion buffer + inhibitors + DTT that contains 3 mM EDTA to lower the salt concentration.

- 26.3 (Optional recovery from pellet)
Centrifuge 20000 x g, 4°C, 00:10:00 . Mix the supernatant with the first soluble fraction and keep on ice On ice . Also, keep the insoluble pellet on ice.

- 27 Take 25 μ L sample from the soluble fraction as an input fraction. Keep on ice. On ice

- 28 For WB: Add equal volume (25 μ L) of 2x SDS loading buffer to input soluble sample.

- 28.1 For WB: Resuspend insoluble pellets in 50 μ L 2x SDS loading buffer.

- 28.2 For WB: Boil at 95 °C for 00:03:00 . Then
 20000 x g, Room temperature, 00:01:00 . Store at -20 °C .

- 29 (Optional DNA analysis): Take another 25 μ L input sample from the soluble fraction. Store at -20°C until needed.



- 29.1 (Optional DNA analysis): When analyzing DNA fragment size, add  1.25 μL 20 mg/ml Proteinase K to 25 μL input sample then mix well.
- 29.2 (Optional DNA analysis): Incubate at  55 $^{\circ}\text{C}$ for  02:00:00 .
- 29.3 (Optional DNA analysis): Add  5.25 μL 6x loading buffer, purple . Analyze the DNA fragment size on 1.5% agarose gel with Gel Red.

Beads-Ab preparation (~30min)

- 30 Transfer  50 μL Dynabeads-Protein G to a 1.5-ml Protein LoBind tube.

Note

The volume is per sample; adjust the buffer volume as needed. Beads for up to 10 samples can be collected in a single tube as a bulk.

- 31 Place on a magnetic stand. Discard the supernatant, then resuspend in  1 mL Chromatin digestion buffer (no DTT, no inhibitors, no benzonase) .
 00:01:00 incubation Then place on a magnetic stand and discard the supernatant.

- 32 Resuspend beads in  100 μL Chromatin digestion buffer (no DTT, no inhibitors, no benzonase) .

Note

The volume is per sample; adjust the buffer volume as needed.

- 33 Add  15 μL Rpb1 CTD (4H8) Mouse mAb (CST, Cat# 2629S) .

Note

The volume is per sample; adjust the antibody volume as needed.

- 34 Incubate at  Room temperature for  00:30:00 with rotation .



35 Wash beads twice with
 1 mL Chromatin digestion buffer (no DTT, no inhibitors, no benzonase) .

36 Resuspend in
 50 μ L Chromatin digestion buffer (no DTT, no inhibitors, no benzonase) .

Note

The volume is per sample; adjust the buffer volume as needed.

37 Place on ice until use.  On ice

Immunoprecipitation (1.5h)

38 Complete beads-Ab preparation (see the section above).

39 Add  50 μ L Ab-conjugated beads to each 1 ml soluble chromatin fraction.

Note

When the optional recovery from pellet is performed, add

 50 μ L Ab-conjugated beads to 2 ml pooled soluble chromatin fraction in a 5-ml Protein LoBind tube.

40 Incubate at  4 °C for  01:30:00 with rotation .

Wash and Elution (1.5h)

41 Wash buffer preparation: Transfer
 25 mL Chromatin digestion buffer (no DTT, no inhibitors, no benzonase) to a 50-ml tube and add 12.5 μ l 1 M DTT. Invert to mix, then keep on ice.  On ice

42 Freshly prepare 100 mM ammonium bicarbonate: Dissolve
 316 mg ammonium bicarbonate in  40 mL H₂O . Keep on ice.  On ice



- 43 Place the sample tubes on a magnetic stand.
- 44 [Optional] Transfer 50 µl supernatant to a new 1.5-ml tube as a flow-through fraction.

- 45 Discard the supernatant and resuspend beads in  1 mL ice-cold chromatin digestion buffer (+DTT) . Incubate at  4 °C for  00:01:00 with rotation . Place tubes on a magnet stand, then discard the supernatant.

Note

Do not let the beads dry out during wash steps.

- 46 Repeat wash 3 more times.  [go to step #45](#)

- 47 [Recommended for mass spectrometry] Wash beads with  1 mL freshly prepared, ice-cold 100 mM ammonium bicarbonate .

Note

This wash step will get rid of the detergent from the samples.

- 48 [Recommended for mass spectrometry] Repeat wash again with 100 mM ammonium bicarbonate for a total of two washes.  [go to step #47](#)

- 49 At the last wash, transfer the samples to a new 1.5 ml Protein LoBind tube.

- 50 After the last rotation for washing, briefly spin down then place samples on a magnet stand and bring to the bench at  Room temperature .

- 51 Discard the supernatant carefully. Remove the remaining supernatant using P-10 tip from the bottom of the tube.

**Note**

Do not let the beads dry out during wash steps.

- 52 Resuspend beads in  60 μ L Mild Elution buffer at  Room temperature .

Note

The majority of co-immunoprecipitated proteins are eluted in this buffer. However, the beads partially retain the antibodies and its antigen (RPB1). This helps reduce the amount of co-eluted antibody.

- 53 Incubate at  Room temperature for  00:15:00 with 1,100 rpm on a ThermoMixer.

Note

Make sure to use Protein LoBind tubes and close the lid tightly to avoid leaking. Alternatively, incubate the samples on the bench and flick the tube to mix every 5 minutes.

- 54 Briefly spin down, then place on a magnetic stand. Transfer  60 μ L supernatant to a new 1.5-ml Protein LoBind tube.

- 55 Repeat elution from the beads  [go to step #52](#) then pool the two elution fractions. Briefly vortex to mix a pooled sample and spin down. Total elution volume =  120 μ L .

- 56 Prepare two aliquots of  45 μ L elution fraction in new 1.5-ml Protein LoBind tubes. Snap freeze in liquid nitrogen and store the samples at  -80 $^{\circ}$ C . These samples are ready for sending to a proteomics core for analysis using mass spectrometry.

- 57 For WB/Gel analysis: Prepare 5x SDS loading buffer (+b-ME): Add  100 μ L 2-Mercaptoethanol to  900 μ L 6x SDS loading buffer . Mix by inverting tubes and spin down.

- 57.1 For WB/Gel analysis: Mix  30 μ L elution fraction with  7.5 μ L 5x SDS loading buffer .

- 57.2 For WB/Gel analysis: Resuspend the remaining bead in  50 μ L 2x SDS loading buffer .



- 57.3 For WB/Gel analysis: [Optional] Add equal volume of 2x SDS loading buffer to the flow-through fraction.
- 57.4 For WB/Gel analysis: Boil at 95 °C for 00:03:00 . Then 20000 x g, Room temperature, 00:01:00 . Store at -20 °C .
- 57.5 For WB/Gel analysis: Test each fraction for western blotting and Oriole gel staining analysis to validate purification. When thawing samples for WB/gel analysis, Boil at 95 °C for 00:03:00 . Then 20000 x g, Room temperature, 00:01:00 .

Citations

Aoi Y, Iravani L, Mroczek IC, Gold S, Howard BC, Shilatifard A. SPT5 regulates RNA polymerase II stability via Cullin 3-ARMC5 recognition.

<https://doi.org/10.1126/sciadv.adt5885>

Aoi Y, Shah AP, Ganesan S, Soliman SHA, Cho BK, Goo YA, Kelleher NL, Shilatifard A. SPT6 functions in transcriptional pause/release via PAF1C recruitment.

<https://doi.org/10.1016/j.molcel.2022.06.037>

Aoi Y, Takahashi YH, Shah AP, Iwanaszko M, Rendleman EJ, Khan NH, Cho BK, Goo YA, Ganesan S, Kelleher NL, Shilatifard A. SPT5 stabilization of promoter-proximal RNA polymerase II.

<https://doi.org/10.1016/j.molcel.2021.08.006>

Acknowledgements

We thank the Shilatifard Lab members for contributing to development and improvement of the protocol, I. C. Mroczek for critical reading, and S. G. Mackintosh for proteomic sample preparation advice.