

Jan 09, 2026

PEI transfection protocol for imaging of transient expression of target protein

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn6r7ql5d/v1>

Julia Heiby¹

¹FLI Leibniz Institute on Aging



Julia C. Heiby

FLI Leibniz Institute on Aging

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn6r7ql5d/v1>

Protocol Citation: Julia Heiby 2026. PEI transfection protocol for imaging of transient expression of target protein . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn6r7ql5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 21, 2024



Last Modified: January 09, 2026

Protocol Integer ID: 106099

Keywords: pei transfection protocol for imaging, u2os cells with pei, pei transfection protocol, transient expression of target protein, u2os cell, fused protein, target protein, transient expression of gpf, pei, transient expression, imaging

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://www.protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://www.protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

How to transfect U2OS cells with PEI, for transient expression of GPF-fused protein.



Materials

- U2OS cells (RRID: CVCL_0042)
- 1.5ml tubes (Eppendorf, #0030 120.086)
- 12-well plates (Lab solute, #7696791)
- Coverslips (Carl Roth, #YX03.1)
- Centifuge 5810R (#5811000015, Eppendorf)
- Safety cabinet (Thermo, #Safe S2020 1.2)
- Automated cell counter (Biorad, #TC20)
- Trypan blue solution for cell counting (Sigma, #T8154)
- Counting Slides for cells (Biorad, #145-0011)
- Serologicla Pipettes (Greiner, cellstar)
- 50ml Falcon tubes (Corning, #352070)
- 15ml falcon tube (Corning, #352096)
- Incubator (Thermo, #BBD 6220, CO2 Incubator), settings: 5% CO2, 95%rH, 37°C

- Milli-Q water system (Merck, Advantage A10)
- Poly-D-Lysine (Gibco, #A38904-01)
- PEI (polyethylenimine, MW 25 kDa, Polyscience, #23966)
- DMEM (Sigma, #D6429)
- Opti-MEM (Gibco, #31985-047)
- Formaldehyde (Carl Roth, # CP10.1)
- DAPI (4',6-Diamidino-2-Phenylindole, Dihydrochloride, Thermo Fisher #D1306)
- Permafluor mounting medium (#041300, Menzel)



Day 1

- 1 From a healthy 80% confluent U2OS cells (RRID:CVCL_0042), seed 50 000 cells into a 12-well plate with cover glasses (either $\varnothing$ 12mm)
- 2 Prepare PEI stock solution as the following
 - dissolve 100mg of PEI (polyethylenimine, MW 25 kDa, Polyscience, 23966) in 90ml warm sterile ddH₂O in a glass beaker
 - add HCl dropwise (37%) to pH 2.0 (a few drops should be enough)
 - stir until everything is dissolved (~10min)
 - adjust to pH 7.0 with NaOH (1M)
 - adjust final volume to 100ml
 - filter and store aliquots at -20°C for long term storage
- 3 Prepare the following transfection solution in a 1.5ml tube:
1:3 (μ g DNA: μ g PEI) in 100 μ L total of DMEM or Opti-MEM (without serum and antibiotics!)
 - + 1 μ g DNA (plasmid containing candidate of choice, optionally fused to a fluorescent protein like GFP or RFP)
 - + 3 μ g PEI (prepare 1mg/ml stock solution in sterile ddH₂O)Add up to 100 μ L DMEM (Sigma, D6429) or Opti-MEM (Gibco, 31985-047)
 - 3.1 Vortex quickly
 - 3.2 Incubate transfection solution 15min at RT
- 4 Add the ~100 μ L transfection solution to the well (drop-wise) and gently swirl (depending on the candidate protein, 20-30% of the transfection solution is enough to transfect the cells, reducing the risk of inducing artifacts of overexpression)
- 5 Place cells back to incubator at 37°C

Day 2

- 6 Replace media with complete DMEM media.
(Incubation time with PEI might affect the degree of expression of the candidate protein, and might vary depending on the candidate)



Imaging after 12h or up to 48h

- 7 Time of expression after PEI removal might vary depending on the candidate protein. Usually, after 12h to 48h the cells are ready to be imaged

Fixing cells

- 8
- Wash cells three times with PBS
 - Fix with 4% formaldehyde (v/v) in PBS for 10 min (prepared freshly!)
 - Incubated 5 min with DAPI (4',6-Diamidino-2-Phenylindole, Dihydrochloride, Thermo Fisher D1306, 0.02 µg/µl in PBS) at RT (note: do not perform this step on ice!)
 - Wash 3 times with PBS
 - Mount coverslips in Permafluor mounting medium using glass slides (041300, Menzel)
 - Let dry at room temperature overnight
 - Store samples at 4°C in the dark until further analysis by microscopy
note: If lysotracker was added prior fixation, the fluorescence will decrease over time. Imaging is best after the mounting media has dried.