

Oct 09, 2024

Permanent specimen preparation by protargol staining

DOI

dx.doi.org/10.17504/protocols.io.q26g71or9gwz/v1

Marek Valt¹, Michael Kotyk¹

¹Charles University, Faculty of Science, Department of Zoology, Prague, Czechia

SUM-K



Marek Valt

Charles University, Faculty of Science, Department of Zoolog...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g71or9gwz/v1>

Protocol Citation: Marek Valt, Michael Kotyk 2024. Permanent specimen preparation by protargol staining. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g71or9gwz/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 09, 2024

Last Modified: October 09, 2024

Protocol Integer ID: 109479

Keywords: staining protargol, small flagellate, permanent specimen preparation, protargol, specimen


Abstract

Protargol staining according to the protocol by Bodian (Bodian 1936) and later modified by Nie (Nie 1950), further optimized for staining of small flagellates.

Troubleshooting

Protargol staining optimized for a small flagellates, used for staining of *Solarion arienae*.

2d 1h 15m

- 1 Mix 1 μm of culture sediment with a drop of sterile egg albumin, on a coverslip.
For *Solarion arienae*, egg albumin was not diluted to achieve better adherence to the cells to coverslip.
- 2 Fix in Bouin-Hollande's solution overnight. (0,175 mol/l picric acid, 0,138 mol/l Cu^{+2} acetate, 4% formaldehyde, and 5% acetic acid in water solution) 
- 3 Next day, transfer coverslips throught a graded ethanol series (70%, 50%).
- 4 Wash in distilled water.
- 5 Bleach in 0,5% KMnO_4 solution for 5 minutes. 5m
- 6 Wash in distilled water (5x).
- 7 Treat coverslips with 5% oxalic acid for 5 minutes. 5m
- 8 Wash again in distilled water (5x).
- 9 Stain in 1% protargol solution, for 48 hours at 37°C in a beaker with copper wire pieces placed between the coverslips. 2d

For *Solarion*, protargol produced by Bayer, I. G. Farbenindustrie Actinengesellschaft (out of business since 1952) was used.
- 10 Wash in distilled water (2x).
- 11 Treat with freshly prepared reducing solution (1% hydroquinone and 5% Na_2SO_3), for 10 minutes. 10m



- 12 Wash in distilled water (5x).
- 13 Tone with 1% AuCl_3 for 5 minutes. 5m
- 14 Wash in distilled water (2x).
- 15 Treat with 2% oxalic acid, for 5 minutes. 5m
- 16 Wash in distilled water (5x).
- 17 Treat with 5% $\text{Na}_2\text{S}_2\text{O}_3$ for 10 minutes 10m
- 18 The final wash is done under a constant stream of tap water for 20 minutes. 20m
- 19 Dehydrate the coverslips in an ethanol series (50%, 70%, 80%, 96%, 100%) and 3 times in xylene (5 minutes each).
- 20 Dehydrate in xylene, 3 times (5 minutes each). 15m
- 21 Finally, the stained, dehydrated coverslips are mounted on glass slides with DPX mounting medium (Sigma–Aldrich).