



Jan 05, 2022

Pharyngeal Pumping Assay

DOI

<https://dx.doi.org/10.17504/protocols.io.b3hiqj4e>

Thomas J. O'Brien¹

¹Imperial College London

Behavioural Genomics



Thomas J O'Brien

Imperial College London

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.b3hiqj4e>

Protocol Citation: Thomas J. O'Brien 2022. Pharyngeal Pumping Assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b3hiqj4e>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 05, 2022

Last Modified: January 05, 2022

Protocol Integer ID: 56586

Keywords: pharyngeal pumping assay the pharynx, proper feeding rate in worm, pharyngeal pumping assay, precise timing of pharyngeal movement, pharyngeal movement, pharynx, neuromuscular pump, pumps per minute, proper feeding rate, many times worm, counting grinder movement, pump, worm, grinder movement, movement of the grinder, alimentary tract, neuron, anterior end of the alimentary tract

Abstract

The pharynx is a neuromuscular pump found at the anterior end of the alimentary tract, consisting of 20 muscles and 20 neurons. A proper feeding rate in worms is coordinated by the precise timing of pharyngeal movements, with one complete cycle of synchronous contraction and relaxation of the corpus and terminal bulb termed a "pump". A simple way to measure *C. elegans* feeding is to count how many times worms pump in a minute (pumps per minute). Movement of the grinder (in the terminal bulb) can easily be observed using a stereomicroscope, and because cycles of contraction/relaxation are synchronised along the pharynx, pumps per minute can be measured simply by counting grinder movements.



Materials

Consumables per worm strain:

- 4 × 90 mm petri dish
- 5 × 60 mm petri dish

Reagents for bleaching:

- Sodium hypochlorite, 5% Chlorine (Fisher Scientific, Catalog #419550010)
- 1 M Sodium Hydroxide (Merck Millipore, Catalog #1091371000)

Reagents for preparing NGM agar:

	A	B	C
	Pre-autoclave	For 1L:	For 500mL:
	Sodium chloride	3g	1.5g
	BioAgar	17g	8.5g
	Bacto Peptone	2.5g	1.25g
	Sterile water	975mL	487.5mL

	A	B	C
	Post-autoclave	For 1L:	For: 500mL
	5mg/mL cholesterol (store at 4°C and away from light)	1mL	500uL
	1M CaCl ₂	1mL	500uL
	1M MgSO ₄	1mL	500uL
	1M KPO ₄	25mL	12.5mL



Pick L4 worms for bleaching (Day 1)

- 1 Using an eyelash pick, pick 10 x L4s of each worm strain to be assayed onto 2 × 90 mm plates (pre-seeded with OP50) and incubate at  20 °C

Pour 60 mm NGM agar plates (Day 2)

- 2 Prepare and pour the following number of NGM plates per worm strain to be studied:
 - 5 × 60 mm ( 15 mL agar per plate)
 - 2 × 90 mm ( 35 mL agar per plate)

Therefore make  150 mL NGM per strain

To ensure accuracy when weighing out media components it is recommended that you do not prepare agar in volumes less than 250 mL at a time!

- 2.1 Prepare NGM agar and autoclave

Protocol for preparing NGM: dx.doi.org/10.17504/protocols.io.bvh7n39n

- 2.2 In the laminar flow hood, pour 5 × 60 mm and 2 × 90 mm plates per worm strain

Protocol

NAME

Plate Pouring

CREATED BY

Priota Islam

Preview

- 2.3 Once plates have dried, store agar-side up in an airtight container in the cold room

Dry/Seed 90 mm plates and bleach worms (Day 5)

- 3 Dry 2 × 90 mm plates per worm strain in the cabinet dryer (setting 2) for 1.5 hours



- 3.1 Seed the 90 mm plates with OP50 and leave to dry overnight at room temperature
- 4 Bleach synchronise worms (picked in Step 1) and store eggs/larvae in a 15 mL falcon tube on the rotator that is constantly spinning at  20 °C until feeding

Protocol

NAME

Bleach synchronisation of *C. elegans*

CREATED BY

Ida Barlow

Preview

Refeed L1s (Day 8)

- 5 Refeed bleach synchronised L1s
- 5.1 At 16:50, spin L1s bleached in step 4 using centrifuge program 1 (2500 rpm, 2 mins)
- 5.2 Using a plastic pipette remove the supernatant, taking care to not disturb the worm pellet, and leave ~1 mL of M9 in the tube
- 5.3 Gently flick the falcon tube to resuspend the L1s
- 5.4 Using a glass pipette, drop 1 small droplet around the edges (off food) of the pre-seeded 90 mm plates (prepared in Step 3)

It is recommended that you check the density of worms in the droplet down a microscope and add another if necessary



- 5.5 Place plates agar-side down in the  20 °C incubator to allow droplets to dry (~5-10 mins), then invert plates and allow to grow for 2.5 days

Dry and seed 60 mm plates (Day 10)

- 6 In the morning, take the 60 mm plates from the cold room, and weigh three random plates without their lids
- 7 Dry plates in in cabinet dryer (setting 2) with lids off until they have lost 3-5% of their weight (approx 2 hours)
- 8 Aseptically prepare 1:10 dilution of OP50 in M9 in a 15 mL falcon tube:
 1 mL OP50
 9 mL M9
- 9 Working by a flame, dispense  30 µL of the dilute OP50 into the centre of each 60 mm plate and allow droplets to dry at room temp
- Ensure that you do not touch the agar with the pipette tip and that plates are kept on a flat surface so you get a nice, even circle of food once dried*
- 10 Once the food has dried, leave plates agar-side up at room temp overnight

Assaying rate of pharyngeal pumping (Day 11)

- 11 Label 5 × 60 mm plates seeded with  30 µL OP50 (around the edge, so not to obscure the bacterial lawn) per worm strain
- 12 At 09:00, use an eyelash-pick to transfer 10 x young adult worms (prepared in Step 5) onto the edge of each 60 mm plate (away from food)
- This yields a total of 50 worms per strain*
- 13 Incubate plates, agar-side up, at  20 °C for 1 hour
- This allows worms to crawl to the bacterial lawn and start foraging again*



14 Working one plate at a time, place a 60 mm plate upon a dissecting stereo microscope at the highest magnification possible and bring the worms/bacterial lawn to focus

14.1 Working on a single worm at a time, identify the grinder (in the terminal bulb of the head) of each worm and ensure this is clearly in focus

See:

http://www.wormbook.org/chapters/www_measurepharyngeal/measurepharyngeal.html
for detailed images/diagrams of pharynx and terminal bulb

14.2 Use a tally counter to count the grinder movements of a worm over a 20 second period

Repeat this 3 times per worm and record the result as an average of the three counts

14.3 Go clockwise around the bacterial lawn, and repeat the steps above for each worm on the 60 mm plate. Then repeat for each 60 mm plate remaining in the incubator

15 Calculate pharyngeal pumps per minute of each strain:

pumps per minute = pumps per 20 seconds * 3