

Sep 23, 2025

Version 1

Phenol-chloroform DNA extraction protocol for biodiversity V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6z3rzgqe/v1>

Mateus Chamone Burgarelli¹

¹UFMG

DNA extraction Animals



Mateus Chamone MCB Burgarelli

UFMG

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6z3rzgqe/v1>

Protocol Citation: Mateus Chamone Burgarelli 2025. Phenol-chloroform DNA extraction protocol for biodiversity. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l6z3rzgqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 21, 2025

Last Modified: May 18, 2026

Protocol Integer ID: 227824

Keywords: chloroform dna extraction protocol for biodiversity, chloroform dna extraction protocol, dna extraction, extracting dna, biodiversity sample type, biodiversity, phenol, ml tube

Abstract

This protocol is optimized for DNA extraction in 1.5 mL tubes, but it can be adapted for larger tubes without depending strictly on volume ratios. It was developed for extracting DNA from a wide range of biodiversity sample types.

Before start

Preheat the Lysis Solution in a water bath at 55 °C. Mammalian blood: Preheat Tris NH₄Cl at 37 °C.

DIGESTION

- 1 Transfer a small piece of tissue or blood into 1.5 mL tubes.
Muscle tissue (10–25 mg); Liver (5–10 mg).
Mammalian blood (20–100 μ L); Avian blood (2–10 μ L).
Avian eggshell membrane (25–40 mg – target blood vessels).
Coagulated blood: Heat the sample at **55 °C for 10 minutes**, then immediately place it on ice. Invert or gently shake to disperse the clot.
- 1.1 Rigid tissues: Add Milli-Q H₂O and incubate in a **55 °C water bath for 10–20 min**. Discard the supernatant and proceed to the next step.
- 2 Thoroughly mince the tissue using a micro-scissor and/or pestle.
- 2.1 Mammalian blood: Add 1 mL pre-warmed Tris NH₄Cl (37 °C). Gently mix by inversion, incubate for 5 minutes at room temperature, and centrifuge at 7,000 rpm for 5 minutes. Discard the supernatant and proceed.
- 3 Add 1 mL saline solution (0.85% NaCl).
For blood: Preferably use 1 mL PBS instead of saline.
- 4 Briefly vortex three times and centrifuge at **7,000 rpm for 5 minutes**.
- 5 Discard the supernatant and retain the pellet.
- 6 Repeat steps 3–5.
- 7 Add 400 μ L Lysis Solution 1 (100 mM Tris-HCl, 10 mM EDTA, 100 mM NaCl, 1% SDS), pre-warmed at 55 °C.
- 8 Add 5–15 μ L Proteinase K (20 mg/mL).
- 9 Briefly vortex five times.



- 10 Incubate in a water bath or Thermomixer at **55 °C** with moderate rotation (~550 rpm) for 2–12h, or until complete digestion.
If digestion is incomplete after overnight incubation, add an additional 10 µL Proteinase K and incubate for at least 2 more hours.
Rigid tissues: digestion time may be extended up to 24 h.

CLEAN UP AND PRECIPITATION (ABSOLUTE ETHANOL)

- 11 Briefly vortex the lysed sample.
- 12 Add 10 µL RNase A (10 mg/mL).
- 13 Incubate for 3 minutes at room temperature.
- 14 Add 120 µL Milli-Q H₂O.
- 15 In a fume hood, add 1 volume Phenol:Chloroform:Isoamyl Alcohol (25:24:1).
- 16 Mix vigorously for 30 seconds.
- 17 Centrifuge at 14,000 rpm for 3 minutes.
- 18 Transfer the supernatant to a new tube.
For tissues rich in structural proteins or incomplete lysis: repeat steps 5–8.
- 19 Add 1 volume Chloroform:Isoamyl Alcohol (24:1).
- 20 Mix vigorously for 30 seconds.
- 21 Centrifuge at 14,000 rpm for 3 minutes.



- 22 Transfer 400 μ L of the supernatant to a new tube.
If the aqueous phase is turbid or yellowish, repeat steps 9–12.
- 23 Change gloves.
- 24 Outside the fume hood, add 100 μ L Ammonium Acetate 10 M (final concentration 2 M).
- 25 Add 1 mL ice-cold Absolute Ethanol, then mix well by inversion (20 times).
- 26 Store the sample at **-80 °C for at least 1 h, or overnight at -20 °C**, to stimulate DNA precipitation.
DNA can be indefinitely stored in ethanolic solutions at 0 °C or -20 °C. Invert the tube a few times before proceeding.
- 27 Centrifuge at 14,000 rpm for 20 minutes at 0 °C.
- 28 Carefully discard the supernatant with a pipette without disturbing the pellet.
- 29 Add 1 mL 70% Ethanol and gently invert the tube 3 times.
- 30 Centrifuge at 14,000 rpm for 5 minutes at 4 °C.
Recommended: Repeat steps 18–20.
- 31 Carefully remove the supernatant with a pipette without disturbing the pellet, and let it air dry.
If necessary, use pipettes with progressively smaller volumes until all ethanol is removed.
Do not let the pellet dry completely – it should remain moist before elution.
- 32 Dissolve the pellet in an appropriate volume of Low TE or ultrapure water (20–100 μ L).
Preheating Low TE aliquots at 55 °C is recommended for better elution.
- 33 Centrifuge at 14,000 rpm for 1 minute.
- 34 Incubate in a Water Bath or Thermomixer with minimal rotation at 55 °C for 1–2 h.



To remove ethanol traces, keep tubes open for 10 min at 55 °C.
Alternatively, incubate overnight at room temperature after a brief pre-heating at 55 °C.

- 35 Store the samples in the refrigerator or freezer.

CLEAN-UP AND PRECIPITATION (ISOPROPANOL)

- 36 Briefly vortex the lysed sample.
- 37 Add 10 µL RNase A (10 mg/mL).
- 38 Incubate for 3 minutes at room temperature.
- 39 Add 300 µL Milli-Q H₂O.
- 40 In a fume hood, add 1 volume Phenol:Chloroform:Isoamyl Alcohol (25:24:1).
- 41 Mix vigorously for 30 seconds.
- 42 Centrifuge at 14,000 rpm for 3 minutes.
- 43 Transfer the supernatant to a new tube.
For tissues rich in structural proteins or incomplete lysis: repeat steps 5–8.
- 44 Add 1 volume Chloroform:Isoamyl Alcohol (24:1).
- 45 Mix vigorously for 30 seconds.
- 46 Centrifuge at 14,000 rpm for 3 minutes.

- 47 Transfer the supernatant to a new tube.
If the aqueous phase is turbid or yellowish, repeat steps 9–12.
- 48 Change gloves.
- 49 Outside the fume hood, add Ammonium Acetate 10 M to a final concentration of 2–2.5 M.
Tip: For a final concentration of 2 M, add 1/4 of the recovered volume of 10 M A.A.
- 50 Add 0.7–1 volume of room-temperature Isopropanol, then mix well by inversion (20 times).
- 51 Incubate for 3 minutes at room temperature, inverting occasionally.
- 52 Centrifuge at 14,000 rpm for 25 minutes at 4 °C.
- 53 Carefully discard the supernatant with a pipette without disturbing the pellet.
- 54 Add 1 mL 70% Ethanol at room temperature and gently invert the tube 3 times.
- 55 Centrifuge at 14,000 rpm for 5–10 minutes at 4 °C.
- 56 Repeat steps 18–20.
- 57 Carefully remove the supernatant with a pipette without disturbing the pellet, and let it air dry.
If necessary, use pipettes with progressively smaller volumes until all ethanol is removed.
Do not let the pellet dry completely – it should remain moist before elution.
- 58 Dissolve the pellet in an appropriate volume of Low TE or ultrapure water (20–100 µL).
Preheating Low TE aliquots at 55 °C is recommended for better elution.
- 59 Centrifuge at 14,000 rpm for 1 minute.



- 60 Incubate in a Water Bath or Thermomixer with minimal rotation at 55 °C for 1–2 h.
To remove ethanol traces, keep tubes open for 10 min at 55 °C.
Alternatively, incubate overnight at room temperature after a brief pre-heating at 55 °C.
- 61 Store the samples in the refrigerator or freezer.

Protocol references

Green, Michael R., and Joseph Sambrook. "Isolation of high-molecular-weight DNA using organic solvents." Cold Spring Harbor Protocols 2017.4 (2017): pdb-prot093450.

Green, Michael R., and Joseph Sambrook. "Precipitation of DNA with isopropanol." Cold Spring Harbor Protocols 2017.8 (2017): pdb-prot093385.

Green, Michael R., and Joseph Sambrook. "Precipitation of DNA with ethanol." Cold Spring Harbor Protocols 2016.12 (2016): pdb-prot093377.

Stanzick, Kira J., et al. "DNA extraction from clotted blood in genotyping quality." Biotechniques 74.1 (2023): 23–29.

Ashley Jones, Neeraj Purushotham, Jamila Nasim, Benjamin Schwessinger 2021. DNA clean-up and size selection for long-read sequencing. protocols.io <https://dx.doi.org/10.17504/protocols.io.bwkdcpsc6>

Adair Borges 2023. Virome DNA extraction with phenol-chloroform. protocols.io <https://dx.doi.org/10.17504/protocols.io.ewov1oxdklr2/v1>

<https://www.qiagen.com/us/knowledge-and-support/knowledge-hub/bench-guide/dna/handling-dna/isopropanol-precipitation-of-dna>

<https://www.qiagen.com/us/resources/faq?id=5d591b8b-968a-4a17-849f-9d0f719b40af&lang=en>

Dneasy Blood & Tissue Handbook. Qiagen 2023.

https://www.neb.com/en/protocols/2018/10/23/protocol-for-extraction-and-purification-of-genomic-dna-from-cells-t3010?srsItd=AfmBOorawHfZQkcuQ8uViQS7AQtpFM4VvMNumrRIRRDthMqrQv_4xLL

ChatGPT (<https://chat.openai.com/>), acesso em 24 de março de 2025