



May 27, 2025

pHrodo labeled synaptosome phagocytic assay

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpqnrldzp/v1>

Lucas Blasco-Agell^{1,2}, Meritxell Pons-Espinal^{1,2}, Jara Montero-Muñoz^{1,2}, Antonella Consiglio^{1,2}, Miquel Vila³

¹Department of Pathology and Experimental Therapeutics, Bellvitge University Hospital- IDIBELL, 08908 Hospitalet de Llobregat, Spain;

²Institute of Biomedicine of the University of Barcelona (IBUB), Carrer Baldori Reixac 15-21, Barcelona 08028, Spain;

³VHIR-CIBERNED-ASAP

Vilalab Public



Miquel Vila

VHIR-CIBERNED-ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpqnrldzp/v1>



Protocol Citation: Lucas Blasco-Agell, Meritxell Pons-Espinal, Jara Montero-Muñoz, Antonella Consiglio, Miquel Vila 2025. pHrodo labeled synaptosome phagocytic assay. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gpqnrldzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 26, 2025

Last Modified: May 27, 2025

Protocol Integer ID: 218937

Keywords: synaptosome phagocytic assay phrodo

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020505

Abstract

pHrodo labeled synaptosome phagocytic assay

- 1 Appropriate amount of human synaptosomes (SYNs) are thawed and centrifuged at full speed for 3-4 minutes at 4 °C.
- 2 Pellet is resuspended and well-mixed in cold 100 mM of sodium bicarbonate, pH 8.5 to a concentration of 1 mg/mL.
- 3 SYNs are then incubated with amine-reactive pHrodo Red Succinimidyl Ester (SE) dye (Thermo Scientific) for 60 minutes at RT in a twist shaker (Labnet International).
- 4 Conjugated SYNs are rinsed with DPBS and centrifuged at full speed for 7 times. 10 µg/mL of SYNs were added to hMG plated on a CellCarrierTM-96 Ultra Microplate (Perkin Elmer[®])
- 5 Imaged using 20x objective of a Zeiss SP05 confocal microscopy, maintaining constant levels of temperature (37°C) and CO₂ (5%).
- 6 Four to five images are taken every minute for 5 hours using the phasecontrast and red fluorescence mode.
- 7 6 hours LPS pre-stimulated cells are employed as positive phagocytic controls.
- 8 Phagocytic Index (PI) is quantified using ImageJ, extracting background from the red fluorescent pHrodo.
- 9 Total fluorescence per image is divided by the total number of cells per image.
- 10 PI values for every hour of the time lapse are plot and the area under the curve (AUC) is quantified to perform statistical analysis.