




May 24, 2023

Version 1

PHYTOMap in Arabidopsis root tips V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzbp4xvx1/v1>



Tatsuya Nobori^{1,2}, Joseph Ecker^{1,2}

¹The Salk Institute for Biological Studies; ²Howard Hughes Medical Institute

Nobori Lab @TSL



Tatsuya Nobori

The Sainsbury Laboratory

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzbp4xvx1/v1>

Protocol Citation: Tatsuya Nobori, Joseph Ecker 2023. PHYTOMap in Arabidopsis root tips. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzbp4xvx1/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 28, 2023


Last Modified: May 24, 2023

Protocol Integer ID: 81174

Keywords: In situ hybridization, plant biology, spatial transcriptomics, phytomap in arabidopsis root tip, cell type marker genes in arabidopsis root, arabidopsis root tip, seq datasets in complex plant tissue, targeted observation of gene expression map, spatial analysis of gene expression, spatial mapping of marker gene, gene expression map, complex plant tissue, gene expression, cell rna, marker gene, mount plant tissue, multiplexed fluorescence in situ hybridization method, plant hybridization, phytomap, cell wall digestion enzyme solution, major cell type, rna, single cell, individual cell, gene, cell, complex responses of individual cell

Abstract

Retrieving the complex responses of individual cells in the native three-dimensional tissue context is crucial for a complete understanding of tissue functions. Here, we present PHYTOMap (Plant HYbridization-based Targeted Observation of gene expression Map), a multiplexed fluorescence *in situ* hybridization method that enables single-cell and spatial analysis of gene expression in whole-mount plant tissue in a transgene-free manner and at low cost. We applied PHYTOMap to simultaneously analyze 28 cell type marker genes in *Arabidopsis* roots and successfully identified major cell types, demonstrating that our method can substantially accelerate the spatial mapping of marker genes defined in single-cell RNA-seq datasets in complex plant tissue.



PHYTOMap principles. a, Target mRNA molecules are hybridized by DNA probes (SNAIL probes) that harbor mRNA species-specific barcode sequences (pink bars). The barcode-containing DNA probes are circularized by ligation and amplified *in situ* by rolling circle amplification (RCA). b, Amplified DNA barcodes are detected by sequence-by-hybridization. Different fluorescent probes target four DNA barcodes for each imaging round. After the imaging, fluorescent probes are stripped away, and another set of four genes is targeted.

[Reagents]

DPBST

0.1% (vol/vol) Tween-20 in 1x DPBS

DPBSTR

DPBST + 1:100 SUPERaseIN

Cell wall digestion enzyme solution (CWDES; 10x stock)

250 mg macerozyme, 250 mg cellulase, 500 mg pectinase in 50 mL Nuclease-free water. Filter sterilize (0.22 μ m filter) and store aliquots of 1 mL at -20°C.



Proteinase K buffer

	A	B
	Reagent	Amount
	1M Tris-HCl (pH 8.0)	5 mL
	0.5 M EDTA (pH 8.0)	5 mL
	Nuclease Free Water	up to 50 mL

FAA

	A	B
	Reagent	Amount
	32% formaldehyde	450 µL
	Acetic acid	50 µL
	Ethanol	500 µL

Guidelines

See the following paper for more information on PHYTOMap.

PHYTOMap: Multiplexed single-cell 3D spatial gene expression analysis in plant tissue

Tatsuya Nobori, Marina Oliva, Ryan Lister, and Joseph R. Ecker, bioRxiv, 2022

DOI: [10.1101/2022.07.28.501915](https://doi.org/10.1101/2022.07.28.501915)

- To avoid RNase contamination, special precautions are necessary when handling samples in this protocol. It is advisable to allocate a dedicated area and equipment specifically for RNA work, and to clean them with commercial RNase and DNase inactivating agents before wiping them with ethanol.
- Probes have to be accurately designed. Detailed information on probe design is available in the publication above.

Materials

[Materials]

- Poly-D-Lysine coated dish (MatTek, P35GC-1.5-14-C)
- T4 DNA Ligase (Thermo Scientific, EL0011)
- EquiPhi29 DNA Polymerase (Thermo Scientific, A39391)
- SUPERaseIn RNase Inhibitor (Invitrogen, AM2696)
- Aminoallyl dUTP (AnaSpec, AS-83203)
- Dulbecco's Phosphate Buffered Saline (DPBS) (Sigma, D8662)
- BSA, molecular biology grade (New England Biolabs, B9000S)
- dNTPs (New England Biolabs, N0447S)
- Fluorescent Brightener 28 disodium salt solution (Sigma, 910090)
- Formaldehyde Solution for Molecular Biology, 36.5-38% in Water (Sigma, F8775)
- Triton-X (Sigma, 93443)
- Proteinase K (Invitrogen, 25530049)
- Nuclease-Free Water (Invitrogen, AM9937)
- BS(PEG)9 (Thermo Scientific, 21582)
- 20 ×SSC buffer (Sigma-Aldrich, S6639)
- Ribonucleoside vanadyl complex (RVC) (New England Biolabs, S1402S)
- Formamide (Sigma, F9037)
- Tris, pH 8.0, RNase-free (Invitrogen, AM9855G)
- EDTA, pH 8.0, RNase-free (Invitrogen, AM9260G)
- Gel Slick™ Solution (Lonza, 50640)
- Cellulase (Yaklut, YAKL0013)
- Macerozyme (Yakult, YAKL0021)
- Pectinase (Fisher Scientific, ICN19897901)

[Reagents]

DPBST

0.1% (vol/vol) Tween-20 in 1x DPBS

DPBSTR

DPBST + 1:100 SUPERaseIN

Cell wall digestion enzyme solution (CWDES; 10x stock)

250 mg macerozyme, 250 mg cellulase, 500 mg pectinase in 50 mL Nuclease-free water. Filter sterilize (0.22 µm filter) and store aliquots of 1 mL at -20°C.

Proteinase K buffer



	A	B
	Reagent	Amount
	1M Tris-HCl (pH 8.0)	5 mL
	0.5 M EDTA (pH 8.0)	5 mL
	Nuclease Free Water	up to 50 mL

FAA

	A	B
	Reagent	Amount
	32% formaldehyde	450 µL
	Acetic acid	50 µL
	Ethanol	500 µL

Safety warnings

- ! See safety data sheets for proper chemical handling, precautionary measures and waste disposal. Obey all local regulations/guidelines for handling and disposal of used reagents and solutions containing reagents mixed in.

Formamide:

Handle with proper attire including gloves and eye protection. Work under fume hood when handling solution and dispose of waste appropriately.

Suspected of causing cancer.

May damage fertility or the unborn child.

May cause damage to organs (Blood) through prolonged or repeated exposure if swallowed.

Formaldehyde:

Handle with proper attire including gloves and eye protection. Work under fume hood when handling solution and dispose of waste appropriately.

May cause cancer.

Toxic if swallowed, in contact with skin or if inhaled.

Causes severe skin burns and eye damage.

May cause an allergic skin reaction.

May cause respiratory irritation.

Suspected of causing genetic defects.

Causes damage to organs (Eyes).

Sampling

- 1 Cut Arabidopsis root tips (approximately 1 cm) with razor blade and place them on a Poly-D-Lysine coated dish. The tissue should adhere to the dish well.

2m

Tissue Fixation

- 2 Incubate the tissue in  100 μ L of FAA (see MATERIALS) at  Room temperature for  01:00:00 .

Note

If most of the root tips detach from the dish, it is necessary to optimize the mounting process for the samples.

- 3 Wash the tissue successively for  00:10:00 at  Room temperature with  100 μ L 70% EtOH, 90% EtOH, and twice with 100% EtOH.

10m

- 4 Wash the tissue twice with  100 μ L 100% MeOH for  00:10:00 and leave the tissue in MeOH after the second wash at  -20 $^{\circ}$ C  Overnight .

40m

Tissue Permeabilization

- 5 Rehydrate the tissue by  00:05:00 successive washes with  100 μ L 75%, 50%, 25% MeOH in DPBST.

5m

- 6 Incubate the tissues in  100 μ L 1x CWDES (see MATERIALS)  On ice for  00:05:00 .

5m

	A	B
	Reagent	Amount



	10x CWDES	10 µL
	SUPERase IN	1 µL
	DPBST	89 µL

1x CWDES

- 7 Remove the CWDES and add  100 µL fresh & cold 1x CWDES. Then, incubate the tissue at  Room temperature for  00:30:00 . 30m
- 8 Wash the tissue twice with  100 µL DPBSTR (see MATERIALS).
- 9 Fix the tissue by incubating in  100 µL 10% (v/v) formaldehyde in DPBST at  Room temperature for  00:30:00 . 30m
- 10 Wash the tissue twice with  100 µL DPBSTR at  Room temperature ( 00:01:00 each). 1m
- 11 Incubate the tissue in  100 µL digestion solution at  37 °C for  00:30:00 . 30m

A	B
Reagent	Amount
Proteinase K buffer (see MATERIALS)	99 µL
Proteinase K	1 µL

Digestion solution

- 12 Wash the tissue twice with  100 µL DPBSTR at  Room temperature ( 00:01:00 each). 1m



- 13 Fix the tissue by incubating in  100 μL 10% (v/v) formaldehyde in DPBST at  Room temperature for  00:30:00 . 30m
- 14 Wash the tissue twice with  100 μL DPBSTR at  Room temperature ( 00:05:00 each). 5m

Gene Specific Probe Hybridization

- 15 Mix gene specific probes at the concentration of 5 nM per oligo.
- 16 Heat the probe mixture at  90 $^{\circ}\text{C}$ for  00:03:00 and let it cool down at  Room temperature . 3m
- 17 Incubate the tissue in the hybridization mixture at  40 $^{\circ}\text{C}$ for  03:00:00 or  Overnight 6h

A	B
Reagent	Amount
20xSSC	10 μL
Formamide	30 μL
10% Triton-X	10 μL
200 mM RVC	10 μL
SUPERase IN	1 μL
Probe mix (500 nM per oligo)	2 μL
Nuclease Free Water	37 μL

Hybridization mixture

Note

The design of probes is described in our [paper](#).

18 Wash the tissue twice with  100 µL DPBSTR at  37 °C for  00:30:00 . 30m

19 Wash the tissue twice with  100 µL 4xSSC in DPBSTR at  37 °C for  00:30:00 . 30m

	A	B
	Reagent	Amount
	20xSSC	20 µL
	DPBST	79 µL
	SUPERase IN	1 µL

4xSSC in DPBSTR

20 Rinse the tissue with  100 µL DPBSTR at  Room temperature .

Ligation

21 Incubate the tissue in  100 µL ligation mixture WITHOUT ligase  On ice for  00:05:00 . 5m

	A	B
	Reagent	Amount
	10x ligation buffer	10 µL
	BSA (2mg/ml)	0.5 µL

A	B
SUPERase IN	1 µL
Nuclease Free Water	88.5 µL

Ligation mixture without ligase

- 22 Incubate the tissue in  100 µL ligation mixture WITH ligase at  Room temperature  Overnight .

A	B
Reagent	Amount
10x ligation buffer	10 µL
BSA (20 mg/ml)	0.5 µL
SUPERase IN	1 µL
T4 DNA ligase	2 µL
Nuclease Free Water	86.5 µL

Ligation mixture

Rolling Circle Amplification (RCA)

- 23 Wash the tissue with  100 µL DPBSTR for  00:05:00 .

5m

- 24 Incubate the tissue in  100 µL RCA mixture WITHOUT equiPhi29 DNA polymerase  On ice for  00:05:00 .

5m

A	B
Reagent	Amount
10x equiPhi29 DNA polymerase buffer	10 µL
10 mM dNTP	2.5 µL

A	B
4 mM aminoallyl-dUTP	0.5 µL
BSA (20 mg/ml)	0.5 µL
DTT (100 mM)	1 µL
SUPERase IN	1 µL
Nuclease Free Water	84.5 µL

RCA mixture without equiPhi29 DNA polymerase

- 25 Incubate the tissue in  100 µL RCA mixture WITH equiPhi29 DNA polymerase at  37 °C  Overnight .

5m

A	B
Reagent	Amount
10x equiPhi29 DNA polymerase buffer	10 µL
10 mM dNTP	2.5 µL
4 mM aminoallyl-dUTP	0.5 µL
BSA (20 mg/ml)	0.5 µL
DTT (100 mM)	1 µL
SUPERase IN	1 µL
equiPhi29 DNA polymerase	5 µL
Nuclease Free Water	79.5 µL

RCA mixture with equiPhi29 DNA polymerase

- 26 Wash the tissue twice with  100 µL DPBSTR for  00:10:00 .

10m



Post-Amplification Fixation

- 27 Incubate the tissue in  100 μ L BS(PEG9) solution at  Room temperature for  01:00:00 .

1h

A	B
Reagent	Amount
DPBST	98 μ L
BS(PEG9) stock	2 μ L

BS(PEG9) solution

Note

BS(PEG9) stock is made by adding 465 μ L DMSO into a vial of 100 mg BS(PEG9) and is stored desiccated at -20°C.

- 28 Aspire BS(PEG9) solution and incubate the tissue in  100 μ L 1 M Tris-HCl pH 8.0 at  Room temperature for  00:30:00 .

30m

- 29 Rinse the tissue in  100 μ L DPBST.

Note

If multiple rounds of imaging are not necessary, gel embedding may be skipped. However, it is important to note that without gel embedding, there is a higher risk of tissue detachment or movement between imaging rounds. Therefore, it is recommended to carefully handle the tissue.

Gel Embedding

30 Incubate the tissue in  100 μL monomer solution  On ice for  00:30:00 .

30m

A	B
Reagent	Amount
20% acrylamide	20 μL
2% bis-acrylamide	10 μL
2xSSC	70 μL

Monomer solution

31 Aspirate the monomer solution and add  50 μL gelling solution. Place a Gel Slick-coated glass coverslip on top of the tissue, and carefully aspirate any excess gelling solution. Incubate the tissue at  Room temperature for  01:00:00 —
 02:00:00 until the gel solidifies.

3h

A	B
Reagent	Amount
Monomer solution	48.9 μL
10% APS stock	1 μL
TEMED	0.1 μL

Gelling solution

Note

Coating a coverslip with Gel Slick:
 Add Gel Slick Solution onto the coverslip, then wipe gently with a Kimwipe to spread the Gel Slick Solution.

- 32 Carefully remove the coverslip with forceps, and wash the tissue with  100 μL DPBST for  00:05:00 . 5m
- 33 Incubate the tissue in ClearSee until imaging.

Target Detection with Sequence-By-Hybridization

- 34 Wash the tissue twice with  100 μL 2xSSC for  00:01:00 . 1m
- 35 Incubate the tissue in  100 μL Bridge probe mixture at  Room temperature for  01:00:00 . 1h

A	B
Reagent	Amount
2x hybridization buffer	50 μL
Bridge probes (10uM stock)	1 μL each (typically 4 probes)
Nuclease free water	up to 100 μL

Bridge probe mixture

- 36 Wash the tissue twice with  100 μL 2xSSC for  00:01:00 . 1m
- 37 Incubate the tissue in  100 μL Detection probe mixture at  Room temperature for  01:00:00 . 1h

A	B
Reagent	Amount
2x hybridization buffer	50 μL



A	B
Detection probes (10uM stock)	1 µL each (4 probes)
Calcofluor White	1 µL
Nuclease free water	up to 100 µL

Detection probe mixture

38 Wash the tissue twice with  100 µL 2xSSC for  00:01:00 .

1m

39 Wash the tissue with  100 µL ClearSee, and keep it in ClearSee for imaging.

40 Imaging with a confocal microscope.

41 After the imaging, strip the bridge/detection probes by incubating the tissue in  100 µL stripping solution (65% formamide in 2xSSC) at  30 °C for  00:30:00 .

30m

A	B
Reagent	Amount
Formamide	65 µL
20xSSC	10 µL
Nuclease free water	up to 100 µL

Stripping solution

42 Go to **Step 35** for the next round of imaging.