

Oct 20, 2019

## Pichia pastoris transformation through electroporation

DOI

[dx.doi.org/10.17504/protocols.io.8heht3e](https://dx.doi.org/10.17504/protocols.io.8heht3e)

Igem Dusseldorf<sup>1</sup>

<sup>1</sup>Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.8heht3e](https://dx.doi.org/10.17504/protocols.io.8heht3e)

**Protocol Citation:** Igem Dusseldorf 2019. Pichia pastoris transformation through electroporation. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8heht3e>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 20, 2019

**Last Modified:** October 20, 2019

**Protocol Integer ID:** 28934



## Materials

linearized vector DNA with electrocompetent *Pichia* cells

1 mL of 1:1 1 M sorbitol

YPD Liquid Medium or Plates (500 ml):

- 5 g yeast extract
- 10 g peptone
- 10 agar (if preparing plates)
- Fill to 450 mL with water

→ autoclave

YPDS + antibiotic Plates (1 liter):

- 10 g yeast extract
- 20 g peptone
- 182.2 g sorbitol
- Fill to 800 mL with water
- Mix until dissolved
- Transfer 400 mL two 1 L bottles
- pre-filled with 10 g BactoAgar
- Fill each bottle to 450 mL with water

→ autoclave

- Cool bottles to ~60°C, add 50 mL of 10x D
- Mix well, then aliquot in 4 bottles of 250 mL



## Prior

- 1 *Prior to transformation plasmid DNA containing the gene(s) of interest was linearized in restriction digests using an enzyme that only cuts once following the manufacturers manual !*

## Transformation

- 2 Mix approx. 150 ng of linearized vector DNA with electrocompetent Pichia cells
- 3 Transfer mixture to ice-cold 0.2 cm electroporation cuvette
- 4 Incubate cuvette for 2 min
- 5 Electroporate using a BioRad MicroPulser with a charging voltage of 1.5 kV
- 6 **Immediately** add 1 mL of 1:1 1 M sorbitol and YPD (v/v) to recover the cells
- 7 Transfer mixture to a plastic culture tube and incubated at 30°C for 1 hour at 225 rpm.
- 8 Centrifuge culture tube at 3.000 x g for 5 min to pellet the cells
- 9 Resuspend pellet in 200 µL YPD media
- 10 Plate out cells on YPDS plates containing desired antibiotic and incubated at 30°C for 3 days