


Aug 20, 2025

Picrosirius red staining for histology

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov11oxpvr2/v1>



Samson Dowland¹

¹Charles Perkins Centre Histology Facility, The University of Sydney, NSW, 2006, Australia

The University of Sydney



Samson Dowland

The University of Sydney

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov11oxpvr2/v1>

Protocol Citation: Samson Dowland 2025. Picrosirius red staining for histology. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov11oxpvr2/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 10, 2025

Last Modified: August 20, 2025

Protocol Integer ID: 224421

Keywords: histology, histochemistry, collagen, microscopy, pathology, histopathology, polarised, psr dye stains collagen fibre, histology picrosirius red, collagen in tissue section, iii collagen fibre, using polarised light microscopy, collagen, polarised light microscopy, histochemical technique for the identification, tissue component, picrosirius red, tissue section, fibrosi, tissue remodelling, histochemical technique, immunohistochemistry

Abstract

Picrosirius red (PSR, aka Sirius red) staining is a histochemical technique for the identification of collagen and its specificity makes it suitable for robust quantitative assessments of collagen in tissue sections. Such digital image analysis is particularly valuable for the evaluation of fibrosis and tissue remodelling.

The PSR dye stains collagen fibres red, while also making them birefringent. There is some non-specific red staining of tissue components, which may include nuclei, keratohyalin granules and mucus, however these are not rendered birefringent so do not impact quantification using polarised light microscopy. The non-specific nuclear staining can be prevented by completing the optional heteropolyacid pre-treatment step in this protocol.

Some authors use PSR to distinguish between type I and type III collagen fibres, but recent analysis has shown that this is not always reliable, with immunohistochemistry providing a more definitive distinction.



Protocol materials

⊗ Phosphotungstic acid hydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P4006**

⊗ Ethanol **POCD Scientific Catalog #ETHABS20M**

⊗ Xylene **POCD Scientific Catalog #XYL5P**

⊗ DPX Mountant for histology **Merck MilliporeSigma (Sigma-Aldrich) Catalog #06522**

⊗ Picric acid **May & Baker LTD**

⊗ Sirius red F3BA **BDH Chemicals Ltd Poole England Catalog #34149**

⊗ Iron (III) chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F2877**

⊗ Hydrochloric Acid **Banksia Scientific Company (AJAX FineChem) Catalog #AJA256**

⊗ Haematoxylin **Fronine Laboratory Supplies Catalog #JL555S**

⊗ Ethanol **Hurst Scientific Catalog #EIMS95-5L**

⊗ Ethanol **POCD Scientific Catalog #ETHABS70WV5**

⊗ Glacial acetic acid **Ajax Finechem (ThermoFisher Scientific) Catalog #AJA1-500ML**



Deparaffinisation and rehydration

32m

1 Dewax and rehydrate tissue sections

32m

Heteropolyacid pre-treatment (optional)

3m 30s

2 Incubate in 0.2% phosphotungstic acid (PTA) solution

2m

PTA solution: 0.2 g PTA in 100 ml distilled water



Phosphotungstic acid hydrate **Merck MilliporeSigma (Sigma-Aldrich)** Catalog #P4006

Note

This PTA pre-treatment prevents non-specific nuclear staining in PSR, thereby enhancing the distinction between collagen and cellular staining. Elimination of the non-specific staining is particularly beneficial if performing automated digital image analysis, improving rigour of the method.

Citation

Mukade Y, Kobayashi S, Nishijima Y, Kimura K, Watanabe A, Ikota H, Shirabe K, Yokoo H, Saio M (1970) . Phosphotungstic Acid-treated Picrosirius Red Staining Improves Whole-slide Quantitative Analysis of Collagen in Histological Specimens.

<https://doi.org/10.1369/00221554221141140>

[LINK](#)

3 Rinse with distilled water

30s

4 Rinse with distilled water

30s



5 Rinse with distilled water

30s

Nuclear staining (optional)

20m

6 Stain in Weigert's Haematoxylin (combine solution A & B immediately before use)

10m

Solution A: 1 g Haematoxylin in 100 ml of 95% Ethanol

Solution B: 1.6 g Ferric chloride hexahydrate (Iron III chloride; $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$) + 95 ml distilled water + 1 ml concentrated HCl

 Haematoxylin **Fronine Laboratory Supplies Catalog #JL555S**

 Ethanol **Hurst Scientific Catalog #EIMS95-5L**

 Iron (III) chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F2877**

 Hydrochloric Acid **Banksia Scientific Company (AJAX FineChem) Catalog #AJA256**

7 Wash with running tap water to remove excess stain

5m



8 Differentiate with acid alcohol (5-10 dips) to clear the slides



[1% Hydrochloric Acid (HCl) in 70% Ethanol]

 Hydrochloric Acid **Banksia Scientific Company (AJAX FineChem) Catalog #AJA256**

 Ethanol **POCD Scientific Catalog #ETHABS70WV5**

9 Wash with running tap water

5m



Note

It is important to ensure that the slide is completely cleared of excess Weigert's Haematoxylin at this point. If there is carryover into the PSR solution, it becomes contaminated with the ferric chloride, resulting in non-specific red background staining.



Collagen staining

10 Stain with picrosirius red solution

1h

Solution: 0.1 g Sirius Red F3BA + 100 ml saturated aqueous picric acid

[M] 1.3 Mass Percent in distilled water

⊗ Sirius red F3BA **BDH Chemicals Ltd Poole England Catalog #34149** CI: 35780

⊗ Picric acid **May & Baker LTD**

Sirius red is also known as "Direct Red 80"

Solution can be re-used many times. Lasts 20 years.

Safety information

Picric acid is an unstable explosive when dry, where it becomes highly sensitive to shock, heat and friction. Ensure it remains hydrated at all times.

Picric acid reacts with metals to form unstable picrate salts which are more explosive. Do not use metal containers, lids or racks.

Picric acid is toxic if swallowed, inhaled or absorbed through the skin.

Consult the MSDS and follow all recommended safety precautions.

Differentiation of PSR

30s

11 Differentiate the PSR staining with 2 changes of a 0.5% acetic acid solution, 5-10 dips each.

30s



The timing of this differentiation (or number of dips) should be optimised for each specific tissue. Use a microscope to check that collagen is stained red and cells are stained yellow.

Solution: 0.5 ml glacial acetic acid + 99.5 ml distilled water

⊗ Glacial acetic acid **Ajax Finechem (ThermoFisher Scientific) Catalog #AJA1-500ML**

12 Drain acetic acid solution from the slides as much as possible

Dehydration, clearing and mounting

11m 30s

13 100% ethanol

30s



Note: If you want to maintain the yellow cytoplasmic staining, ensure the dehydration is quick.

10 dips per ethanol step is usually sufficient.

 Ethanol **POCD Scientific Catalog #ETHABS20M**

14 100% ethanol

30s

 Ethanol **POCD Scientific Catalog #ETHABS20M**

15 100% ethanol

30s

Note: ensure whole slide and slide rack is fully immersed in ethanol to remove all traces of water, preventing water contamination of the xylene

 Ethanol **POCD Scientific Catalog #ETHABS20M**

16 100% Xylene

5m

 Xylene **POCD Scientific Catalog #XYL5P**

17 100% Xylene

5m

 Xylene **POCD Scientific Catalog #XYL5P**

18 Coverslip with a resinous mounting medium, such as DPX

 DPX Mountant for histology **Merck MilliporeSigma (Sigma-Aldrich) Catalog #06522**

Results


19

Expected result


Collagen: Red
Bone: Red
Muscle: Yellow
Epithelium: Yellow
Erythrocytes: Yellow
Nuclei: Brown - Black

Note - collagen and bone stained specifically with Sirius Red are birefringent, enabling identification and quantification under polarised light microscopy.


In the below images, the optional heteropolyacid pre-treatment steps were not used.




Skin tissue stained with PSR, without nuclear counterstain. Imaged with regular brightfield microscopy.



Skin tissue stained with PSR, without haematoxylin nuclear counterstain. Imaged with polarised light microscopy. Note that while keratin stains red with PSR, it does not exhibit birefringence as it is not collagenous.



Skin tissue stained with PSR, with Weigert's haematoxylin nuclear counterstain. Imaged with regular brightfield microscopy.



Skin tissue stained with PSR, with Weigert's haematoxylin nuclear counterstain. Imaged with polarised light microscopy. Note that while keratin stains red with PSR, it does not exhibit birefringence as it is not collagenous.

Protocol references

López De Padilla, C. M., Coenen, M. J., Tovar, A., De la Vega, R. E., Evans, C. H., & Müller, S. A. (2021). Picrosirius Red Staining: Revisiting Its Application to the Qualitative and Quantitative Assessment of Collagen Type I and Type III in Tendon. *The journal of histochemistry and cytochemistry : official journal of the Histochemistry Society*, 69(10), 633–643. <https://doi.org/10.1369/00221554211046777>

Citations

Step 2

Mukade Y, Kobayashi S, Nishijima Y, Kimura K, Watanabe A, Ikota H, Shirabe K, Yokoo H, Saio M. Phosphotungstic Acid-treated Picrosirius Red Staining Improves Whole-slide Quantitative Analysis of Collagen in Histological Specimens.

<https://doi.org/10.1369/00221554221141140>