

Mar 29, 2024

Version 2

Piggybac-mediated stable expression of NGN2 in iPSCs for differentiation into excitatory glutamatergic neurons V.2

 [Cell reports](#)

DOI

dx.doi.org/10.17504/protocols.io.e6nvwj54dlmk/v2

Dan Dou^{1,2}, C. Alexander Boecker³, Erika L.F. Holzbaur^{1,2}

¹Department of Physiology, Perelman School of Medicine, University of Pennsylvania, Philadelphia, PA 19104, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, USA;

³Department of Neurology, University Medical Center Goettingen, 37077 Goettingen, Germany

Protocols



Dan Dou

University of Pennsylvania

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwj54dlmk/v2>



External link: <https://doi.org/10.1016/j.celrep.2023.112448>

Protocol Citation: Dan Dou, C. Alexander Boecker, Erika L.F. Holzbaur 2024. Piggybac-mediated stable expression of NGN2 in iPSCs for differentiation into excitatory glutamatergic neurons. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.e6nvwj54dlmk/v2> Version created by **Dan Dou**

Manuscript citation:

Pantazis, C.B., Yang, A., Lara, E., McDonough, J.A., Blauwendraat, C., Peng, L., Oguro, H., Zou, J., Sebesta, D., Pratt, G., et al. (2022). A reference induced pluripotent stem cell line for large-scale collaborative studies. *BioRxiv* 2021.12.15.472643.

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 29, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 97519

Keywords: iPSC, Differentiation, iNeuron, Piggybac, NGN2, ASAPCRN, inducible ngn2 in human ipsc, glutamatergic neuron, stable expression of ngn2, employing piggybac transfection, ineuron differentiation from human ipsc, excitatory glutamatergic neuron, piggybac transfection, inducible ngn2, ngn2, stable integration of ngn2, mediated stable expression, doxycycline, ineuron differentiation

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000350

Abstract

We adapted a previously-described method (Pantazis et al., 2022) for employing Piggybac transfection to stably express doxycycline-inducible NGN2 in human iPSCs. After stable integration of NGN2, proceed to differentiate iPSCs using protocol "iNeuron differentiation from human iPSCs."

Attachments



[549-1145.pdf](#)

106KB



Guidelines

Citations:









- Pantazis, C.B., Yang, A., Lara, E., McDonough, J.A., Blauwendraat, C., Peng, L., Oguro, H., Zou, J., Sebesta, D., Pratt, G., et al. (2022). A reference induced pluripotent stem cell line for large-scale collaborative studies. *BioRxiv* 2021.12.15.472643.

Materials

Materials


- 10 cm cell culture dish
- 6-well cell culture dish
- Cryovials

Reagents

-  Growth Factor Reduced (GFR) Matrigel® **Corning Catalog #354230**
-  Essential 8™ Medium **Gibco - Thermo Fisher Scientific Catalog #A1517001**
-  Accutase® solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A6964**
-  Y-27632 2HCl **Selleckchem Catalog #S1049**
-  Opti-MEM® I Reduced Serum Medium **Thermo Fisher Catalog #31985070**
-  Lipofectamine™ Stem Transfection Reagent **Thermo Fisher Scientific Catalog #STEM00008**
-  PB-TO-hNGN2 **addgene Catalog #172115** [RRID:Addgene_172115](#)
- piggyBac™ transposase vector (Transposagen/Hera BioLabs) **#SPB-D10**
-  KnockOut™ Serum Replacement **Thermo Fisher Catalog #10828010**
- DMSO (CATALOG)

Troubleshooting

Safety warnings

 Wear proper PPE when transferring cryovials to liquid N₂.



Piggybac-mediated stable expression of NGN2 in iPSCs for differentiation into excitatory glutamatergic neurons

3d 6h


- 1 Culture iPSCs in a 10 cm dish coated with Growth Factor Reduced Matrigel (Corning) and feed daily with Essential 8 media (ThermoFisher).
- 2 Passage iPSCs with warm Accutase into Essential 8 media with [M] 10 micromolar (μM) ROCK inhibitor. Plate 800,000 iPSCs into one Matrigel-coated well of a 6-well plate.
- 3 3 - 6 hours after plating, cells should be healthy and attached. Perform transfection using Lipofectamine Stem and a 2:1 ratio of donor plasmid to transposase:

	A	B
	OptiMEM	200 μL
	PB-TO-hNGN2-puro-BFP plasmid	0.75 μg
	EF1 α -transposase plasmid	0.37 μg
	Lipofectamine Stem	4 μL

- 4 Check for transfection efficiency (BFP-labeled cells) on the next day using fluorescence microscopy.
- 4.1 Passage iPSCs with Accutase to a 10 cm dish when cells are confluent enough for splitting.

Note

Continue to feed iPSCs daily with Essential 8 media without ROCK inhibitor, and confirm division of stably-expressing transfected cells (should observe local clusters of BFP-fluorescent cells).

- 5  72:00:00 after transfection, select for transfected iPSCs with [M] 0.5 Mass Percent puromycin.

3d



5.1 Confirm purity of surviving transfected cells with fluorescence microscopy. When population is pure, withdraw puromycin.



6 Cryopreserve selected iPSCs with



	A	B
	Essential 8 media	70%
	Knockout serum replacement	20%
	DMSO	10%
	ROCK inhibitor (Supplement)	10 μ M

6.1 Proceed to culture and induction to neuronal fate using doxycycline (see "Protocol: iNeuron differentiation from human iPSCs").