

Sep 12, 2025

Plant extraction for LC-MS Analysis

 Forked from [Lyophilized plant extraction for LC-MS Analysis](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbkj5qgpk/v1>

Tito Damiani¹, Kateřina Kučerová²

¹IOCB Prague; ²IOCB

Pluskal Lab



Fede B

IOCB

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbkj5qgpk/v1>

Protocol Citation: Tito Damiani, Kateřina Kučerová 2025. Plant extraction for LC-MS Analysis. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.81wgbkj5qgpk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 16, 2025

Last Modified: September 12, 2025

Protocol Integer ID: 220309

Keywords: Plant, Metabolomics, LC-MS, lyophilization, plant extraction for lc, extracting metabolite, plant extraction, plant material for lc, lyophilized plant material, ms analysis, transfer to lc

Abstract

This protocol describes a rapid and efficient method for extracting metabolites from lyophilized plant material for LC-MS analysis. Samples are homogenized, extracted with ethanol–water, dried, and reconstituted in water–acetonitrile before final centrifugation and transfer to LC-MS vials. The method is optimized for small-scale use and ensures compatibility with LC-MS workflows.

Attachments



[Workflow.png](#)

4.4MB



Materials

Solvents

- DMSO
- CH₃CN:H₂O (50:50)

Experimental Supplies

- 2.0 ml Eppendorf tubes
- LC-MS vials
- Stainless Steel grinding beads
- Dry ice

Labels

- For weighted material
- For LC-MS vials

Instruments

- Analytical balance
- Analytical mill
- Thermal mixer
- TissueLyser
- Pipette (10–1000 µL) + tips
- Centrifuge
- Vortex
- LC-MS system



Plant collection, transportation, lyophilisation

- 1 Collect a portion of fresh plant material (leaf, root, stem, flower, fruit) into a 50 mL Falcon tube
- 2 Close the tubes tightly
- 3 Put the tubes in dry ice immediately for transportation to the lab
- 4 Unscrew the tubes
- 5 Cover the top and sides with one piece of thin net fabric
- 6 Attach the fabric to the tube with a rubber band
- 7 Insert the tubes into the lyophiliser flasks
- 8 Run the lyophilisation for 3 days (follow lyophiliser instructions)
- 9 Insert all the lyophilised material into a lab grinder
- 10 Grind the material to obtain a powder
- 11 Store the material in dry room in a tightly closed Falcon tube.

Plant extraction

4h 45m

- 12 Weight **25 ± 2.5 mg** of lyophilized plant material into a 2.0 Eppendorf tube with round bottom

20m



- 13 Add a metal bead to each tube 30s
- 14 Homogenize in TissueLyser **30 sec, 25 rpm** 30s
- 15 Add **2 000 µL** of DMSO (HPLC grade) and vortex to mix 10m
- 16 Incubate **3 min** at **40 °C** 3m
- 17 Homogenise in TissueLyser **60 sec, 25 rpm** 1m
- 18 Centrifuge **5 min, 14 000 rpm** 5m
- 19 Transfer **10 µL** of supernatant into LC-MS vials with insert and add **90 µL** of CH₃CN-H₂O solution (50:50, v/v) 20m