

Jun 12, 2019

(3) Plant Infiltration for *Mimulus in Planta* Transformation



In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.3rtgm6n



¹University of Connecticut

Mimulus



Andrea Sweigart

University of Georgia

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: https://dx.doi.org/10.17504/protocols.io.3rtgm6n

External link:

http://mimubase.org/FTP/Protocols/Stable_Transformation/Mimulus%20in%20planta%20transformation.pdf





Protocol Citation: Yaowu Yuan 2019. Plant Infiltration for Mimulus in Planta Transformation. protocols.io

https://dx.doi.org/10.17504/protocols.io.3rtgm6n

Manuscript citation:

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 04, 2019

Last Modified: June 12, 2019

Protocol Integer ID: 24083

Keywords: Plant infiltration, stable transformation, Mimulus, Mimulus in planta, mimulus in planta transformation, planta transformation this protocol, planta transformation, mimulus

Abstract

This protocol is part of a **collection for** *Mimulus in planta* **transformation**.

Attachments



Mimulus in planta tr...

39KB

Troubleshooting

Safety warnings



See SDS (Safety Data Sheet) for hazards and warnings.

Before start

Ideally want 8-12 young healthy plants in square pots that have branched and have lots of buds for infiltration (When plants are young, cut/pinch the plants so that they will be short and branched).

1 Trim off mature buds from the plants. 2 Gently peel back the leaves around the young buds. 3 Spray the young buds 2-3 times with the Agro (spray right before they will be put in the vacuum). 4 Vacuum infiltrate 1-2 plants at a time. pull a vacuum until 29 inches Hg, seal the vacuum andturn off the pump, wait 2 minutes, then release the vacuum. 4.1 Pull a vacuum until 29 inches Hg. 4.2 Seal the vacuum and turn off the pump. 4.3 Wait 00:02:00 , then release the vacuum. 5 Mark the plant with a colored sword. 6 Cover the plants with plastic tent overnight. 7 Take off plastic and put in tomato cage. 8 Pollinate for 2 weeks. 9 Collect the seeds. 10 Plant in flats using a salt shaker, about 15-20 seed capsules per flat.



- 11 When the seeds start to germinate, begin spraying with 1:1000 Basta.
- 12 Continue spraying with Basta every other day.
- 13 Transplant the healthy looking, resistant plants.