

Jan 31, 2024 Version 2

## Plant RNA purification using TRIzol (TRI reagent) V.2

DOI

[dx.doi.org/10.17504/protocols.io.36wgq4enyvk5/v2](https://dx.doi.org/10.17504/protocols.io.36wgq4enyvk5/v2)

Diep R Ganguly<sup>1</sup>

<sup>1</sup>University of Pennsylvania

Pogson Group



Diep R Ganguly

University of Pennsylvania, The Australian National Universi...

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.36wgq4enyvk5/v2](https://dx.doi.org/10.17504/protocols.io.36wgq4enyvk5/v2)

External link: <https://doi.org/10.1073/pnas.2121362119>

**Protocol Citation:** Diep R Ganguly 2024. Plant RNA purification using TRIzol (TRI reagent). **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.36wgq4enyvk5/v2> Version created by **Diep R Ganguly**

**Manuscript citation:**



**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 30, 2024

**Last Modified:** January 31, 2024

**Protocol Integer ID:** 94427

**Keywords:** plant rna purification, extraction of total rna, using trizol, total rna, rna, dna nuclease treatment, plant tissue, extraction, purification

## Abstract

Extraction of total RNA from plant tissue using TRIzol (or TRI reagent) followed by DNA nuclease treatment.

## Attachments



[TRIreagent\\_SDS.pdf](#)

196KB



[chloroform\\_SDS.pdf](#)

134KB

## Guidelines

TRI reagent and chloroform are hazardous - handle with care, operate in fume hood, wear vinyl gloves and safety glasses.

## Materials

- TRIzol (or TRI reagent)
- Chloroform (or chloroform : IAA [24:1])
- Isopropanol
- 80% Ethanol
- Nuclease-free H<sub>2</sub>O (e.g. 0.01% DEPC-treated H<sub>2</sub>O) or Tris-EDTA (10 mM Tris-Cl, pH 6.5, 0.1 mM EDTA)
- 2 mL safe-lock Eppendorf tubes
- 1.5 mL microcentrifuge tubes
- Tissue lyser or mortar and pestle
- RNase-Free DNase Set (Qiagen)
- RNA loading dye, 2X (NEB)



## Safety warnings

- ! TRI reagent - hazardous
- Chloroform - hazardous

Ensure you read SDS documents (attached) and organise appropriate waste vessels (fume hood).

## Before start

Ensure benches and equipment are RNase free.

For RNase-Free DNase Set (Qiagen): Prepare DNase I stock solution by dissolving the lyophilized DNase I (1500 Kunitz units) in 550 µl of the RNase-free water. Mix gently by inverting the vial. Divide into single use aliquots and store at -20 °C. Thawed aliquots can be stored at 2–8°C for up to 6 weeks. Do not refreeze the aliquots after thawing.



## RNA purification

- 1 Collect 50-100 mg of plant tissue and freeze immediately in liquid N<sub>2</sub>.
- 2 Grind tissue to fine powder under liquid N<sub>2</sub> using tissue lyser or mortar + pestle, then immediately add 1 mL TRI reagent (1 mL per 100 mg tissue).

Note, achieving a fine grind is critical to high yields of intact RNA.
- 3 Invert each tube by hand ~20x and incubate at room temperature for 5 minutes (DO NOT vortex samples as it may result in RNA degradation).
- 4 Add 1/5 volume of pre-mixed chloroform:isoamyl alcohol (24:1), cap tubes, shake vigorously (by hand) for 15 seconds (solution should become cloudy), then incubate at room temperature for 3 minutes.
- 5 Centrifuge at 14,000 rcf for 10 minutes at 4°C.
- 6 Transfer the upper aqueous phase to a new microfuge tube (approx. 400-600 µL).
- 7 Repeat steps 4 and 5 (approx. 300-400 µL).

Note, if you are observing buffer and salt carryover in your purified RNA (high 230 nm absorbance), reduce volume of upper-phase recovered.
- 8 (Optional) If the expected RNA concentration is ≤10 µg/mL, add 1/10 volume of 3M NaOAc (pH 5.5) and/or 100 µg/mL GlycoBlue (or glycogen).
- 9 Add equal volume of 100% isopropanol, then mix by inverting tubes ~20x by hand.
- 10 Incubate at -20°C for 1 hour. Alternatively, incubate overnight to capture small RNAs.
- 11 Centrifuge samples at 14,000 rcf for 20 minutes at 4°C.
- 12 Remove the supernatant, you should observe a white pellet.



- 13 Add 1 mL of 80 % ethanol and invert tube ~10x.
- 14 Centrifuge samples at 10,000 rcf for 5 minutes at room temperature.
- 15 Remove supernatant, carefully since the pellet often becomes dislodged at this step.
- 16 Air-dry pellet at room temperature for 5 minutes.
- 17 Resuspend pellet in RNase-free water (e.g. 0.01% DEPC-treated water) or Tris-EDTA (10 mM Tris-Cl, pH 6.5, 0.1 mM EDTA).

## DNA nuclease treatment and ethanol precipitation

- 18 Make up volume of RNA solution to 87.5  $\mu$ L with nuclease-free water. This can be performed with an aliquot or total sample from the previous step.
- 19 Add 10  $\mu$ L Buffer RDD and 2.5  $\mu$ L DNase I stock solution (Qiagen RNase Free DNase Set) and mix with gentle pipetting.
- 20 Incubate at room temperature for 5-10 minutes.
- 21 Add 500  $\mu$ L of 100% ethanol.
- 22 (Optional) Add 3  $\mu$ L glycogen or GlycoBlue, and 10  $\mu$ L NaOAc (pH 5.5) to aid RNA precipitation.
- 23 Mix by gentle inversion.
- 24 Incubate samples for at least 1 hour at -20  $^{\circ}$ C (overnight, if purifying small RNAs).





- 25 Centrifuge at 14,000 rcf at 4 °C for 20 minutes.
- 26 Remove supernatant and rinse pellet with 1 mL of 80% ethanol.
- 27 Centrifuge at 10,000 rcf for 5 minutes at 4 °C.
- 28 Remove supernatant without disturbing pellet and air-dry for 2 minutes.
- 29 Resuspend pellet in RNase-free water (e.g. 0.01% DEPC-treated water) or Tris-EDTA (10 mM Tris-Cl, pH 6.5, 0.1 mM EDTA).

## Quality control

- 30 Take 50-100 ng aliquot of RNA and mix 1:1 with 2X RNA loading dye (NEB).
- 31 Incubate RNA at 65 °C for 5 minutes.
- 32 Load and run samples on a 1% agarose TBE gel.
- 33 Nanodrop RNA to check for purity and quantity.