

Jan 14, 2026

Plasma isolation for cell-free DNA analysis from dogs using EDTA blood tubes

 [Cancer Research Communications](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov11bzyvr2/v1>

Patricia Filippesen Favaro¹, Muhammed Murtaza¹

¹Center for Precision Medicine, University of Wisconsin-Madison

ctdnalab



Muhammed Murtaza

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov11bzyvr2/v1>

External link: <https://doi.org/10.1158/2767-9764.CRC-25-0373>

Protocol Citation: Patricia Filippesen Favaro, Muhammed Murtaza 2026. Plasma isolation for cell-free DNA analysis from dogs using EDTA blood tubes. protocols.io <https://dx.doi.org/10.17504/protocols.io.ewov11bzyvr2/v1>



Manuscript citation:

Favaro PF, Wang Y, McDonald BR, Wang TX, Jadhav S, Cheng HH, Marcinak CT, Pan X, Murtaza M (2026) Comparative Analysis of Cell-Free DNA Fragmentation Patterns in Canines with Sarcoma and Tumor-Free Canines and Humans. Cancer Research Communications 6(2). doi: [10.1158/2767-9764.CRC-25-0373](https://doi.org/10.1158/2767-9764.CRC-25-0373)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 07, 2025

Last Modified: January 14, 2026

Protocol Integer ID: 231793

Keywords: cell-free DNA, canine, blood, plasma, liquid biopsy, free dna analysis from dog, edta blood tubes laboratory protocol for blood processing, using edta blood tubes laboratory protocol, plasma isolation for cell, tubes laboratory protocol for blood processing, free dna analysis, using blood sample, blood sample, plasma isolation, tubes laboratory protocol, edta tube, blood processing, laboratory, dna

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Laboratory protocol for blood processing and plasma isolation for cell-free DNA analysis from dogs using blood samples collected in EDTA tubes.



Guidelines

1. After blood collection, turn the tube gently 10 times to homogenize the blood and place at room temperature, in upright position, until processing.
2. It is critical that blood samples are centrifuged twice as described in the protocol, within 1-3 hours of collection to minimize the possibility of white blood cells (WBCs) lysis that can affect ability to measure circulating tumor DNA (ctDNA).
3. Collect buffy coat in a separate 2 mL sterile screw-capped DNase-free microcentrifuge tube.
4. Collect plasma aliquots in separate 2 mL sterile screw-capped DNase-free microcentrifuge tubes.

Procedure

- 1 Collect blood in K2-EDTA tubes with 6 mL or 3 mL capacity, based on expected blood volume and size of the dog
Note: The tube should be filled during collection to ensure the correct proportion of additive and blood.
- 2 Gently invert the tube 10 times, and leave it on upright position at room temperature prior to centrifugation.
- 3 Within 1-2 hours, centrifuge the tubes at 820 x g for 10 minutes at room temperature, with the centrifuge brake off (gentle slowdown speed).
- 4 Collect the supernatant (plasma), transferring 1 mL aliquots to microtubes, taking care to not disturb the whitish band (which is the buffy coat).
Label the tubes: #1 #2 #3 (...) in the order of aliquoting the plasma.
Use these in Step 6.
- 5 Collect the buffy coat (~700 μ L) in a separate microtube (cryotube, screw-capped): start collecting from slightly above (~2 mm) the white band.
Note: It is expected to pipette a small amount of red blood cells during buffy coat collection.
Store microtubes with buffy coat at -80°C.
- 6 Centrifuge the plasma aliquots (from step 4) at 16,000 x g for 10 minutes at room temperature, with the centrifuge brake off (gentle slowdown speed).
Record the time that the centrifugation started (centrifugation #2).
Note: a pellet will likely be formed.
- 7 Collect the supernatant (plasma), transferring 900 μ L of each aliquot to a new microtube (cryotube, screw-capped)
Note: be careful to not disturb the pellet.
Store microtubes at -80°C.