

Jan 30, 2025

🌐 Plasmid Cloning - PureLink Kit

DOI

dx.doi.org/10.17504/protocols.io.q26g7m988gwz/v1

Ty Flanagan¹

¹University of Minnesota



Ty Flanagan

University of Minnesota

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7m988gwz/v1>

Protocol Citation: Ty Flanagan 2025. Plasmid Cloning - PureLink Kit. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7m988gwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 29, 2025

Last Modified: January 30, 2025

Protocol Integer ID: 119307

Keywords: purelink quick plasmid miniprep kit, plasmid cloning, purelink kit this protocol, plasmid ext, purelink kit, cloning, goolge drive

Abstract

This protocol is for plasmid cloning using *E. coli* and the PureLink Quick Plasmid Miniprep Kit. See the Plasmid Ext - PureLink Quick Plasmid Miniprep Kit Protocol.pdf on the Goolge Drive.

Troubleshooting



Grow modified *E. coli* overnight

12m

- 1 Prepare working space for bacterial work.
 - 70% ethanol spray down.
 - Bunsen burner.
 - Gloves.
 - Follow aseptic technique.
- 2 Prepare 2 tubes with  2.5 mL LB liquid broth with 50 ug/mL of ampicillin.

5m

3m

Note

Use a non-sealing tube, or at the very end, before incubation, keep the tube cap slightly unscrewed.

The 5 mL non-sealing tubes are recommended.

- 3 Inoculate the media with one "tip-worth" of frozen bacterial stock. Use a 1000 uL pipette tip to physically scrape the frozen stock to collect the bacteria in the pipette tip's tip. Mix the media with the pipette tip.
- 4 Incubate the culture at  37 °C in the shaking incubator overnight.

2m

2m

Note

Keep volumes low relative to tube size as to prevent spilling. Also, make sure that the tube is supported as to not move much.
Do not overtighten screw-on caps.



Follow PureLink Quick Plasmid Miniprep Kit Protocol

58m

- 5 Begin heating  200 µL of TE Buffer in a 1.5 mL tube at  70 °C using a dry bath.
- 6 Centrifuge the cultures and decant the media.
 12000 x g, 00:03:00

5m

5m



- 7 Add  250 μL of Resuspension Buffer (R3) and completely resuspend the bacteria pellet using the pipette. 5m

Note

If using a newer kit, make sure that RNase A has been added to the R3 solution.

- 8 Add  250 μL of Lysis Buffer (L7). Gently mix by inversion 3 times. Incubate at room temp for  00:05:00 . 7m

Note

Do not vortex.

- 9 Add 350 μL of Precipitation Buffer (N4). Quickly after addition, mix by inversion or vigorous shaking. 4m

Note

Do not vortex.

- 10 Centrifuge. 10m

 15000 x g, 00:10:00

- 11 Transfer the supernatant into a spin column with a 2 mL wash tube. Centrifuge and discard the flowthrough. Reuse the same wash tube until noted. 5m

 12000 x g, 00:01:00

- 12 Wash the spin column with  500 μL of Wash Buffer (W10). Incubate at room temp for  00:01:00 . Centrifuge and discard flowthrough. 6m

 12000 x g, 00:01:00

Note

If using a newer kit, make sure that ethanol has been added to the W10 solution.

13 Add  700 μL of Wash Buffer (W9). Centrifuge and discard flowthrough.

6m

 12000 x g, 00:01:00

Centrifuge and place spin column into new 1.5 mL tube.

 12000 x g, 00:01:00

Note

If using a newer kit, make sure that ethanol has been added to the W9 solution.

14 Add  75 μL of the preheated TE Buffer (TE). Incubate at room temp for  00:01:00

2m

15 Centrifuge and discard the column.

3m

 12000 x g, 00:02:00

Note

Store DNA at  4 °C for short term and  -20 °C for long term. Freeze-thaw cycles damage DNA products, but frozen is better for long term storage. ~4 weeks or 30 days is general guidance on shelf life or refrigerator DNA.